**B-1000 Series**

# INSTRUCTION MANUAL

Model
B-1000METBF
B-1000METDK

Ver. 2.5    2023

---

## Table of Contents

1.	Warning	3
2.	Safety Information	3
3.	Package content	4
3.1	B-1000METBF Transmitted and reflected light	4
3.2	B-1000METBF Reflected light	5
3.3	B-1000METDK Transmitted and reflected light	6
3.4	B-1000METDK Reflected light	7
4.	Unpacking	8
5.	Intended use	8
6.	Symbols and conventions	8
7.	Instrument description	9
7.1	B-1000METBF Transmitted and reflected light	9
7.2	B-1000METBF Reflected light	10
7.3	B-1000METDK Transmitted and reflected light	12
7.4	B-1000METDK Reflected light	13
7.5	Motorised version	15
8.	Assembling	16
8.1	B-1000METBF	16
8.2	B-1000METDK	19
8.2.1	Use of nosepiece adapter rings	19
8.3	Motorised version	20
9.	Transmitted light observation procedures	21
10.	Use of the microscope (transmitted light)	22
10.1	General switch on	22
10.2	Control keyboard	22
10.3	Brightness adjustment	22
10.4	Adjust the observation head	22
10.5	Adjust the interpupillary distance	23
10.6	Diopter adjustment	23
10.7	Use of eyeshields	23
10.8	Light path selection	24
10.9	Coarse focus tension adjustment	24
10.10	Focus stop lever	24
10.11	Stage	25
10.12	Centering the condenser	25
10.13	Effect of field diaphragm	25
10.14	Aperture diaphragm	25
10.15	Motorised version	26
10.15.1	Focusing	26
10.15.2	Stage	26
11.	Brightfield reflected light observation procedures	27
12.	Use of the microscope (reflected light)	28
12.1	General switch on	28
12.2	Brightness adjustment	28
12.3	Switching Brightfield/Darkfield observation	28
12.4	Use of ND filter	29
12.5	Changing the stage height	29
12.6	Centering the field diaphragm (FS)	29
12.7	Effects of the field diaphragm	30
12.8	Centering the aperture diaphragm (AS)	30
12.9	Use of color filters (only B-1000METBF)	30
13.	Use of the microscope in Polarized light (PO)	31
13.1	Use of polarized light	31
13.2	Use of connection plate	32
14.	Microphotography	33
14.1	Use of C-mount cameras	33
14.2	Use of reflex cameras	33
15.	Maintenance	34
16.	Troubleshooting	35
	Equipment disposal	37

---

## 1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

## 2. Safety Information

### Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

### 3. Package content

#### 3.1 B-1000METBF Transmitted and reflected light

- ① Frame
- ② Reflected light illuminator
- ③ Observation head
- ④ Eyepieces
- ⑤ Stage
- ⑥ Objectives
- ⑦ Condenser
- ⑧ Polarizer + Analyzer
- ⑨ Color Filters
- ⑩ Power supply
- ⑪ Allen wrench
- ⑫ Dust cover
- ⑬ Power cord

### 3.2 B-1000METBF Reflected light

- ① Frame
- ② Reflected light illuminator
- ③ Observation head
- ④ Eyepieces
- ⑤ Stage
- ⑥ Objectives
- ⑦ Polarizer + Analyzer
- ⑧ Color Filters
- ⑨ Allen wrench
- ⑩ Dust cover
- ⑪ Power cord

### 3.3 B-1000METDK Transmitted and reflected light

- |                               |                               |
|-------------------------------|-------------------------------|
| ① Frame                       | ⑧ Polarizer + Analyzer        |
| ② Reflected light illuminator | ⑨ Power supply                |
| ③ Observation head            | ⑩ Allen wrench                |
| ④ Eyepieces                   | ⑪ Dust cover                  |
| ⑤ Stage                       | ⑫ Power cord                  |
| ⑥ Objectives                  | ⑬ Adapter rings for nosepiece |
| ⑦ Condenser                   |                               |

### 3.4 B-1000METDK Reflected light

- ① Frame
- ② Reflected light illuminator
- ③ Observation head
- ④ Eyepieces
- ⑤ Stage
- ⑥ Objectives
- ⑦ Polarizer + Analyzer
- ⑧ Allen wrench
- ⑨ Dust cover
- ⑩ Power cord
- ⑪ Adapter rings for nosepiece

---

## 4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

## 5. Intended use

### Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

### IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

## 6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

### CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.


### ELECTRICAL SHOCK


This symbol indicates a risk of electrical shock.

## 7. Instrument description


### 7.1 B-1000METBF Transmitted and reflected light

7.2 B-1000METBF Reflected light


**B-1000METBF Both Versions - Opposite side**

7.3 B-1000METDK Transmitted and reflected light


7.4 B-1000METDK Reflected light

**B-1000METDK Both Versions - Opposite side**

## 7.5 Motorised version


Only the parts related to the motorized version are highlighted.

## 8. Assembling

### 8.1 B-1000METBF

1. Insert the reflected light illuminator ① on the frame and lock the locking knob ② with the provided Allen wrench. (Fig. 1)

2. Remove the cover ③ on the back side of the illuminator. (Fig. 2)


3. Insert the LED housing ④ and lock the locking screws ⑤. (Fig. 3-4)

4. Connect the illuminator cable ⑥ to the connector ⑦ placed on the back side of the power supply. (Fig. 5)

5. Connect the power supply to the connector ⑧ placed on the back side of the power supply. (Fig. 6)

6. Insert the optical head above the illuminator and tighten the locking screw ⑨ with the provided Allen wrench. (Fig. 7)

  - **Hold the head with one hand during the locking in order to avoid that the head falls.**


7. Insert both eyepieces into the tubes of the optical head. (Fig. 8)

8. Screw each objective into the thread of the nosepiece, clockwise with increasing magnification. (Fig. 9)

9. Mount the stage: lower the support using the coarse focus knob, then place the stage and firmly tighten the lock screw ⑩. (Fig. 10)

**Only for transmitted and reflected light version:**

10. Insert the condenser under the stage: position until it is well inserted into its holder (under the condenser there is a pin that must fully enter the holder guide). (Fig. 11)
11. Lock the condenser fixing knob ⑪.

12. Insert the power supply jack in the socket placed in the rear side of the of the microscope. (Fig. 12)

## 8.2 B-1000METDK

- The reflected light illuminator for the B-1000METDK is installed prior to shipment from the factory and requires no special installation procedures.

Repeat steps from 2. to 12. of chapter 8.1 for the microscope assembling.

### 8.2.1 Use of nosepiece adapter rings

The nosepiece can allocate objectives for brightfield/darkfield observation (BD). However, it is possible to install objectives for brightfield using the adapter rings.

1. Screw the adapter ring into the empty hole of the nosepiece. (Fig. 13-14)


2. Screw the objective into the adapter ring. (Fig. 15)
3. Perform the observation as usual.

- The brightfield objective cannot be used for darkfield, therefore switching the mirror cube selector to “DF” could bring to a completely dark image.
- **Using the adapter rings it could happen an imperfect parfocality between brightfield and brightfield/darkfield objectives. This is not a defect.**

### 8.3 Motorised version

1. Assemble the stage in the same way as the manual version.
  - Check the perfect alignment of the rear part of the stage with the rear arm of the frame.
  - A non perfect alignment could lead to an incorrect functioning of the system. (Fig. 16)

2. Connect the cable ① from the stage to the frame and tight the locking screws of the connectors ②. (Fig. 17)


3. Connect the provided cables: ③ 12V power supply for the motorised parts; ④ 6V microscope power supply; ⑤ serial cable; ⑥ PS/2 mouse. (Fig. 18)
- **Connect power cables as the last step.**


## 9. Transmitted light observation procedures

## 10. Use of the microscope (transmitted light)

### 10.1 General switch on


To activate the transmitted light illuminator, put the main switch ①, located on the left side of the stand, to the position “1”. (Fig. 19)

### 10.2 Control keyboard

B-1000 transmitted light illumination can be controlled through the keyboard placed on the left of the stand. (Fig. 20)

- **ON-OFF** (②): press this key (after switching the main switch on “1”) to turn ON or OFF the microscope LED.
- **BOOST** (③): press this button in order to increase the brightness (useful for high-magnification objectives or very opaque specimens).
- **Do not enable boost mode while observing with low magnification objectives (4x, 10x) with fully open diaphragm: the high brightness may hurt user’s eyes.**
- **AUTO OFF** (④): if you want the illuminator to switch off automatically, press this button until 15, 30 or 60 minutes delay is set. After this period of time, the light will turn off. You have to press the ON-OFF button to turn it on again.

### 10.3 Brightness adjustment

Use the brightness adjustment dial ⑤ on the left side of the microscope to increase or decrease the light intensity on the specimen. (Fig. 21)

### 10.4 Adjust the observation head

Loosen the locking screw ①, turn the observation head to a comfortable position for observation, and then lock the locking screw again. (Fig. 22)

### 10.5 Adjust the interpupillary distance

Observing with both eyes, hold the two eyepiece prism assemblies. Rotate them around their common axis until the fields of view coincide.

- **The graduation on the interpupillary distance indicator ②, pointed by the spot “.” on the eyepiece holder, shows the distance between the operator’s eyes. (Fig. 23)**

The range of the interpupillary distance is 48-75 mm.

Fig. 23

### 10.6 Diopter adjustment

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
  2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the dioptic adjustment ring ③ to compensate. (Fig. 24)
- **The adjustment range is  $\pm 5$  diopter. The number indicated on the adjustment ring graduation should correspond to the operator’s dioptic correction.**

Fig. 24

### 10.7 Use of eyeshields

- **Use with eyeglasses**

Fold rubber eyeshields with both hands. Folded eyeshields avoid scratching the lenses of eyeglasses. (Fig. 25)

Fig. 25

- **Use without eyeglasses**

Raise eyeshields and observe at the microscope placing eyes to the shields, avoiding external light to disturb the observation. (Fig. 26)

Fig. 26

## 10.8 Light path selection

- The observation head is equipped with an optical path selector that allows the light to be distributed to the eyepieces and to the photo / TV port.
- Move the selector ① to one of the three possible positions to distribute the light. (Fig. 27)

POSITION	LIGHT
IN	100% EYEPIECES
MIDDLE	50% EYEPIECES / 50% TV
OUT	100% TV


Fig. 27

## 10.9 Coarse focus tension adjustment

The tension of the coarse focusing knob is factory preset.

- To modify the tension according to personal's needs, rotate the ring ②. (Fig. 28)
- Clockwise rotation increases the tension.
  - If the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.


Fig. 28

## 10.10 Focus stop lever

The upper limit knob has two functions: prevent the contact between slide and objective and acts as "focus memory".

- After focussing the specimen, pull the lever ① toward the front of the microscope and lock it. (Fig. 29)
- In this way the focus upper limit is set.
- Now one can lower the stage with coarse focus knob, replace the specimen and raise again the stage up to the upper limit: specimen will be in approximate focus and will need a fine adjustment to get the proper focus.
- Fine focus movement is not affected by the coarse focus lock.**
  - To unlock, move the lever in the opposite direction to the one used for the locking.**
- Two stoppers ② are inserted on the frame. DO NOT REMOVE THE TWO STOPPERS.**


Fig. 29

### 10.11 Stage

Stage accepts standard slides 26 x 76 mm, thickness 1,2 mm or metallurgical samples with thickness max. 45 mm, placed on a standard support 26 x 76 mm. (Fig. 30)

- **The height of the stage is pre-set at the factory to accommodate slides. To change the height of the stage, please refer to the specific chapter.**

It is possible to place two slides side by side on the stage.

- **Open the spring arm of the slide holder ① and place from the front the slide on the stage.**
- **Gently release the spring arm of the slide holder.**
- **A sudden release of the spring arm could cause the falling of the slide.**

### 10.12 Centering the condenser


1. Place the specimen on the stage, insert 10x objective into the light path and focus.
2. Insert the front lens of the swing-out condenser ①. (Fig. 31)
3. Rotate the field diaphragm ring ② clockwise, to fully close the diaphragm.
4. Rotate the condenser height adjustment knob ③ to focus the edges of the diaphragm.
5. Rotate the two centering screws ④ to bring the bright spot in the centre of the field of view.
6. Gradually open the diaphragm. The condenser is centred when the diaphragm image is symmetrical to the field of view.
7. In normal use, open the diaphragm until it circumscribes the field of view.

### 10.13 Effect of field diaphragm

Field diaphragm adjusts the illuminated area to obtain a high contrast image.

Set the diaphragm according to the objective in use until it circumscribes the field of view, in order to eliminate unnecessary light to eyepieces. (Fig. 32)

### 10.14 Aperture diaphragm

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image.
- With low contrast specimens set the numerical aperture value ⑤ (printed on the condenser ring) to about 70%-80% of the objective's N.A. (Fig. 33). If necessary, remove one eyepiece and, looking into empty eyepiece sleeve, adjust the condenser's ring in order to obtain an image like the one in Fig. 34.

**Example: with objective PLAN 40x / 0.65 set the scale to 0.65 x 0.8 = 0.52**


Fig. 34

## 10.15 Motorised version

### 10.15.1 Focusing

- Focus motor is activated through the mouse wheel. Front or rear rotation of the mouse wheel raises or lowers the stage. (Fig. 35)
- 1. By moving the mouse wheel without pressing it, the microscope moves in “fine” mode along the Z axis.
- 2. By simultaneously moving and pressing the mouse wheel the microscope moves along the Z-axis in accelerated mode (“coarse” mode), facilitating specimen change or oil positioning.
- **NOTE: Rotations in accelerated mode are ‘discretized’:** a single step of rotation quickly moves the stage along the Z-axis by approximately 4 mm.
- **NOTE: If after the first rotation, you press and turn the mouse wheel again while the stage is moving, there will be no effect. To obtain a second “step” of the stage, you must wait until the first step is completed.**


Fig. 35

### 10.15.2 Stage

1. Stage is moved through the mouse movement. A mouse movement to the front or to the back ③ causes a stage movement of the stage along the Y axis, while a left or right movement ④ causes a stage movement of the stage along the X axis. (Fig. 36)
2. It is always possible, however, operate on the translation knobs of the stage for a manual movement.



Fig. 36

## 11. Brightfield reflected light observation procedures

## 12. Use of the microscope (reflected light)

### 12.1 General switch on

To activate the reflected light illuminator, put the main switch ①, located on the back side of the power supply, to the position "1". (Fig. 37)

### 12.2 Brightness adjustment

Use the brightness adjustment dial ② on the front panel of the reflected light power supply to increase or decrease the light intensity on the specimen. (Fig. 38)

### 12.3 Switching Brightfield/Darkfield observation

- The microscope allows to work in a brightfield (BF) or darkfield (DF).
  - For observation in darkfield it is necessary to use specific objectives called "BD".
1. Move the selector in "BF" position for a brightfield reflected illumination. (Fig. 39)
  2. Move the selector in "DF" position for darkfield reflected illumination.
- Darkfield objectives can also work either brightfield.

## 12.4 Use of ND filter

- Darkfield observation requires a huge amount of light for an optimal observation.
  - In order not to be dazzled when switching from DF to BF position, a neutral density filter is connected to the BF/DF selector. (Fig. 40)
1. The filter is automatically engaged/disengaged simply moving the BF/DF selector.
  - When the selector is in “BF” position the ND filter is engaged.
  - When the selector is in “DF” position the ND filter is disengaged.
  - User can disable this connection by unscrewing the fixing screw on the left side.
2. Using the provided Allen wrench fully unscrew the screw placed in the “LOCK” hole ①. (Fig. 41)
  3. After this operation, the ND filter and the BF/DF selector are disconnected and can be operated separately.

Fig. 40

Fig. 41

## 12.5 Changing the stage height

- In case of thick metallurgical specimen it could happen that the total stroke of the coarse focus knob cannot compensate the height of the specimen. This makes the focusing of the specimen impossible.
  - In this case the stage can be lowered to accept thick specimens allowing a correct focusing range.
1. Loosen the fixing screw ② and lower the stage to the desired position. (Fig. 42)
  2. Lock the fixing screw.
  3. Start working normally.


Fig. 42

## 12.6 Centering the field diaphragm (FS)

1. Move the selector on the reflected light illuminator in the position “BF”. (Fig. 39)
2. Place the specimen on the stage, insert 10x objective into the light path and focus.
3. Completely pull out the field diaphragm selector ③ to fully close the diaphragm. (Fig. 43)
4. Using the provided Allen wrench use the two centering screws ④ to bring the bright spot in the centre of the field of view.
5. Gradually open the diaphragm. The condenser is centred when the diaphragm image is symmetrical to the field of view. (Fig. 44)
6. In normal use, open the diaphragm until it circumscribes the field of view.


Fig. 43

## 12.7 Effects of the field diaphragm

Field diaphragm adjusts the illuminated area to obtain a high contrast image.

Set the diaphragm according to the objective in use until it circumscribes the field of view, in order to eliminate unnecessary light to eyepieces.


Fig. 44

## 12.8 Centering the aperture diaphragm (AS)

1. Move the selector on the reflected light illuminator in the position "BF". (Fig. 39)
2. Completely pull out the aperture diaphragm selector ⑤ to fully close the diaphragm. (Fig. 45)
3. Remove one eyepiece.
4. While looking into the empty eyepiece sleeve, use the provided Allen wrench into the centering screws ⑥ to bring the bright spot in the centre of the field of view. (Fig. 45)
5. The illuminator is centred when the diaphragm image is symmetrical to the field of view.
  - The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image.
  - With low contrast specimens move the diaphragm lever to about 70%-80% of the objective's N.A. If necessary, remove on eyepiece and, looking into empty sleeve, adjust the lever in order to obtain an image like the one in Fig. 34.


Fig. 45

## 12.9 Use of color filters (only B-1000METBF)

1. Insert the desired filter for the observation in the slot placed in the illuminator's side. (Fig. 46)
2. First click stop position leaves filter not inserted, while second click stop position inserts filter into the light path.

FILTER	APPLICATION
Blue	Converts the source colour temperature to daylight
Green	Increase the contrast in black and white photography
Yellow	Contrast filter for the observation of wafers and semiconductors
Frosted white	Reduces lighting irregularities by creating a homogeneous field of view

Fig. 46

## 13. Use of the microscope in Polarized light (PO)

### 13.1 Use of polarized light

1. Move the selector on the reflected light illuminator in the position "BF". (Fig. 39)
  - **Polarized light observation cannot be performed when the BF/DF selector is in "DF" position.**
2. Insert the polarizer "PO" with the inscription facing the user (before remove the empty slider).
3. Move the polarizer "PO" into the light path until the slider clicks into its position. (Fig. 47)

4. Insert analyser "AN" with the inscription facing up (first removing the covering plates on the illuminator).
5. Move the rotatable analyser "AN" into the light path until the slider clicks into its position. (Fig. 48-49)


6. Put a flat mirror on the stage and get into approximate focus the mirror surface.
7. While observing into the eyepieces, rotate the scale of the analyser ① until you achieve the darkest image. (Fig. 50)
  - Now the light extinction (so-called "crossed Nicol position") is achieved.
8. Put a specimen on the stage, focus the specimen and begin observation.

### 13.2 Use of connection plate

- Analyzer and polarizer can be coupled to facilitate the simultaneous insertion of the two sliders.
1. Screw the fixing screw ① into the holes of the analyzer and polarizer to secure the coupling plate ②. (Fig. 51)
  2. Move the assembled unit analyzer + polarizer using the knob ③.


## 14. Microphotography

### 14.1 Use of C-mount cameras

1. Loosen the clamping screw ① on the trinocular port and remove the dust cap ②. (Fig. 52)

2. Screw the C-mount adapter ③ to the camera ④ and insert the round dovetail of the C-mount into the empty hole of the trinocular port, then tighten the clamping screw ①. (Fig. 53)

### 14.2 Use of reflex cameras

1. Insert the Reflex adapter ① into the relay tube to the microscope ②.
  2. Screw the "T2" ring ③ (not provided) to the reflex adapter.
  3. Connect the Reflex camera ④ to the "T2" ring just installed. (Fig. 54)
  4. Mount the other end of the relay tube ② into the empty hole of the trinocular port, then tighten the clamping screw. (Fig. 52)
- "T2" ring is not provided along with the microscope, but is commercially available.
  - While shooting dark specimens, darken eyepieces and viewfinder with a dark cloth to minimize the diffused light.
  - To calculate the magnification of the camera: objective magnification \* camera magnification \* lens magnification.
  - **If using an SLR camera, mirror movement may cause the camera to vibrate.**
  - **We suggest lifting the mirror, using long exposure times and a remote cord.**

---

## 15. Maintenance

### Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

### To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

### Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

### Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

### For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

## 16. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
<b>I. Optical Section:</b>		
LED does not light	Power cord/supply is unplugged	Connect into the power outlet
LED operates, but field of view remains dark	Brightness is too low	Set brightness to a proper level
	DF selector on the reflected light illuminator is inserted but no DF objectives are available	Move the selector to BF position
	Colour filter slider is in a wrong position	Move to a click stop
Field of view is obscured or not evenly illuminated	Light path selector knob is set to the camera position	Move the selector to the eye position
	Nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Condenser is not properly attached	Re-attach it
	Nosepiece is not properly attached	Push the slide dovetail all the way until it is stopped
	Condenser is not properly centred	Centre the condenser
	Field diaphragm is stopped down too far	Open the field diaphragm until it circumscribes the field of view
	DF selector on the reflected light illuminator is inserted but no DF objectives are available	Move the selector to BF position
Dirt or dust is visible in the field of view	Dirt/dust on the specimen	Clean thoroughly
	Dirt/dust on the to surface of the condenser	
	Dirt/dust on the eyepieces	
Image looks double	Aperture diaphragm is stopped down too far	Open aperture diaphragm
	The field diaphragm (both for reflected light and transmitted light) is not well centred	Set the diaphragm according to Koehler settings
Image quality is poor <ul style="list-style-type: none"> <li>• Image is not good;</li> <li>• Contrast is poor;</li> <li>• Details are indistinct;</li> <li>• Image glares</li> </ul>	Nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture diaphragm is too closed or too open	Adjust aperture diaphragm
	For reflected light observation, the specimen must not have any coverglass	Use specimens without coverglass
	Focus is not even	Slide holder is not flat. Move the specimen to a flat position
One side of the image is unfocused	Nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Nosepiece is not correctly mounted	Push slide dovetail all the way until it is stopped
	Stage is not properly attached	Re-attach it
	Specimen is not correctly mounted on the stage	Place specimen correctly on the stage and secure it with slide holder

Image appears to waver	Nosepiece is not corrected mounted	Push slide dovetail all the way until it is stopped
	Objective is not correctly engaged in light path	Make sure that revolving nosepiece clicks into place correctly
	Condenser/reflected field diaphragm are not properly centered	Center the condenser/field diaphragm properly
Field of view becomes only slightly brighter when the voltage is raised	Condenser/reflected field diaphragm are not properly centered.	Center the condenser/field diaphragm properly
	Condenser is lowered too far	Adjust the condenser height position
<b>II. Mechanical Section:</b>		
Coarse adjustment knob is hard to turn	Tension adjustment ring is tightened excessively	Loose ring.
	You are trying to raise stage while focus-lock lever is kept locked	Unlock focus-lock lever
Stage drifts down by itself or focus is lost during observation	Tension adjustment ring is too loose	Tighten ring
Coarse adjustment will not go all the way up	Focus-lock lever is locked at a too low height	Unlock focus-lock lever
Coarse adjustment will not go all the way down	Condenser holder is too low	Raise condenser holder
Image shifts when you touch stage	Stage is not properly mounted	Clamp stage
Specimen stops midway on the X axis movement	Specimen is not correctly positioned	Place specimen correctly
<b>III. Electrical Section:</b>		
LED doesn't turn on	Power supply not connected	Check for proper connection
Brightness is not enough	Brightness setting is too low	Adjust brightness
Light blinks	Power supply not well connected	Check for proper connection
<b>IV. Observation tube:</b>		
Field of view of one eye does not match that of the other	Interpupillary distance is incorrect	Adjust interpupillary distance
	Incorrect diopter adjustment	Adjust diopter
	Your view is not accustomed to microscope observation	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope
<b>V. Microphotography:</b>		
Image edge is unfocused	To a certain extent it is due to achromatic objectives features	To minimize the problem, set the aperture diaphragm in a proper position
Bright spots appear on the image	Stray light entering in the microscope through eyepieces or camera viewfinder	Cover eyepieces and viewfinder with a dark cloth

---

## Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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**OPTIKA® S.r.l.**