ACCESSORIES Series

INSTRUCTION MANUAL

Model

CLD-01

Ver. 2.3 2024

Table of Contents

1.	Warning	3
2.	Safety Information	3
3.	Package content	4
4.	Unpacking	5
5.	Intended use	5
6.	Symbols and conventions	5
7.	Instrument description	6
8.	Assembling	7
8.1	Light guide assembling (CL-11.1 or CL-12)	7
8.2	Connecting the power supply	7
8.3	Connecting polarizing filters	8
9.	Use of the illuminator	9
9.1	CLD-01	9
9.2	CL-11.1	9
10.	Maintenance	10
11.	Troubleshooting	10
	Equipment disposal	11

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

2. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

3. Package content

① LED illuminator

② Power supply + power cord

③ Allen wrench

* CL-11.1 or CL-12 light guides are optional

4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

5. Intended use

Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION


This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.


7. Instrument description


8. Assembling

8.1 Light guide assembling (CL-11.1 or CL-12)

1. Insert the light guide connector ① in the light exit of the illuminator ②. (Fig. 1)

2. Lock the fixing screws ③. (Fig. 2-3)

8.2 Connecting the power supply

1. Plug the power cord on the back of the illuminator. (Fig. 4)

8.3 Installing polarizing filters

1. Insert the polarizing filter on top of the focussing lens and lock the fixing screw. (Fig. 5-6)

2. Adjust the polarizing effect rotating each filter according to the analyzer installed on the microscope head. (Fig. 6)

9. Use of the illuminator

9.1 CLD-01

1. Use ON/OFF switch ① to turn the illuminator on and to adjust the light intensity. (Fig. 7)
- Set the desired light intensity according the specimen under examination.

9.2 CL-11.1

1. Bend the self-supporting arms ② to aim the light beam. (Fig. 8)

2. Each arm can be oriented in a separate way, in order to achieve a proper illumination. (Fig. 9)

10. Maintenance

Device environment

This device is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the device

- The device should always be kept vertically when moving it and be careful so that no moving parts fall out.
- Never mishandle or impose unnecessary force on the device.
- Never attempt to service the device yourself.
- After use, turn off the light immediately, cover the whole microscope with the provided dust cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in “OFF” position.
- Users should observe all safety regulations of the region.
- The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the device to Optika for maintenance, please use the original packaging.


11. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical System		
LED does not light	Power cord is unplugged	Plug power cord into the power outlet
	Main switch is in “OFF” position	Put the main switch to “ON”
Light too dim on the sample	Intensity set too low	Adjust the intensity knob

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.