OPTIGEM Series

INSTRUCTION MANUAL

Model
OPTIGEM-10
OPTIGEM-20

Ver. 1.2 2023

Table of contents

1.	Warning	3
2.	Safety Information	3
3.	Package content	4
3.1	OPTIGEM-10	4
3.2	OPTIGEM-20	5
4.	Unpacking	6
5.	Intended use	6
6.	Symbols and conventions	6
7.	Instrument description	7
7.1	OPTIGEM-10	7
7.2	OPTIGEM-20	9
8.	Assembling	11
9.	Use of the microscope	13
9.1	Adjusting interpupillary distance	13
9.2	Focusing	13
9.3	Adjusting the tension of the focus knob	13
9.4	Diopter compensation	13
9.5	Magnification	14
9.6	Switching on illumination	14
9.7	Use of illumination	15
9.8	Use of the stand in horizontal position	16
10.	Use of optional accessories	17
10.1	Darkfield analysis	17
10.2	Polarized light analysis	18
10.3	Polarized light + Darkfield analysis	18
10.4	Immersion analysis with translating holder	19
10.5	Immersion analysis with translating holder and vacuum pick up	20
10.6	Immersion and polarization analysis with translating holder	20
11.	Microphotography	22
11.1	Installing the C-mount adapter	22
11.2	Use of reflex cameras	22
12.	Maintenance	23
13.	Troubleshooting	24
	Equipment disposal	25

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

2. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

3. Package content

3.1 OPTIGEM-10

① Microscope stand

② Binocular microscope body

③ Eyepieces

④ Jewel clip

⑤ Dust cover

⑥ Power supply

3.2 OPTIGEM-20

- ① Microscope stand
- ② Trinocular microscope body
- ③ Eyepieces

- ④ Jewel clip
- ⑤ Dust cover
- ⑥ Power supply

4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

5. Intended use

Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.


ELECTRICAL SHOCK


This symbol indicates a risk of electrical shock.

7. Instrument description


7.1 OPTIGEM-10


Opposite side

7.2 OPTIGEM-20

Opposite side

8. Assembling

1. Remove the base from its packaging and place it on a flat surface. The base is already assembled from the factory and does not require any further assembly procedure other than that of mounting the microscope body.
2. Insert the microscope body in the holder. (Fig. 1)

3. Lock the fixing knob ①. (Fig. 2)

4. Remove the dust caps from the eyepiece sleeves and insert the eyepieces. (Fig. 3)

5. Connect the jack of the power supply to the socket in the rear part of the microscope base. (Fig. 4)


6. Insert the jewel clip in one of the three threaded holes on the stage. (Fig. 5)

Installing the photo port (OPTIGEM-20)

1. Loosen the fixing screws ① of the supplied photo port and remove the existing photo port. (Fig. 6)
2. Insert the ring of the photo port ② making sure to leave the fixing screw of the photo adapter ③ facing right. (Fig. 7)
3. Lock the fixing screws ①. (Fig. 6)


9. Use of the microscope

9.1 Adjusting interpupillary distance

Hold the right and left eyepiece tube with both hands and adjust the interpupillary distance by moving the two parts until one circle of light can be seen. (Fig. 8)

- If two circles appear, the interpupillary distance is too big.
- If two overlapped circles appear, the interpupillary distance is too small.

9.2 Focusing

Put the sample to be observed on the stage plate and focus the sample using the focusing knobs ①. (Fig. 9)


9.3 Adjusting the tension of the focus knob

- **This adjustment allows to increase or decrease the tension of the knob by avoiding an involuntary descent of the microscope body under its own weight. Adjust the tension just above the point where the focus is stable.**

Grab the knobs ② with both hands and, while holding the left knob, turn the right knob.

- The tension increases or decreases depending on the direction of rotation of the right knob. (Fig. 10)

9.4 Diopter compensation

- **This compensation makes possible for users wearing glasses to adjust the microscope to their eyes and use the microscope without glasses.**

1. Put the zoom down to the lowest magnification using the zoom knob ③ and focus the specimen with the focusing knobs ④. (Fig. 11)
2. Put the zoom to the maximum magnification and repeat the focusing.
3. Return to the lowest magnification: the specimen will be out of focus.

4. Adjust the diopter compensation ring of the right eyepiece ⑤ (Fig. 12) until the image of the right eyepiece is clear and sharp. Repeat the procedure for the left eyepiece.
5. Check the focus of the image for the whole zoom range. It should be perfectly parfocal (focus is maintained during the change of magnification).

9.5 Magnification

Select the desired magnification by rotating the zoom knob.

- The microscope body is equipped with a “magnification lock” function that allows to obtain a precise setting of the desired magnification. Magnification lock can be activated or deactivated by operating with the provided Allen wrench in the hole ⑥ placed in front of the microscope body. (Fig. 13)

9.6 Switching on illumination

1. Set the main switch ① (located on the left side of the base) to the “I” position to power the base. (Fig. 14)

2. Turn the knob ② (on the right side of the base) to switch the lighting on/off and to increase/decrease the light intensity. (Fig. 15)

9.7 Use of illumination

The OPTIGEM microscope comes with an integrated LED illumination system that consists of three separate light sources:


1. Incident illumination: through a flexible arm attached to the microscope stage.
2. Transmitted illumination: from the hole of the microscope stage.
3. Darkfield illumination: a LED ring inside the stage hole.

By turning the switch ③ on the left side of the base (Fig. 16), the different observation methods can be selected.


The switch has four available position (Fig. 17):

- a) Transmitted light illumination
- b) Darkfield illumination
- c) Darkfield illumination combined with Incident illumination
- d) Incident light illumination

The transmitted light illumination has an additional option: the slider ④ (Fig. 18) allows observation in bright field (a) or “velvet” dark field (b) (Fig. 19) by simply moving the slider to the left or to the right.

9.8 Use of the stand in horizontal position

The OPTIGEM microscope stand can easily be tilted in order to observe objects in a horizontal position.

1. Loosen the head fixing knob (Fig. 2) and rotate the head by 180°.
2. Fix the screw again.
3. Loosen the stand tilting knob ① (Fig. 20) and rotate the stand by 180° (make sure that the power cable has not been twisted too many turns around the base).
4. Lean the stand backwards until it gets stopped by the support ② and the microscope is in the position of figure 21.

10. Use of optional accessories

A wide range of accessories are available for the OPTIGEM microscope: gemology clips, an iris aperture diaphragm, a polarization analysis kit, an optical unit to switch the dark field from a “sharp” (gem exhibition) to a “soft” (diamond analysis) mode and a translating cell holder with vacuum pick-up and quartz immersion cells complete and enhance the instrument usability.

All the accessories can be individually used or combined in a stack, in any order, both when the microscope is being used in horizontal position or when it is in vertical position.

Below some possible configurations for common uses are illustrated.

10.1 Darkfield analysis

- **ST-201 is needed.**


1. Place rods and spacers as shown in Fig. 22.

2. Put the ST-201 on the microscope stage, aligning the holes in the stage with the holes in the ST-201.
3. Put one rod from the bottom of the stage, as shown. (Fig. 23)

4. Screw one spacer on the rod, from the top (Fig. 24). Repeat the same steps for the other rod.

5. Switch on the LED ring illuminator (position “b”, Fig. 17).
 - If you want it is possible to put your gem directly on the diaphragm, after closing it to a proper diameter. (Fig. 25)


Fig. 25

10.2 Polarized light analysis

- **ST-202 is needed.**
1. Repeat steps 1-4 of chapter 10.1 for the installation of ST-202. (Fig. 26)

Fig. 26

2. Screw the analyzer on the bottom of the microscope body. (Fig. 27)
3. Holding your sample with the jewel clip, you can rotate both the polarizer and analyzer to perform polarization analysis.

Fig. 27

10.3 Polarized light + Darkfield analysis

- **ST-201 + ST-202 are needed.**
1. Install ST-201 as described in chapter 10.1.
 2. Screw another spacer on both rods.
 3. Stack ST-202 and fix it with another spacer. (Fig. 28)
 4. Screw the analyzer on the bottom of the head. (Fig. 27)
 5. Holding your sample with the jewel clip, you can rotate both the polarizer and analyzer to perform polarization analysis.



Fig. 28


10.4 Immersion analysis with translating holder

- ST-203 + ST-204 + ST-207 are needed.


1. Install ST-207 as described in chapter 10.1. (Fig. 29)


2. Stack translating holder ST-204 and fix it with another spacer. (Fig. 30)


3. Rotate the microscope into horizontal position. (Fig. 31)

4. Put the immersion cell (ST-203) on the holder's plate. (Fig. 32)


5. Install the jewel clip on the vertical rod, in order to hold the sample inside the immersion cell. (Fig. 33)


10.5 Immersion analysis with translating holder and vacuum pick up

- **ST-203 + ST-204 + ST-207 + ST-205 are needed.**
1. Install the accessories as described in chapter 10.4.
 2. Insert the adapter for the vacuum pick-up pen on the vertical rod. (Fig. 34)
-
3. Switch on the vacuum pump by pressing the black button. (Fig. 35)
 4. Pick up the gem with vacuum by placing it in contact with the rubber point of the pick-up pen.


10.6 Immersion and polarization analysis with translating holder

- **ST-202 + ST-203 + ST-204 are needed.**
1. Install ST-202 as described in chapter 10.1.
 2. Screw the analyzer on the bottom of the head. (Fig. 27)
 3. Stack the translating holder and fix it with the spacers. (Fig. 36)


-
4. Now rotate the microscope into horizontal position and put the immersion cell ST-203 on the holder's plate. (Fig. 37)

11. Microphotography

11.1 Installing the C-mount adapter

1. Loosen the clamping screw ① on the trinocular port and remove the dust cap ②. (Fig. 38)

2. Screw the C-mount adapter ③ to the camera ④ and insert the round dovetail of the C-mount into the empty hole of the trinocular port (Fig. 39), then tighten the clamping screw ①.

11.2 Use of reflex cameras

1. Insert the Reflex adapter ① into the relay tube ②.
 2. Screw the "T2" ring ③ (not provided) to the reflex adapter.
 3. Connect the Reflex camera ④ to the "T2" ring just installed. (Fig. 40)
 4. Mount the other end of the relay tube ② into the empty hole of the trinocular port, then tighten the clamping screw. (Fig. 38)
- "T2" ring is not provided along with the microscope, but is commercially available.
 - While shooting dark specimens, darken eyepieces and viewfinder with a dark cloth to minimize the diffused light.
 - To calculate the magnification of the camera: objective magnification * camera magnification * lens magnification.
- **When using an SLR camera, mirror movement may cause the camera to vibrate.**
 - **We suggest lifting the mirror, using long exposure times and a remote cord.**

12. Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 75 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off- position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.


13. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
The illumination is ON, but the field of view is dark.	The plug is not connected to the illumination	Connect the cable
	The brightness is too low	Adjust to a proper setting
The edge of the field of view is vignetted or the brightness is asymmetric.	The incident illuminator is not correctly oriented	Change the angle of the incident illuminator
Dust and stains can be seen in the field of view.	There are stains and dust on the specimen	Clean the specimen
	There are stains and dust on the eyepiece	Clean the eyepiece
Poor image quality: • The image is not sharp • The contrast is not high • The details are not clear • Image glares	The lenses (additional lens, objective, eyepieces) are dirt	Thoroughly clean all the optical system
One side of the image is out of focus.	The specimen is out of place (tilted)	Place the specimen flat on the stage.
II. Mechanical Section:		
The focus knob is hard to turn.	The tension adjustment collar is too tight	Loosen the tension adjustment collar
The focus is unstable.	The tension adjustment collar is too loose	Tighten the tension adjustment collar
III. Electric Section:		
The LED doesn't turn on.	No power supply	Check the power cord connection
The brightness is not enough	The brightness adjustment is low	Adjust the brightness
The light blinks	The power cord is poorly connected	Check the power cord
IV. Viewing tube assembly:		
The field of view of the two eyes is different	The interpupillary distance is not correct	Adjust the interpupillary distance
	The diopter correction is not right	Adjust the diopter correction
	The viewing technique is not correct, and the operator is straining the eyesight	When look into the eyepieces, do not stare at the specimen but look at the whole field of view. Periodically, move the eyes away to look at a distant object, then back into the eyepieces
V. Microphotography and video:		
The image is unfocused	Incorrect focusing	Adjusting the focus system as in the present manual
Bright patches appear on the image	Stray light is entering the microscope through the eyepieces and through the camera viewfinder	Cover the eyepieces and the viewfinder with a dark cloth

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.