OPTIKA

M I C R O S C O P E S

I T A L Y

Ver. 1.1.0

OPTIGEM 3-4

 **OPERATION MANUAL**

 **GUIDA UTENTE**

OPTIKA MICROSCOPES - ITALY

www.optikamicroscopes.com - info@optikamicroscopes.com

1.0	INTRODUCTION	page 3
2.0	GENERAL CHARACTERISTICS	page 3
3.0	SAFETY GUIDELINES	page 4
3.1	Electrical safety	
4.0	WARNING/DANGER SIGNALS USED IN THIS MANUAL	page 4
5.0	THE MICROSCOPE AND ITS COMPONENTS	page 5
6.0	UNPACKING AND ASSEMBLY	page 6
6.1	Environmental operating conditions	
6.2	Unpacking	
6.3	Installation	
6.4	Installation of the photo/video adapter (for OPTIGEM-4)	
7.0	USING THE MICROSCOPE	page 11
7.1	Illumination system	
7.2	Adjustment of the stereo head	
8.0	MAINTENANCE AND CLEANING	page 15
9.0	TROUBLESHOOTING	page 18
10.0	ACCESSORIES	page 19
11.0	DIMENSIONS	page 19
12.0	RECOVERY AND RECYCLING	page 20

This microscope is a precision scientific instrument that has been developed for a professional use. It has been built according to high optical and mechanical standards to withstand many years of daily laboratory use.

This is a classical gemmological microscope with tilting stand. The darkfield condenser employs a typical optical configuration based on an halogen illumination and a reflection system. A flexible arm with fluorescent tube simulates the daylight colour temperature in order to provide an accurate colour evaluation.

The advanced illumination system, using a combination of transmitted and incident light, allows a thorough analysis of the object surface. The microscope stand can be inclined up to 45° relative to the vertical line, to offer a higher user comfort. The base can be rotated 360°, which also provides versatility.

The stereoscopic head is equipped with a continuous zoom optical system (from 0.67x up to 4.5x) that guarantees an excellent field depth. The trinocular port on the head (OPTIGEM-4) offers the possibility to connect a photo/video camera to view and elaborate images on a PC.

2.0 GENERAL CHARACTERISTICS

- ✓ Zoom Head 0.67x up to 4.5x (binocular and trinocular).
- ✓ Eyepieces: WF10x/22mm, high-point
- ✓ 360° Head rotation
- ✓ Trans-illumination: bright and darkfield, adjustable intensity, both corresponding to the sun colour temperature (6000 K).
- ✓ Iris diaphragm assembled on the darkfield condenser
- ✓ Epi-illumination: fluorescent tube mounted on a flexible arm
- ✓ Clamps for Gemmological objects
- ✓ 45°-tilting stand and 360° rotating base
- ✓ Models: OPTIGEM-3 with binocular head;
 OPTIGEM-4 with trinocular head.

3.0 SAFETY GUIDELINES

This manual contains important information on safety concerning the installation, use and maintenance of the gemmological Optika microscopes OPTIGEM-3 and OPTIGEM-4. We recommend to read carefully this manual before using this instrument. For a safe use, it's essential to read and abide by the instructions set out in this manual.

The OPTIKA products have been designed for a safe use in normal operating conditions. The microscope and accessories described hereby have been made and tested according to industrial safety standards for laboratory instruments.

An incorrect use can cause injuries or damages to the instrument.

Keep this manual always ready for a quick consultation whenever necessary.

3.1 Electrical safety

- Before connecting the power supply, make sure that the line voltage in your geographical area matches the operating voltage of this equipment, and that the lamp switch is in the off position.
- Users should observe all safety regulations of their country. This equipment is supplied with CE safety label. However, users do have full responsibility for using it safely.

4.0 WARNING/DANGER SIGNALS USED IN THIS MANUAL

Users need to be aware of safety issues when using this microscope. The warning/danger signals – which have been used in this manual - are shown below.

DANGER

Follow the instructions denoted by this symbol to avoid potential serious risks for people.

WARNING

Warning sign, indicating that an incorrect use of the instrument can cause injuries or damages

WARNING

Possibility of an electric shock

HOT!

Warning: surface with high temperatures. Avoid direct contact.

NOTE

Technical notes or use recommendations.

The picture below shows the OPTIGEM-3 and OPTIGEM-4 microscope along with its components:

This microscope must be used in a dust-free environment. Keep it away from solvents, chemical vapours and humidity. Do not expose it to high temperature variations, directly to the sun or to vibrations, which all could jeopardize the instrument's performance.

6.1 Operating environmental conditions

- Temperature: 10 - 36°C (50 - 96.8°F)
- Relative humidity: 0 - 80% up to 30°C (86°F)

6.2 Unpacking

Check the packing to make sure that nothing is missing. Make a list of all accessories so as to have a quick reference when ordering spare parts or contacting the customer support. No accessories or other parts have to be left inside the packing. Store the original packing away, to be used for future shipments of the microscope and/or its accessories.

Do not touch optical surfaces such as lenses, filters or glasses. Grease traces or other residual substances may compromise the quality of the final image and quickly corrupt the optical surfaces.

6.3 Installation

- Remove all contents from the original packing and place it in a safe and dry place.
- Place the main body of the microscope on a flat and stable surface.
- **Do not connect the power supply at this stage.**


- Remove the front plug which protects the objectives of the optical head

- Insert the head in the focus holder and tighten the screw on the left side



- Remove the plugs from the eyepiece slots and insert the eyepieces. The slots are equipped with a safety screw to fix the eyepieces.

6.0 UNPACKING AND ASSEMBLY

- Insert the clip for holding gemological objects by pressing it down into one of the two dedicated holes.
- Connect the power supply to the socket located on the base.


The power supply cable must be connected only to sockets equipped with a suitable earthing.

- Let a skilled technician check your electrical system to ensure it is suited.
- If other accessories do not have to be installed, the microscope is now ready for use.

6.4 Installation of the photo/video adapter (for OPTIGEM-4)

- Remove the protection cover from the Photo/video output.


C-MOUNT Adapter

Use the following adapters to connect cameras with C-Mount standard connection:

- ST-147 Adapter for cameras (1/3" sensors)
- ST-147.1 Adapter for cameras (1/2" sensors)


- Screw the adapter part equipped with optics on the head photo/video output.

- Screw the adapter ring on the camera C-mount connection.

- Place the camera on the top of the adapter and tighten the fixing screw.

- Check the camera user guide for further information.

REFLEX Adapter


Use the following adapter to connect SRL cameras:

ST-146 Adapter for REFLEX camera

- Insert the photographic eyepiece in the lower part of the adapter.

- Screw the adapter on the head photo/video output.


- Screw the other adapter extremity on the SRL camera T2 ring – which can be provided by any supplier of photographic material.

6.0 UNPACKING AND ASSEMBLY

- Place the camera on the top of the adapter and tighten the fixing screw.

7.0 USING THE MICROSCOPE

After installing the necessary components, the microscope is ready for use.

OPTIGEM is a laboratory instrument designed to last for a long time. Always handle it carefully and avoid vibrations, hits and abrupt movements.

Always disconnect the power supply when the microscope is inactive for a long time, when cleaning it or when performing maintenance operations.



DO NOT DISASSEMBLE THE INSTRUMENT

Do not disassemble this instrument. This would invalidate the warranty and may cause failures.


7.1 Illumination system

- Push the main switch button located sideways on the base.
- Rotate the potentiometer for transmitted light adjustment till a suited intensity level has been reached.
- Thanks to the knob located on the object plate, a black plate can be inserted in order to allow a dark-field observation of the sample.
- The Iris diaphragm allows an accurate control of the transmitted light with the darkfield. Adjust its opening by means of the dedicated lever.

- The incident light is activated by the red button located on the lamp. Adjust the directional arm to find the proper illumination angle.

7.2 Adjustment of the stereo head

Interpupillary distance

Hold the eyepiece tubes with both hands and adjust the distance between them until a single circle can be observed (while observing with both eyes). If two separate circles appear, reduce the interpupillary distance; if two overlapping circles appear, increase such distance.

Focusing and tension adjustment

Hold the sample with the clamp and set the zoom at the minimum magnification. Proceed to focusing the sample image by means of the pair of lateral knobs.

The focusing tension can be adjusted by turning the right focusing knob while holding the left one still.

Dioptric compensation

This compensation allows people wearing glasses to adjust the microscope so as to use it without them.

- Set the compensation on the eyepieces at zero.

- Set the zoom at the minimum magnification and focus the sample by observing it only with your right eye.
- Then observe it only with your left eye and – without touching the focusing knobs – adjust the left eyepiece dioptic compensation until obtaining a clear and sharp image with this eye as well.
- Set the zoom at the maximum magnification and then adjust the focusing knobs to get a sharp sample image. The focusing will remain constant with all magnifications (parafocal zoom).

Magnification and working distance

The total magnification depends on the eyepieces, on the zoom selected and on possible additional lenses used. By means of the lateral zooming knobs, the magnification can be modified by a factor of 6.7 with a perfect parafocality and sample centring.

Change the eyepieces and/or add an appropriate magnifying objective lens to increase the total magnification. The total magnification can be calculated with the following equation:

Total magnification = Eyepiece magnification x Zoom magnification x additional lens magnification

The accessories employed to modify the total magnification on the OPTIGEM-3 and OPTIGEM-4 are listed below:


ST-141	Eyepieces (pair) WF10x/22mm
ST-142	Eyepieces (pair) WF15x/16mm
ST-143	Eyepieces (pair) WF20x/12mm
ST-144	Eyepieces (pair) WF25x/9mm
ST-145	Micrometric eyepiece WF10x/22mm
SA01.5	Additional lens 1,5x
SA02	Additional lens 2x

The normal working distance for the standard configuration (without additional lenses) is 100 mm.

NOTE

Using an additional lens with a magnification higher than 1x means reducing the working distance. It may be necessary to lower the head focus holder, fixing it on the dedicated holes.

WARNING

- Always disconnect the power supply when the microscope is inactive for a long time, when cleaning it or when performing maintenance operations.
- Make sure that all electrical components are protected from humidity and liquids.

NOTE

- After using the microscope, always protect it with the provided dust cover.
- Keep the optical surfaces clean to ensure the best performance.

Protection from water

Do not use or store this instrument in humid environments. Humidity favours an oxidation process of the mechanical parts and the development of fungus and moulds on the optical surfaces. Humidity must not exceed 80%.

Cleaning

To clean the mechanical parts, use a dry non-abrasive cloth. In case of stains particularly difficult to remove, moisten the cloth with a little bit of ethanol.

NOTE

The lens cleaning has to be carried out very carefully so as not to damage the surface.

To remove dust, the use of compressed air is recommended (for example, spray cans such as those used to remove dust from computer keyboards). Alternatively a buff (such as the one normally provided with glasses) can be suited for this task.

NOTE

Avoid the use of solvents, which may jeopardize the system quality.

Moving the microscope

The microscope should not be moved often; when it is necessary ensure that all parts are properly protected. It would be better to disassemble the instrument and place it in its original packing.

Electric parts of the microscope

WARNING

Except when changing the lamps, do not operate on the electric parts. In case of problems, get in touch with your dealer.

Replacing the bright field lamp


WARNING

- Always disconnect the power supply cable from the microscope when not using it for a long time, when cleaning it and when performing any maintenance operations.
- Use only **6V 30W halogen lamps**
- Turn the stand upside-down so as to be able to access the bottom part of the illumination system.

- Loosen the lamp-holder box screw and remove the old lamp.


- Insert a new lamp. Beware not to touch with your hands the lamp glass surface.
- Close the lamp-holder box.

Replacing the lamp for incident light

- Always disconnect the power supply cable from the microscope when not using it for a long time, when cleaning it and when performing any maintenance operations.
- Use only a **fluorescent lamp PL 7W 6400 K**
- Turn the screw to open the fluorescence lamp box.
- Remove the old lamp by pulling it .

- Insert a new lamp by pushing it deep into the connector.

Replacing the fuses

- Always disconnect the power supply cable from the microscope when not using it for a long time, when cleaning it and when performing any maintenance operations.
- Use only **T3.15A 250V 20mm x 5mm fuses.**

- By means of a flat head screwdriver, loosen the fuse holder and remove it.

8.0 MAINTENANCE AND CLEANING

- Remove the fuse from the fuse holder and replace it with a new one.

- Put the fuse holder back in the microscope base. If necessary, replace the other fuse as well (there are two of them in the base)

9.0 TROUBLESHOOTING

- *The illumination system doesn't turn on.*
 - a) Check if the power cable is connected and the main switch on the position 1.
 - b) Check both fuses located in the microscope base.
- *A dark field image of the sample can not be obtained*
 - Check if the black plate has been inserted by means of the small lever located on the object plate near the illumination hole.
 - Get the sample close to the illumination hole by moving the clamp and use the iris diaphragm to find the ideal darkfield light angle.
 - Increase the transmitted light intensity by means of the adjustment knob located in the lower part of the microscope base
- *The image is not sharp for all zoom values.*
 - a) Carry out the dioptic compensation adjustment as described in the chapter dedicated to the use of the microscope
 - b) Focus the sample with the zoom set at the maximum magnification, then set it at the minimum magnification and check that the parfocality is kept constant.
- *The camera doesn't show any images*
 - a) Put the photo/video output selector on the optical head on the IN position.


The accessories' and spare parts' codes for OPTIGEM-3 and OPTIGEM-4 are reported in the following chart:

ST-141	Eyepieces (pair) WF10x/22mm
ST-142	Eyepieces (pair) WF15x/16mm
ST-143	Eyepieces (pair) WF20x/12mm
ST-144	Eyepieces (pair) WF25x/9mm
ST-145	Micrometric eyepiece WF10x/22mm
SA01.5	Additional lens 1,5x
SA02	Additional lens 2x
ST-231	Polarisation analysis set
ST-207	Diamonds proportion analysis set
ST-147	Camera adapter (1/3" sensor)
ST-147.1	Camera adapter (1/2" sensor)
ST-146	REFLEX camera adapter


11.0 DIMENSIONS

[measures are shown in mm]

Art.13 Dlsg 25 July 2005 N°151. “According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal.”

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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