B-1000 Series

INSTRUCTION MANUAL

Model
B-1000FL HBO

Ver. 4.0 2019

Table of Contents

1. Warning	4
2. Symbols and conventions	4
3. Safety Information	4
4. Intended use	4
5. Overview	5
5.1 Manual version	5
5.2 Motorised version	7
6. Unpacking	8
7. Assembling	8
7.1 Manual version	8
7.2 Motorised version	9
7.3 Assembling the microscope	10
7.4 Only for motorised version	14
8. Brightfield observation procedures (transmitted light)	15
9. Use of the microscope in Brightfield (transmitted light)	16
9.1 General switch on	16
9.2 Control keyboard	16
9.3 Brightness adjustment	16
9.4 Adjust the observation head	17
9.5 Adjust the interpupillary distance	17
9.6 Diopter adjustment	17
9.7 Use of eyeshields	18
9.8 Light path selection	18
9.9 Coarse focus tension adjustment	18
9.10 Focus stop lever	19
9.11 Stage	19
9.12 Condenser centering	20
9.13 Effect of field diaphragm	20
9.14 Aperture diaphragm	20
9.14 Use of oil immersion objective	21
9.15 Only for motorised version	21
9.15.1 Nosepiece rotation	21
9.15.2 Focusing	22
9.15.3 Stage	22
10. Fluorescence observation procedures	23
11. Use of the microscope in Fluorescence	24
11.1 HBO bulb switch on	24
11.2 Centering of HBO bulb	24
11.3 Use of diaphragms	26
11.3.1 Field diaphragm	26
11.3.2 Aperture diaphragm	26
11.4 Use of the shutter	26
11.5 Use of fluorescence	27
11.6 Use of light excluding plate	27
12. Use of condenser for Brightfield / Darkfield / Phase Contrast (optional)	28
12.1 Brightfield observation (BF)	28
12.2 Darkfield observation (DF)	28
12.3 Phase contrast observation (PH)	29
12.4 Use of green filter	30
13. Simultaneous observation in Phase Contrast + Fluorescence	30
14. Differential Interference Contrast DIC Observation (optional)	31
14.1 Koehler DIC transmitted light	31
14.2 Nomarski DIC transmitted light	32
15. Simultaneous Observation in Fluorescence and DIC	34
15.1 Koehler DIC reflected light	34
15.2 Nomarski DIC reflected light	36
16. Microphotography	39
16.1 Use of "C"mount camera	39
16.2 Use of Reflex camera	39

17. Maintenance	40
18. Troubleshooting	41
Equipment disposal	43

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

2. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

3. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

4. Intended use

Standard models


For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

IVD Models


Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

5. Overview

5.1 Manual version


Opposite side

5.2 Motorised version

Only the parts related to the motors are showed; all the other components of the microscope remain unchanged with respect to the manual version.

Opposite side

6. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

7. Assembling

Once opened the box, the microscope parts are the following:

7.1 Manual version

- | | |
|--------------------------------|--|
| ① Frame | ⑩ Microscope power supply |
| ② Fluorescence epi-illuminator | ⑪ Fluorescence power supply |
| ③ Objectives | ⑫ Power cord for fluorescence power supply |
| ④ Condenser | ⑬ Immersion oil (if 100x is included in the configuration) |
| ⑤ Observation head | ⑭ Allen wrench |
| ⑥ Eyepieces | ⑮ Dust cover |
| ⑦ Stage | |
| ⑧ HBO lamp housing | |
| ⑨ Darkening plate | |

7.2 Motorised version

- ① Frame
- ② Fluorescence epi-illuminator
- ③ Objectives
- ④ Condenser
- ⑤ Observation head
- ⑥ Eyepieces
- ⑦ Stage
- ⑧ HBO lamp housing
- ⑨ Darkening plate
- ⑩ Microscope power supply

- ⑪ Fluorescence power supply
- ⑫ Motor power supply
- ⑬ Serial cable
- ⑭ PS/2 Mouse
- ⑮ Power cord for fluorescence power supply
- ⑯ Immersion oil (if 100x is included in the configuration)
- ⑰ Allen wrench
- ⑱ Dust cover

7.3 Assembling the microscope

1. Insert the fluorescence epi-illuminator over the microscope and tighten the screw with the 2 mm Allen wrench. (Fig. 1)

2. Screw the extension tube on the rear of the attachment, using the 3 provided allen screws. (Fig. 2)

3. Screw the lamp housing on the extension tube, tightening the screws inside the holes. (Fig. 3)

4. Insert the optical head above the attachment, using the 2 mm Allen wrench to tighten the screw. (Fig. 4)


5. Insert eyepieces into the empty eyepiece sleeves. (Fig. 5)

6. Insert the condenser under the stage.
 - Check that it is correctly inserted in its housing (under the condenser there is a plug that must completely fit into the guide of the condenser holder). (Fig. 6)
7. Lock the condenser fixing knob ①.

8. Mount the stage: lower the support using the coarse focus knob, then place the stage and firmly tighten the lock screw ②. (Fig. 7)

9. Screw the objectives into the revolving nosepiece, clockwise, in order of magnification. (Fig. 8)


10. Insert the power supply jack in the socket placed in the rear side of the frame. (Fig. 9)

11. Open the lamp housing using the door lock screw ① and remove the lamp holder. (Fig. 10)

12. Remove the plastic block ② from the lamp holder (or the exhausted one in case of replacement) by loosening the two locking screws ③. (Fig. 11)

13. Insert the mercury bulb ④ (respect the polarity of the bulb), tighten the locking screws and refit the lamp holder inside the lamp housing. (Fig. 12)

14. Connect the cable from the lamp housing to the external fluorescent power supply and then lower the metal fixing tab ①. (Fig. 13)

15. Plug the power cord into the socket ②. (Fig. 14)


**Before connecting the power cord, secure the lamp housing cable to the power supply.
If the power cord is connected first, there may be a risk of electrical shock.**


- Disconnect all electrical cables before installing or replacing the bulb.
- The bulb has an anode and a cathode of different sizes. Respect the polarity during assembly, respecting the bulb dimensions.
- Do not touch the bulb of the lamp with bare hands to leave no traces of grease on the bulb. If this happens, clean the bulb with a soft cloth before turning on the lamp.
- The lamp has an average life of about 200-250 hours: a time counter and a voltage indicator are shown on the lamp power supply. Replace the lamp when the hour count exceeds 250 or if the voltage drops below 4.5A.
- During use, the lamp, the lamp housing and the surrounding environment become hot.
- Before replacing the lamp, switch off the power supply, disconnect all cables and wait for the lamp and the lamp housing to cool.
- After switching on the lamp, wait at least 10-15 minutes before switching it off.
- After switching off the lamp, wait for 5-10 minutes before switching it on again so that the mercury vapors have time to condense.
- The lamp contains ultraviolet radiation that could be harmful to eyes and skin. Always look at the lamp arc through the provided orange screen.

7.4 Only for motorised version

16. Assemble the stage in the same way as the manual version. Check the perfect alignment of the rear part of the stage with the rear arm of the frame. An imperfect alignment could lead to an incorrect functioning of the system. (Fig. 15)

17. Connect the cable ① from the stage to the frame and tight the locking screws of the connectors ②. (Fig. 16)


18. Connect the provided cables: ③ 12V power supply for the motorised parts; ④ 6V microscope power supply; ⑤ serial cable; ⑥ PS/2 mouse. (Fig. 17)

- **Connect power cables as the last step.**


8. Brightfield observation procedures (transmitted light)

9. Use of the microscope in Brightfield (transmitted light)

9.1 General switch on


To activate the transmitted light illuminator turn the main switch ①, located on the left side of the stand, to the position “1”. (Fig. 18)

9.2 Control keyboard

B-1000 illumination can be managed through the keyboard placed on the left of the stand. (Fig. 19)

- **ON-OFF (②)**: press this key (after switching the main switch on “1”) to turn ON or OFF the microscope LED.
- **BOOST (③)**: press this button in order to increase the brightness (useful for high-magnification objectives or very opaque specimens). **Do not enable boost mode while observing with low magnification objectives (4x, 10x) with fully open diaphragm: the high brightness may hurt user's eyes.**
- **AUTO OFF (④)**: if you want the illuminator to switch off automatically, press this button until 15, 30 or 60 minutes delay is set. After this period of time, the light will turn off. You have to press the ON-OFF button to turn it on again.

9.3 Brightness adjustment

Use the brightness adjustment dial ⑤ on the left side of the microscope to increase or decrease the light intensity on the specimen. (Fig. 20)

9.4 Adjust the observation head

Loosen the lock-screw ①, turn the observation head to a comfortable position for observation, and then lock the lock-screw. (Fig. 21)


9.5 Adjust the interpupillary distance

Observing with both eyes, hold the two eyepiece prism assemblies. Rotate them around their common axis until the fields of view coincide.

- **The graduation on the interpupillary distance indicator ②, pointed by the spot “.” on the eyepiece holder, shows the distance between the operator’s eyes. (Fig. 22)**

The range of the interpupillary distance is 48-75 mm.

9.6 Diopter adjustment

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
 2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the dioptic adjustment ring ③ to compensate. (Fig. 23)
- **The adjustment range is ± 5 diopter. The number shown on the adjustment ring graduation should correspond to the operator’s dioptic correction.**

9.7 Use of eyeshields

• Use with eyeglasses

Fold rubber eyeshields with both hands. Folded eyeshields avoid scratching the lenses of eyeglasses. (Fig. 24)

• Use without eyeglasses

Raise eyeshields and observe at the microscope placing eyes to the shields, avoiding external light to disturb the observation. (Fig. 25)

9.8 Light path selection

- The observation head is equipped with an optical path selector that allows the light to be distributed to the eyepieces and to the photo / TV port.
- Move the selector ① to one of the three possible positions to distribute the light. (Fig. 26)

POSITION	LIGHT
IN	100% EYEPIECES
MIDDLE	50% EYEPIECES / 50% TV
OUT	100% TV

9.9 Coarse focus tension adjustment

The tension of the coarse focusing knob is factory preset.


- To modify the tension according to personal's needs, rotate the ring ②. (Fig. 27)
- Clockwise rotation increases the tension.
 - If the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.

9.10 Focus stop lever

The upper limit knob has two functions: prevent the contact between slide and objective and acts as "focus memory".

1. After focussing the specimen, pull the lever ① toward the front of the microscope and lock it. (Fig. 28).
- In this way the focus upper limit is set.
2. Now one can lower the stage with coarse focus knob, replace the specimen and raise again the stage up to the upper limit: specimen will be in approximate focus and will need a fine adjustment to get the proper focus..
- **Fine focus movement is not affected by the coarse focus lock.**
 - **To unlock, move the lever in the opposite direction to the one used for the locking.**
 - **Two stoppers ② are inserted on the frame. DO NOT REMOVE THE TWO STOPPERS.**

9.11 Stage

Stage accepts standard slides 26 x 76 mm, thickness 1,2 mm with coverslide 0,17mm. (Fig. 29)

It is possible to place two slides side by side on the stage.

1. Open the spring arm of the slide holder ① and place from the front the slide on the stage.
 2. Gently release the spring arm of the slide holder.
- **A sudden release of the spring arm could cause the falling of the slide.**

9.12 Condenser centering

1. Place the specimen on the stage, insert 10x objective into the light path and focus.
2. Insert the front lens of the swing-out condenser using the lever ①. (Fig. 30)
3. Rotate the field diaphragm ring ② in the direction showed by the arrow, to fully close the diaphragm.
4. Rotate the condenser height adjustment knob ③ to focus the edges of the diaphragm.
5. Rotate the two centering screws ④ to bring the bright spot in the center of the field of view.
6. Gradually open the diaphragm. The condenser is centered when the diaphragm image is symmetrical to the field of view.
7. In normal use, open the diaphragm until it circumscribes the field of view.

Fig. 30

9.13 Effect of field diaphragm

Field diaphragm adjusts the illuminated area to obtain a high contrast image.

Set the diaphragm according to the objective in use until it circumscribes the field of view, in order to eliminate unnecessary light to eyepieces. (Fig. 31)


Fig. 31

9.14 Aperture diaphragm

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image
- With low contrast specimens set the numerical aperture value ⑤ (printed on the condenser ring) to about 70%-80% of the objective's N.A. (Fig. 32). If necessary, remove one eyepiece and, looking into empty sleeve, adjust the condenser's ring in order to obtain an image like the one in Fig. 33.

Example: with objective PLAN 40x / 0.65 set the scale to $0.65 \times 0.8 = 0,52$

Fig. 32


Fig. 33

9.14 Use of oil immersion objective

1. Focus the specimen with a low power objective.
2. Lower the stage (remembering to lock the coarse upper limit knob).
3. Put a drop of oil (provided) on the area of the specimen to be observed. (Fig. 34)
- **Make sure that there are no oil bubbles. Air bubbles in the oil damage the image quality.**
- To check for bubbles: remove an eyepiece, fully open the aperture diaphragm and observe the objective exit pupil. (The pupil must be round and bright).
- To remove the bubbles, gently move the nosepiece to the right and left to move the immersion objective a few times and allow the air bubbles to move.
4. Insert immersion objective.
5. Return the stage to the upper focusing point and obtain an optimal focus using the fine focus knob.
6. After use, gently remove the oil with a soft paper towel or a lightly moistened optic paper with a mixture of ethyl ether (70%) and absolute alcohol (30%).
- **The immersion oil, if not immediately cleaned, could crystallize creating a glass-like layer. In this situation the observation of the specimen would be difficult if not impossible due to the presence of an additional thickness on the objective.**

9.15 Only for motorised version

9.15.1 Nosepiece rotation

1. To change magnification it is possible operate on the nosepiece movement buttons located on the right side of the frame. (Fig. 35). Orange button ① rotates the nosepiece clockwise, while blue button ② rotates the nosepiece counterclockwise.
2. As an alternative it is possible operate on the right and left mouse buttons.


9.15.2 Focusing

Focus motor is activated through the mouse wheel. Front or rear rotation of the mouse wheel raises or lowers the stage. (Fig. 36)


9.15.3 Stage

1. Stage is moved through the mouse movement. A mouse movement to the front or to the back ③ causes a stage movement of the stage along the Y axis, while a left or right movement ④ causes a stage movement of the stage along the X axis. (Fig. 37)
2. It is always possible, however, operate on the translation knobs of the stage for a manual movement.


10. Fluorescence observation procedures

11. Use of the microscope in Fluorescence

11.1 HBO bulb switch on

1. Turn on the power supply using the main switch. (Fig. 38)
2. The “HBO ON” LED on the front panel of the power supply unit will only light up if the lamp is turned on correctly.
 - Wait until a value close to 4.5 A appears on the current display. If the current drops below 4 A, replace the bulb.
 - You should now wait for at least 5 minutes before aligning and using the bulb.

3. Move the lever to the “O” position to open the shutter. (Fig. 39)


4. Select the filterset by rotating the wheel to the desired position. (Fig. 40)


11.2 Centering of HBO bulb



1. Rotate the revolving nosepiece into an empty position (no objectives) and remove the dust cap, or remove one objective from the nosepiece.
2. Insert into the light path the filter cube “B” (fig. 40) and put a piece of white paper on the stage. (Fig. 42)


3. Operating on the focusing screw of the collector lens ① and on the centering screws ② (Fig. 41) try to obtain the bright spot of the bulb arc. (Fig. 42).

4. Using the focusing screw of the collector lens ① focus the arc's image on the paper. The bright spot must be as sharper as possible. (Fig. 43)
5. Using the centering screws ② on the side of the lamp housing, center the image of the arc. (Fig. 43-44)

6. Using the focusing screw of the collector lens ① enlarge the image until a homogeneous illumination is achieved. (Fig. 45). At this point, insert an objective into the optical path and, looking into the eyepieces, optimize the illumination always using the screws ① and ②..

- After replacing the exhausted bulb, reset the time counter on the power supply by pressing the “Reset” button ③. (Fig. 46)

11.3 Use of diaphragms

The illuminator is provided with centerable field and aperture diaphragms. (Fig. 47)

Centering and use procedure of reflected light diaphragm is the same to the one described for transmitted light diaphragms in chapter 9.11 “Condenser centering”.


11.3.1 Field diaphragm

Operate on the “FS” ① to close the field diaphragm and, using the provided Allen wrench, operate on the centering screws ②. (Fig. 48)

11.3.2 Aperture diaphragm

Operate on the “FS” ③ to close the aperture diaphragm and, using the provided Allen wrench, operate on the centering screws ④. (Fig. 48)

11.4 Use of the shutter

- Open the shutter by moving the lever to the “○” position to work in fluorescence. (Fig. 49)
 - Close the shutter by moving the lever to the “●” position, as the observation should be interrupted for a limited time and the sample should not be subjected to unnecessary illumination during the observation period.
- Frequent switching on and off of the HBO bulb significantly reduces its service life.

11.5 Use of fluorescence

The filter turret is provided with 6 positions.

- B and G filters are installed in the factory, while two other filters (UV and V) are optional.
- It is possible, however, use empty filter cubes and install custom filtersets.

FILTER NAME	EXCITATION FILTER	DICHROIC MIRROR	BARRIER FILTER	APPLICATIONS
B	460-490 nm	500 nm	520LP nm	• FITC: fluorescent antibodies • Achridine orange: DNA, RNA • Auramine
G	510-550 nm	570 nm	590LP nm	• Rhodamine, TRITC: fluorescent antibodies • Propidium Iodide: DNA, RNA • RFP
UV	330-385 nm	400 nm	420LP nm	• AT-selective • Chromosome and nuclear counterstaining • Chromosome banding
V	400-410 nm	455 nm	455LP nm	• Amine reaction

11.6 Use of light excluding plate

- **Microscope is provided with a light excluding plate that can be placed on the stage and prevents flare and reflections coming from the condenser front lens.**

The plate can be used in two different ways.

1. Mode n° 1: place the plate on the stage (under the slide holder) and place the slide directly over the plate. (Fig. 50)
 2. Mode n° 2: lower the condenser and insert the plate between the two layers of the stage. (Fig. 51).
- **In both cases it is possible to move the sample using the stage X-Y translation knobs.**

Fig. 50

Fig. 51

12. Use of condenser for Brightfield / Darkfield / Phase Contrast (optional)

Universal condenser provided with B-1000 allows observation in brightfield, darkfield and phase contrast.

Observation mode	Condenser turret position
Brightfield	BF (Fig. 52)
Darkfield	DF (Fig. 53)
Phase contrast 10x	10/20 (Fig. 54)
Phase contrast 20x	10/20 (Fig. 54)
Phase contrast 40x	40 (Fig. 55)
Phase contrast 100x	100 (Fig. 56)

12.1 Brightfield observation (BF)





1. Rotate the condenser turret to insert the “BF” position.
2. Now repeat the steps described in the procedure “*Brightfield observation procedures (transmitted light)*”.

12.2 Darkfield observation (DF)

1. Rotate the condenser turret to insert the “DF” position.
 2. Open the aperture diaphragm.
 3. Place a specimen on the stage and focus.
 4. Observing into eyepieces raise or lower the condenser until a homogeneous illumination of the specimen can be achieved, thus obtaining a proper darkfield effect.
- **Darkfield requires a huge amount of light. Switching from darkfield to brightfield, one could be dazzled. Do not keep your eyes on the eyepieces when moving the condenser turret from DF to BF.**
 - **“Dry” darkfield observation, that is, without the use of oil, is only possible with objectives with N.A. lower than 0.7.**
 - **When observing in darkfield, it may be necessary to raise the condenser from the normal position to obtain a more homogeneous illumination. This is not a defect.**

12.3 Phase contrast observation (PH)

1. Center the condenser as already described in chapter 9.11.
 2. Rotate the condenser turret to insert the “10/20” position
 3. Insert 10x objective into the light path.
 4. Place a specimen on the stage and focus.
 5. Remove one eyepiece and insert the centering telescope. (Fig. 57)
 6. Rotate the upper part of the centering telescope until the two phase rings (one dark and one bright) visible in the telescope are in focus. (Fig. 57)
 7. Using centering screws on the condenser ① (Fig. 58), center the phase rings to make the bright ring ② be concentric to the dark ring ③. (Fig. 59-60)
 8. Insert 20x objective (do not rotate the condenser turret) and check the centering of the two rings. (Fig. 60)
 9. Repeat the same operation with other objectives to check the ring centering: 40x objective – turret position “40”, 100x objective – turret position “100”.
 10. At the end remove the centering telescope, reinstall the eyepiece and begin observation.
- **With 40x and 100x objectives it may be useful to slightly raise the condenser, to obtain a better projection of the phase rings. This is not a defect.**
 - **With the 4X objective, the condenser could have a dark halo at the periphery of the field of view. This is not to be considered a defect.**

12.4 Use of green filter

- Green filter is used to increase the contrast of the image during phase contrast observation.
- Place the filter on the field lens of the microscope and begin the observation. (Fig. 61)
- For brightfield or darkfield observation it advisable to remove the green filter from the light path.

13. Simultaneous observation in Phase Contrast + Fluorescence

- The microscope allows observation in transmitted light Phase contrast in combination with reflected light Fluorescence.
 - Samples with rapid decay must first be observed in Fluorescence and then in Phase Contrast.
 - The combined observation allows you to easily identify some areas of the sample that emit fluorescence.
1. Turn on the power supply for the HBO fluorescent bulb and wait 5 minutes before the arc stabilizes.
 2. Move the filter selector to an empty position or, if the filter holder is complete, to the position containing the UV filter.
 3. Insert the desired PH lens and rotate the phase contrast condenser turret to the position containing the corresponding phase ring.
 4. Focus the sample.
 5. Adjust the light intensity of the transmitted light.
 6. Move the fluorescence filter selector to the desired position.
 7. To obtain the proper observation of the sample, adjust the light intensity of the transmitted light to modulate the intensity of the fluorescence with respect of the phase contrast.

14. Differential Interference Contrast DIC Observation (optional)

- The microscope allows the observation in Differential Interferential Contrast (DIC) with two different methods: Koehler DIC and Nomarski DIC.
- Samples with rapid decay must first be observed in Fluorescence and then in DIC.
- The Koehler DIC method is the simplest both from the point of view of installation and from the point of view of use, while the Nomarski DIC method provides for a more complex setup.
- Both methods work in transmitted light but can be used in combination with fluorescence observation reflected light and therefore the transmitted light only set-up is different from the combined fluorescence observation set-up.


14.1 Koehler DIC transmitted light

The observation in Koehler DIC in transmitted light requires the kit consisting of the following accessories: Polarizer ①, Analyzer for transmitted light ②, Interferential green filter IF550 ③ and DIC slider ④. (Fig. 62)


1. Place the polarizer on the lens at the base of the microscope.

2. Remove the dummy slider from the nosepiece and insert the analyzer into the hole of the dummy slider, then insert the assembly ⑤ into the slot ⑥. (Fig. 63)
3. Remove the slide from the stage.
4. Rotate the polarizer at the base of the microscope to achieve maximum darkening of the eyepieces.

5. Once the maximum darkening is achieved, remove the analyzer from the dummy slider and insert it into the DIC prism. Now insert the DIC slider ⑦ into the slot ⑥. (Fig. 64)
6. Close the condenser aperture diaphragm a little.


7. Put a specimen on the stage and focus.
 8. Begin the observation by rotating the DIC ⑧ slider knob to obtain a three-dimensional sample effect. (Fig. 65)
- For a better effect on the image you can use the green filter IF550 which must be placed on top of the polarizer.

14.2 Nomarski DIC transmitted light


The observation in Nomarski DIC in transmitted light requires the kit consisting of the following accessories: Universal condenser ① (containing the dedicated DIC prisms according to the objectives in use), Analyzer for transmitted light ②, DIC slider ③. (Fig. 66)

Universal Condenser Controls


- Ⓐ Optical inserts markers
- Ⓑ Aperture diaphragm scale
- Ⓒ Aperture diaphragm lever
- Ⓓ Top lens
- Ⓔ Top lens lever
- Ⓕ Optical inserts centering screws


1. Using the knob ①, insert the polarizer ② embedded in the condenser and loosen the polarizer rotation fixing screw ③. (Fig. 67)

- Ⓖ Polarizer rotation fixing screw
- Ⓗ Polarizer rotation knob
- Ⓘ Polarizer in/out knob
- ⓵ Polarizer slider locking screw
- ⓶ Polarizer
- ⓷ Indicator markers

2. Remove the dummy slider from the nosepiece and insert the analyzer into the dummy slider, then insert the ④ assembly into the slot ⑤. (Fig. 68)

3. Remove the slide from the stage.
4. Turn the polarizer knob ⑧ under the condenser to achieve maximum darkening of the eyepieces, and then tighten the polarizer locking screw ⑥. (Fig. 69)

Fig. 69

5. Once the maximum darkening is found, remove the slider from the nosepiece, remove the analyzer from the dummy slider and insert it into the DIC prism. Now insert the DIC slider ⑥ into the slot ⑤. (Fig. 70)


Fig. 70

6. Rotate the condenser turret ⑦ to insert the DIC prism matching the objective in use. (Fig. 71)
- **The condenser is supplied with magnetic markers. Each marker is specific to the type of insert mounted in the condenser (DIC, PH, DF, etc.).**

Fig. 71

7. Put a specimen on the stage and focus.
8. Begin the observation by turning the knob on the DIC slider ⑧ to obtain a three-dimensional effect of the sample. (Fig. 72)

Fig. 72

15. Simultaneous Observation in Fluorescence and DIC

- The microscope allows observation in transmitted light Differential Interference Contrast DIC in combination with reflected light Fluorescence.
- Samples with rapid decay must first be observed in Fluorescence and then in DIC.
- The combined observation allows you to easily identify some areas of the sample that emit fluorescence.

15.1 Koehler DIC reflected light

Koehler DIC observation combined with fluorescence needs the use of the following accessories: Polarizer ①, Reflected light analyzer ②, Interferential green filter IF550 ③ and DIC slider ④. (Fig. 73)

1. Place the polarizer on the lens at the base of the microscope.

Fig. 73

2. Insert the analyzer into the slot ⑤ placed on the right side of the reflected light illuminator. (Fig. 74)
3. Move the filter holder selector to an empty position or, if the filter holder is complete, to the position containing the UV filter.

Fig. 74

4. Remove the slide from the stage.
5. Put on "0" the scale of the reflected light analyzer. (Fig. 75)

Fig. 75

6. Rotate the polarizer at the base of the microscope to achieve maximum darkening of the eyepieces.
7. Once the maximum darkening is achieved, remove the analyzer from the dummy slider and insert it into the DIC prism. Now insert the DIC slider ⑥ into the slot ⑦. (Fig. 76)
8. Close the condenser aperture diaphragm a little.


Fig. 76

9. Put a specimen on the stage and focus.
10. Begin the observation by rotating the DIC slider knob ⑧ to obtain a three-dimensional sample effect. (Fig. 77)


11. Insert the desired fluorescence filter and open the shutter ⑨. (Fig. 78)
12. Adjust the transmitted light intensity to optimize fluorescence and DIC observation until the best image contrast can be achieved.

15.2 Nomarski DIC reflected light

Nomarski DIC observation combined with fluorescence needs the use of the following accessories: Universal condenser ① (containing the dedicated DIC prisms according to the objectives in use), Analyzer for reflected light ②, DIC slider ③. (Fig. 79)


Universal Condenser Controls


- ① Universal condenser
- ② Analyzer for reflected light
- ③ DIC slider

1. Using the knob ①, insert the polarizer ② embedded in the condenser and loosen the polarizer rotation fixing screw ③. (Fig. 80)

2. Insert the analyzer into the slot ④ placed on the right side of the reflected light illuminator. (Fig. 81)
3. Move the filter holder selector to an empty position or, if the filter holder is complete, to the position containing the UV filter.

- ④ Polarizer rotation fixing screw
- ⑤ Polarizer rotation knob
- ⑥ Polarizer in/out knob
- ⑦ Polarizer slider locking screw
- ⑧ Polarizer
- ⑨ Indicator markers

4. Remove the slide from the stage.
5. Put on "0" the scale of the reflected light analyzer. (Fig. 82)

Fig. 82

6. Turn the polarizer knob ⑥ under the condenser to achieve maximum darkening of the eyepieces, and then tighten the polarizer locking screw ⑦. (Fig. 83)

Fig. 83

7. Once the maximum darkening is found insert the DIC slider ⑧ into the slot ⑨. (Fig. 84)


Fig. 84

8. Rotate the condenser turret ⑩ to insert the DIC prism matching the objective in use. (Fig. 85)
- **The condenser is supplied with magnetic markers. Each marker is specific to the type of insert mounted in the condenser (DIC, PH, DF, etc.).**

Fig. 85

9. Put a specimen on the stage and focus.
10. Begin the observation by rotating the DIC slider knob to ⑧ obtain a three-dimensional sample effect. (Fig. 86)

11. Insert the desired fluorescence filter and open the shutter ⑨. (Fig. 87)
12. Adjust the transmitted light intensity to optimize fluorescence and DIC observation until the best image contrast can be achieved.

16. Microphotography

16.1 Use of “C”mount camera

1. Loosen the clamping screw ① on the trinocular port and remove the dust cap ②. (Fig. 88)

2. Screw the C-mount adapter ③ to the camera ④ and insert the round dovetail of the C-mount into the empty hole of the trinocular port, then tighten the clamping screw ①. (Fig. 89)

16.2 Use of Reflex camera

1. Insert the Reflex adapter ① into the relay tube ②.
2. Screw the “T2” ring ③ (not provided) to the reflex adapter.
3. Connect the Reflex camera ④ to the “T2” ring just installed (Fig. 90).
4. Mount the other end of the relay tube ② into the empty hole of the trinocular port, then tighten the clamping screw. (Fig. 88)
 - “T2” ring is not provided along with the microscope, but is commercially available.
 - While shooting dark specimens, darken eyepieces and viewfinder with a dark cloth to minimize the diffused light.
 - To calculate the magnification of the camera: objective magnification * camera magnification * lens magnification.
 - **If using an SLR camera, mirror movement may cause the camera to vibrate.**
 - **We suggest lifting the mirror, using long exposure times and a remote cord.**

17. Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

18. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
LED operates, but field of view remains dark.	Power supply is unplugged.	Connect
	Brightness is too low	Set brightness to a proper level
	Fluorescence filter selector is not in a click stop	Move the selector to a click stop
	Fluorescence shutter is closed	Open the shutter lo shutter
	Fluorescence filter is not suitable for the specimen	Use a suitable filter
Field of view is obscured or not evenly illuminated	Revolving nosepiece is not correctly engaged.	Make sure that the revolving nosepiece clicks properly into place.
	The turret of the phase contrast condenser is in an incorrect position	Move the turret to a click stop
Dirt or dust is visible in the field of view.	Dirt/dust on the specimen	Clean the specimen
	Dirt/dust on the eyepieces	Clean the eyepieces
Image looks double	Aperture iris diaphragm is stopped down too far.	Open aperture iris diaphragm.
	The condenser is not well centered or it is in a wrong height	Set the condenser according to Koehler settings.
Visibility is poor. • Image is not sharp. • Contrast is poor. • Details are indistinct. • Image glares	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture iris diaphragm is too closed or too open.	Adjust aperture iris diaphragm.
	Dust or dirt on lenses (condenser, objectives, eyepieces and slide)	Clean thoroughly.
	For transmitted light observation, the coverglass thickness must not exceed 0.17mm	Use a coverglass with thickness 0.17mm
	For phase contrast observation, a brightfield objective is used instead a phase contrast one	Use a phase contrast objective
	Phase rings of objective and condenser are not well centered	Operate on centering screws
	Objective in use is not compatible with condenser phase ring	Use a compatible objective
	Focus is not even	Slide holder is not flat. Move the specimen to a flat position.
One side of the image is unfocused	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Slide is mounted not in a flat position (tilted)	Place the specimen in a flat position on the stage
	Poor quality of the glass slide	Use a glass slide with higher quality
II. Mechanical Section:		
Coarse focus knob is hard to turn	Tension adjustment ring is too tight	Loosen tension adjustment ring
Focus is unstable	Tension adjustment ring is too loose	Tighten tension adjustment ring
III. Electrical Section:		
LED doesn't turn on.	Power supply not connected	Check for proper connection
Brightness is not enough	Brightness setting is too low	Adjust brightness
Light blinks	Power supply not well connected	Check for proper connection

IV. Observation tube:		
The field of view of one eye does not match the field of view of the other eye	Interpupillary distance is incorrect.	Adjust interpupillary distance.
	Incorrect diopter adjustment.	Adjust diopter.
	Your view is not accustomed to microscope observation.	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope.
V. Microphotography:		
Image edge is unfocused	To a certain extent it is due to achromatic objectives features	To minimize the problem, set the aperture diaphragm in a proper position
Bright spots appear on the image	Stray light entering in the microscope through eyepieces or camera viewfinder	Cover eyepieces and viewfinder with a dark cloth

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.