B-190 Series

INSTRUCTION MANUAL

Model
B-190TBPL+

Ver. 1.4 2026

Table of contents

1.	Warning	3
2.	Safety Information	3
3.	Package content	4
4.	Unpacking	5
5.	Intended use	5
6.	Symbols and conventions	5
7.	Instrument description	6
8.	Assembling	7
8.1	Assembling the microscope	7
8.2	Polarizing set (optional)	9
9.	Use of the microscope	10
9.1	Switching on the microscope	10
9.2	Light intensity adjustment	10
9.3	Stage	10
9.4	Adjust the interpupillary distance	10
9.5	Diopter adjustment	11
9.6	Motorized focus system	11
9.6.1	Turning on the microscope	11
9.6.2	Use of focus	12
9.6.3	Stage retraction	12
9.7	Use of oil immersion objective	13
9.8	Aperture diaphragm	13
9.9	Use of the polarizer (optional)	13
10.	Use of software and digital head	14
11.	Micrometric Slide M-005	14
12.	Maintenance	15
13.	Troubleshooting	16
	Equipment disposal	17

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

2. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

3. Package content

- ① Frame
- ② Digital observation head
- ③ Eyepieces
- ④ Objectives (4X/10X/40X/100X)
- ⑤ Dust cover
- ⑥ Green filter

- ⑦ Immersion oil
- ⑧ Microscope power supply
- ⑨ Digital imaging unit power supply
- ⑩ USB cable 0,5 m
- ⑪ Digital imaging unit

NOTE: OPTIKA reserves the right to make corrections, modifications, enhancements, improvements and other changes to its products at any time without notice.

4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

5. Intended use

Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION


This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.


7. Instrument description

8. Assembling

8.1 Assembling the microscope

1. Remove the dust cap from the microscope frame and from the bottom of the observation head.
2. Insert the optical head above the stand and tighten the screw. (Fig. 1)
 - **Hold the head with one hand during the locking in order to avoid that the head falls.**
3. Insert both eyepieces into the tubes of the optical head. (Fig. 2)
4. Insert the power supply jack in the socket placed at the rear side of the microscope. (Fig. 3)
5. Fix the rotating part of the junction using the black wing-nut ①. (Fig. 4)




6. Slide the digital imaging unit onto the head holder from top to bottom and slide it all the way down its slot. (Fig. 5)
- To unlock the digital imaging unit proceed with the opposite operation: push toward up and remove it from the holder.

7. Plug one side of the USB cable ② to the digital head and the other side to the digital imaging unit using the USB-A connector ③. (Fig. 6-7)
8. Plug the power supply cable to the digital imaging unit for battery recharge using the USB-C connector ④. (Fig. 8)

- **NOTE: The tablet comes with a 12V/2A power adapter featuring a European-style plug and a USB-C plug. In regions where a different type of wall outlet is used, the user must use a USB-A/USB-C cable and a Power Delivery adapter with the following maximum specifications:**
Input 100/240V – 50/60Hz – 1.7A
Output 5/20V – 6.2/3.25A (Max 67W)

- The digital imaging unit was set to work only in landscape to facilitate operations with the integrated camera management software.
- To switch between landscape and portrait mode, try turning on the rotation lock option (to disable automatic screen rotation) or by selecting the relevant sub-option from the display orientation menu. Go to Settings > System > Display, then click the arrow down under display orientation to select landscape or portrait mode.


8.2 Polarizing set (optional)


1. Place the polarizer on the light exit ① at the base of the microscope. (Fig. 9)

2. Loosen the head fixing knob ② and remove the head from the microscope frame. (Fig. 10)

3. Insert the analyzer into the hole inside the frame ③. (Fig. 11)
4. Put back the head into its original position and lock the fixing knob.

9. Use of the microscope

9.1 Switching on the microscope

Operate on the main switch ① placed in the rear side of the microscope, moving the selector on "I". (Fig. 12)


Fig. 12

9.2 Light intensity adjustment

Operate on the light intensity dial to increase or decrease the illumination intensity. (Fig. 13)

Fig. 13

9.3 Stage

Stage accepts standard slides 26 x 76 mm, thickness 1.2 mm with coverslide 0.17 mm. (Fig. 14)

1. Open the spring arm of the slide holder ② and place the slide from the front on the stage.
 2. Gently release the spring arm of the slide holder.
- **A sudden release of the spring arm could cause the falling of the slide.**


Fig. 14

9.4 Adjust the interpupillary distance

Hold the right and left parts of the observation head using both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen. (Fig. 15)

- The graduation on the interpupillary distance indicator ③, pointed by the spot "." on the eyepiece holder, shows the distance between the operator's eyes.

The range of the interpupillary distance is 48- 75 mm.

Fig. 15


9.5 Diopter adjustment

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
 2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the diopter adjustment ring ① to compensate. (Fig. 16)
- The adjustment range is ± 5 diopter. The number indicated on the adjustment ring graduation should correspond to the operator's diopter correction.

9.6 Motorized focus system

- **NOTE: The microscope is equipped with a motorized focusing system. Focusing can be achieved ONLY by the left knob.**
- **The right knob is fixed and has no effect.**
- A label is attached to the back of the microscope explaining the functionality of the system. (Fig. 17)
- On the plastic front panel are 3 LEDs ① (Fig. 18). These are the indicators of the different focus modes (Macro "M", Semi-Micro "sm" and Micro "m").
- Each time the microscope is turned on (Chapter 9.1), the system performs an "auto-calibration" of the focus: this is for the system to detect the end limits.

9.6.1 Turning on the microscope

1. When the microscope is turned on, the stage is automatically lowered to the lower limit of its travel and then returned to where the stage was at the time the microscope was turned off the previous time.
- **During this operation the focus knob is inactive. Any rotation of the knob will have NO effect.**
2. When the stage is returned to the initial starting point, the three LEDs flash twice.
3. When the self-calibration operation is completed, all three LEDs ② (Macro mode) will remain lit. (Fig. 19)

9.6.2 Use of focus

- The motorized focusing system involves the use of a single knob ③. (Fig. 20)
- When the focus knob is turned, individual mechanical clicks are felt that correspond to individual motor steps.
- This is an intended effect to allow the operator to “figure out” by how many steps the motor moves during the focusing operation.
- **This cannot be considered a defect.**
- **NOTE: The system has an electronic upper limit switch. When the upper limit point allowed by the limit switch is reached, the system will not allow further ascent of the stage.**

1. Pressing the focus knob briefly (< 1 sec) illuminates the LEDs on the control panel (Fig. 21-22) related to the current system setting.
 - LED “m”: “micro” mode (only left LED lit)
 - LED “sm”: “semi-micro” mode (left LED + middle LED lit)
 - LED “M”: “macro” mode (three LEDs lit)
 - **NOTE: Pressing the button cycles through the three modes: macro --> semi-micro --> micro --> macro.**
 - **NOTE: Each time you turn on the microscope, the system resets to the “Macro” mode.**


Mode	Nominal resolution (single step) (μm)
Macro	250
Semi-micro	8
Micro	1

9.6.3 Stage retraction

The microscope includes a stage retraction function. This function is used to:

- Safely replace the slide.
 - Put oil on the slide without having to rotate the revolver.
1. Press the focus knob for about 2 seconds.
 2. After about 2 seconds, the three LEDs start flashing continuously, waiting for the knob to be released. After release, the stage drops to the lower limit point (about 17 mm).
 - **NOTE: Activating the retract function disables the focus knob. The knob will be reactivated only when the knob is pressed again and the stage returns to its original position. This is an intended effect and should not be considered a defect.**
 3. Once the lower limit is reached, the three LEDs begin flashing, waiting for the slide to be replaced or the oil to be placed in anticipation of using the 100x objective (see Chapter 9.7).
 4. When the desired operation has been performed, press the focus knob again: hold the knob down until the three LEDs stop flashing, then release it. The stage will return to its initial position, the LEDs will briefly flash twice, and then the LEDs that were the indicators of the focus mode in use before the retraction was activated will remain lit.
 5. Now normal operation can be resumed.

9.7 Use of oil immersion objective

1. Focus the specimen with a low power objective.
 2. Lower the stage.
 3. Put a drop of oil (provided) on the area of the specimen to be observed. (Fig. 23)
- **Make sure that there are no oil bubbles. Air bubbles in the oil damage the image quality.**
 - To check for bubbles: remove an eyepiece, fully open the aperture diaphragm and observe the objective exit pupil. (The pupil must be circular and bright).
 - To remove the bubbles, gently move the nosepiece to the right and left to move the immersion objective a few times and allow the air bubbles to move.
4. Insert immersion objective.
 5. Return the stage to the upper focusing point and obtain an optimal focus using the fine focus knob.
 6. After use, gently remove the oil with a soft paper towel or a lightly moistened optic paper with a mixture of ethyl ether (70%) and absolute ethyl alcohol (30%).
- **The immersion oil, if not immediately cleaned, could crystallize creating a glass-like layer. In this situation the observation of the specimen would be difficult (even not impossible) due to the presence of an additional thickness on the objective.**


9.8 Aperture diaphragm

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image. Move the diaphragm ring ① (Fig. 24) on the value corresponding to the objective in use. In this case the optimal setting of the condenser is achieved.
- With low contrast specimens set the numerical aperture to about 70%-80% of the objective's N.A. If necessary, remove on eyepiece and, looking into empty sleeve, adjust the condenser's diaphragm in order to obtain an image like the one in Fig. 25.
- **NOTE: The condenser has an A.N. of 1.25. However, this value is only obtained when the condenser is used "with oil," i.e., by placing a drop of oil on the front lens of the condenser. If the condenser is used "dry," the maximum A.N. that can be obtained is 0.9.**

9.9 Use of the polarizer (optional)

1. Remove the specimen from the stage.
2. Looking inside the eyepieces, rotate the polarizer until the darkest position is achieved.
3. Once the dark is achieved ("extinction" or "Crossed Nicol" position) it is possible to begin the observation.

10. Use of software and digital head

The camera inside the digital head is driven by PROVIEW software.

For the instructions about the use of the software, please refer to the specific instruction manual.

The manual can be accessed by selecting *Help > Help Contents*. A PDF file containing the software instruction manual will open.


11. Micrometric Slide M-005

Micrometric slide, 26x76mm, with 2 scales
(1mm/100div. for biological microscopes / 10mm/100div. for stereomicroscopes)

1 DIV=0.01mm

For biological microscopes calibration

1 DIV=0.1mm

For stereo microscopes calibration

12. Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region.
- The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

13. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
LED operates, but field of view remains dark	Power supply is unplugged	Connect
	Brightness is too low	Set brightness to a proper level
Dirt or dust is visible in the field of view	Dirt/dust on the specimen	Clean the specimen
	Dirt/dust on the eyepieces	Clean the eyepieces
Image looks double	Aperture diaphragm is stopped down too far	Open aperture diaphragm
Visibility is poor <ul style="list-style-type: none"> • Image is not good • Contrast is poor • Details are indistinct • Image glares 	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture diaphragm is too closed or too open	Adjust aperture diaphragm
	Dust or dirt on lenses (condenser, objectives, eyepieces and slide)	Clean thoroughly
	For transmitted light observation, the coverglass thickness must not exceed 0.17mm	Use a coverglass with thickness 0.17mm
	Focus is not even	Slide holder is not flat. Move the specimen to a flat position
One side of the image is out of focus	The nosepiece is not in the center of the light path	Turn the nosepiece to a click stop
	The specimen is out of place (tilted)	Place the specimen flat on the stage
	The optical performance of the sample cover glass is poor	Use a cover glass of better quality
II. Electric section:		
The LED doesn't turn on	No power supply	Check the power cord connection
The brightness is not enough	The brightness adjustment is low	Adjust the brightness
The light blinks	The power cord is poorly connected	Check the power cord
III. Observation tube:		
Field of view of one eye does not match that of the other	Interpupillary distance is incorrect	Adjust interpupillary distance
	Incorrect diopter adjustment	Adjust diopter
	Your view is not accustomed to microscope observation	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.