B-810 Series

INSTRUCTION MANUAL

Model
B-810

Ver. 2.3 2023

Table of Contents

1.	Warning	3
2.	Safety Information	3
3.	Package content	4
4.	Unpacking	5
5.	Intended use	5
6.	Symbols and conventions	5
7.	Instrument description	6
7.1	Microscope body	6
7.2	Brightfield condenser (M-1189 / M-1191)	8
7.3	Brightfield / Darkfield / Phase Contrast condenser (M-1152.NO)	8
7.4	Brightfield / Phase Contrast condenser (M-1124.NO)	9
8.	Assembling	10
9.	Brightfield observation procedures	12
10.	Use of the microscope	13
10.1	Microscope switch on	13
10.2	Light intensity adjustment	13
10.3	Coarse focus tension adjustment	13
10.4	Coarse upper limit lever	14
10.5	Stage	14
10.6	Dioptric adjustment	14
10.7	Adjusting the interpupillary distance	15
10.8	Use of eye shields	15
10.9	Light path selection	15
10.10	Centering the condenser	16
10.11	Effects of the field diaphragm	16
10.12	Aperture diaphragm	16
10.13	Use of oil immersion objective	17
11.	Brightfield / Darkfield / Phase Contrast Condenser	18
11.1	Brightfield Observation (BF)	18
11.2	Darkfield Observation (DF)	18
11.3	Phase Contrast Observation (PH)	19
11.4	Use of the green filter	20
12.	Brightfield / Phase Contrast Condenser	20
12.1	Brightfield Observation (BF)	21
12.2	Phase Contrast Observation (PH)	21
12.3	Use of the green filter	21
13.	Microphotography	22
13.1	Use of C-mount cameras	22
13.2	Use of Reflex cameras	22
14.	Maintenance	23
15.	Troubleshooting	24
	Equipment disposal	26

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

2. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

3. Package content

- | | |
|--|---|
| ① Microscope frame | ⑦ Allen wrench |
| ② Eyepieces | ⑧ Dust cover |
| ③ Objectives | ⑨ Power supply |
| ④ Observation head | ⑩ Green filter (only for PH configuration) |
| ⑤ Condenser (according to the configuration) | ⑪ Centering telescope (only for PH configuration) |
| ⑥ Immersion oil (if 100x objective is included into the configuration) | |

4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

5. Intended use

Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.


ELECTRICAL SHOCK


This symbol indicates a risk of electrical shock.

7. Instrument description

7.1 Microscope body

Opposite side


7.2 Brightfield condenser (M-1189 / M-1191)

7.3 Brightfield / Darkfield / Phase Contrast condenser (M-1152.NO)

7.4 Brightfield / Phase Contrast condenser (M-1124.NO)

8. Assembling


1. Insert the optical head above the stand and tighten the screw. (Fig. 1)
- **Hold the head with one hand during the locking in order to avoid that the head falls.**

2. Insert both eyepieces into the tubes of the optical head. (Fig. 2)

3. Insert the condenser under the stage. Check that it is correctly inserted in its housing (under the condenser there is a plug that must completely enter into the guide of the condenser holder). (Fig. 3)
4. Tighten the condenser fixing screw ①.


5. Screw each objective into the thread of the nosepiece, clockwise with increasing magnification. (Fig. 4)

-
6. Insert the power supply jack in the connector placed at the rear side of the microscope. (Fig. 5)

9. Brightfield observation procedures

10. Use of the microscope


10.1 Microscope switch on

Move the main switch ①, located on the left side of the stand, to position “1” to turn on the microscope. (Fig. 6)

10.2 Light intensity adjustment


Operate on the light intensity adjustment knob ② to increase / decrease the illumination voltage. (Fig. 7)

10.3 Coarse focus tension adjustment

The tension of the coarse focusing knob is factory preset.


1. To modify the tension according to personal's needs, rotate the ring ③. (Fig. 8)
- Clockwise rotation increases the tension.
 - If the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.

10.4 Coarse upper limit lever

The upper limit knob has two functions: prevent the contact between slide and objective and acts as “focus memory”.

1. After focusing the specimen, pull the lever ① toward the front of the microscope and lock it. (Fig. 9).
 - In this way the focus upper limit is set.
 2. Now one can lower the stage with coarse focus knob, replace the specimen and raise again the stage up to the upper limit: specimen will be in approximate focus and will need a fine adjustment to get the proper focus..
 - **Fine focus movement is not affected by the coarse focus lock.**
 - **To unlock, move the lever in the opposite direction to the one used for the locking.**
- **Two stoppers ② are inserted on the frame. DO NOT REMOVE THE TWO STOPPERS.**

10.5 Stage

Stage accepts standard slides 26 x 76 mm, thickness 1,2 mm with coverside 0,17mm.

It is possible to place two slides side by side on the stage.

- **Open the spring arm of the slide holder ③ and place frontally the slides on the stage. (Fig. 10)**
- **Gently release the spring arm of the slide holder.**
- **A sudden release of the spring arm could cause the falling of the slide.**

10.6 Diopter adjustment

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
 2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the diopter adjustment ring ④ to compensate. (Fig. 11)
- **The adjustment range is ± 5 diopter. The number indicated on the adjustment ring graduation should correspond to the operator's diopter correction.**


10.7 Adjusting the interpupillary distance

Observing with both eyes, hold the two eyepiece prism assemblies. Rotate them around their common axis until the fields of view coincide.

- **The graduation on the interpupillary distance indicator ①, pointed by the spot “.” on the eyepiece holder, shows the distance between the operator’s eyes. (Fig. 12)**

The range of the interpupillary distance is 48-75mm.

10.8 Use of eye shields

- **Use with eyeglasses**

Fold rubber eyeshields with both hands. Folded eyeshields avoid scratching the lenses of eyeglasses. (Fig. 13)

- **Use without eyeglasses**

Raise eye shields and observe at the microscope placing eyes to the shields, avoiding external light to disturb the observation. (Fig. 14)

10.9 Light path selection

The observation head is equipped with an optical path selector that allows the light to be distributed to the eyepieces and to the photo / TV port.

1. Fully insert the selector ② to get 100% of the light to the eyepieces. (Fig. 15)
- In this case the photo / TV port is not enabled.
2. Pull out the selector to enable the light also at the photo / TV port. The light split in this case will be 50% to the eyepieces and 50% to the photo / TV port.

10.10 Centering the condenser

1. Place the specimen on the stage, insert 10x objective into the light path and focus.
2. Insert the front lens of the swing-out condenser ①. (Fig. 16)
3. Rotate the field diaphragm ring ② in counterclockwise direction, to fully close the diaphragm.
4. Rotate the condenser height adjustment knob ③ to focus the edges of the diaphragm.
5. Rotate the two centering screws ④ to bring the bright spot in the center of the field of view.
6. Gradually open the diaphragm. The condenser is centered when the diaphragm image is symmetrical to the field of view. (Fig. 17)
7. In normal use, open the diaphragm until it circumscribes the field of view.

Fig. 16

10.11 Effects of the field diaphragm

Field diaphragm adjusts the illuminated area to obtain a high contrast image.

Set the diaphragm according to the objective in use until it circumscribes the field of view, in order to eliminate unnecessary light to eyepieces.


Fig. 17

10.12 Aperture diaphragm

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image.
- With low contrast specimens set the numerical aperture value ⑤ (printed on the condenser ring) to about 70%-80% of the objective's N.A. (Fig. 18) If necessary, remove on eyepiece and, looking into empty sleeve, adjust the condenser's ring in order to obtain an image like the one in fig. 19.

Example: with objective PLAN 40x / 0,65 set the scale to 0.65 x 0.8 = 0,52

Fig. 18


Fig. 19

10.13 Use of oil immersion objective

1. Focus the specimen with a low power objective.
2. Lower the stage (remembering to lock the coarse upper limit lever).
3. Put a drop of oil (provided) on the area of the specimen to be observed. (Fig. 20)
 - **Make sure that there are no oil bubbles. Air bubbles in the oil damage the image quality.**
 - To check for bubbles: remove an eyepiece, fully open the aperture diaphragm and observe the objective exit pupil. (The pupil must be round and bright).
 - To remove the bubbles, gently move the nosepiece to the right and left to move the immersion objective a few times and allow the air bubbles to move.
4. Insert immersion objective.
5. Return the stage to the upper focusing point and obtain an optimal focus using the fine focus knob.
6. After use, gently remove the oil with a soft paper towel or a lightly moistened optic paper with a mixture of ethyl ether (70%) and absolute ethyl alcohol (30%).
 - **The immersion oil, if not immediately cleaned, could crystallize creating a glass-like layer. In this situation the observation of the specimen would be difficult if not impossible due to the presence of an additional thickness on the objective.**

11. Brightfield / Darkfield / Phase Contrast Condenser

The M-1152.NO condenser allows observation in bright field, dark field and phase contrast.

Fig. 20

Fig. 21


Fig. 22


Fig. 23


Fig. 24

Observation mode	Condenser Turret position
Brightfield	BF (Fig. 20)
Darkfield	DF (Fig. 21)
Phase contrast 10x	10/20 (Fig.22)
Phase contrast 20x	10/20 (Fig.22)
Phase contrast 40x	40 (Fig. 23)
Phase contrast 100x	100 (Fig. 24)

11.1 Brightfield Observation (BF)



1. Rotate the condenser turret to insert the "BF" position.
2. Now repeat the steps described in the chapter "*Brightfield observation procedure*".

11.2 Darkfield Observation (DF)

1. Rotate the condenser turret to insert the "DF" position.
 - **By inserting the darkfield ring, the aperture diaphragm opens automatically. This is a desired effect and should not be considered a defect.**
2. Place a specimen on the stage and focus.
3. Observing into eyepieces raise or lower the condenser until a homogeneous illumination of the specimen can be achieved, thus obtaining a proper darkfield effect.
 - **Darkfield requires a huge amount of light. Switching from darkfield to brightfield, one could be dazzled. Do not keep your eyes on the eyepieces when moving the condenser turret from DF to BF.**
 - **"Dry" darkfield observation, that is, without the use of oil, is only possible with objectives with N.A. lower than 0,7.**
 - **When observing in darkfield, it may be necessary to raise the condenser from the normal position to obtain a more homogeneous illumination. This is not a defect.**

11.3 Phase Contrast Observation (PH)

1. Center the condenser as already described in the chapter 110.10.
 - This condenser does not have a swing-out lens, so the operation described in step 2 is not needed.
2. Rotate the condenser turret to insert the "10/20" position.
3. Insert 10x objective into the light path.
 - **By inserting any phase ring, the aperture diaphragm opens automatically. This is a desired effect and should not be considered a defect.**
4. Place a specimen on the stage and focus.
5. Remove one eyepiece and insert the centering telescope. (Fig. 25)
6. Rotate the upper part of the centering telescope until the two phase rings (one dark and one bright) visible in the telescope are in focus. (Fig. 26)
7. Using centering screws on the condenser ① (Fig. 27), center the phase rings to make the bright ring ② be concentric to the dark ring ③. (Fig. 28)
8. Insert 20x objective (do not rotate the condenser turret) and check the centering of the two rings.
9. Repeat the same operation with other objectives to check the ring centering: 40x objective – turret position "40", 100x objective – turret position "100".
10. At the end remove the centering telescope, reinstall the eyepiece and begin observation.
- **With 40x and 100x objectives it may be useful to slightly raise the condenser, to obtain a better projection of the phase rings. This is not a defect.**
11. **With the 4X objective, the condenser could have a dark halo at the periphery of the field of view. This is not to be considered a defect.**

11.4 Use of the green filter

- The green filter is used to increase the contrast of the image during phase contrast observation.
- Place the filter on the field lens of the microscope (Fig. 29) and begin the observation.
- For observation in brightfield or darkfield it is advisable to remove the filter from the optical path.

Fig. 29

12. Brightfield / Phase Contrast Condenser

Slider condenser M-1124.NO allows observation in bright field and phase contrast with objectives 10X, 20X and 40X.

Fig. 31

Fig. 32

Fig. 33

Observation mode	Slider position
Brightfield	O (Fig. 30)
Phase contrast 10x	10/20 (Fig. 31)
Phase contrast 20x	10/20 (Fig. 31)
Phase contrast 40x	40 (Fig. 32)

12.1 Brightfield Observation (BF)

1. Move the condenser slider to the middle position to insert the empty position. (Fig. 33)
2. Now repeat the steps described in the chapter “*Brightfield observation procedure*”.


12.2 Phase Contrast Observation (PH)

1. Center the condenser as already described in the chapter 9.10.
 - This condenser does not have a swing-out lens, so the operation described in step 2 is not needed.
2. Move the condenser slider to the right to insert the phase ring 10 (Fig. 34) or to the left to insert the phase ring 40. (Fig. 35)
3. Insert the corresponding objective in to the light path.
4. Open the aperture diaphragm.
5. Place a specimen on the stage and focus.
6. Remove one eyepiece and insert the centering telescope. (Fig. 25)
7. Rotate the upper part of the centering telescope until the two phase rings (one dark and one bright) visible in the telescope are in focus. (Fig. 25-26)
8. Using the centering screws on the slider ① (Fig. 36), center the rings as already described in chapter 11.3.
9. At the end remove the centering telescope, reinstall the eyepiece and begin observation.
 - **With 40x and 100x objectives it may be useful to slightly raise the condenser, to obtain a better projection of the phase rings. This is not a defect.**
 - **With the 4X objective, the condenser could have a dark halo at the periphery of the field of view. This is not to be considered a defect.**

12.3 Use of the green filter


- The green filter is used to increase the contrast of the image during phase contrast observation.
- Place the filter on the field lens of the microscope (Fig. 29) and begin the observation.
- For observation in brightfield it is advisable to remove the filter from the optical path.

13. Microphotography

13.1 Use of C-mount cameras

1. Loosen the clamping screw ① on the trinocular port and remove the dust cap ②. (Fig. 37)

2. Screw the C-mount adapter ③ to the camera ④ and insert the round dovetail of the C-mount into the empty hole of the trinocular port, then tighten the clamping screw ①. (Fig. 38)

13.2 Use of Reflex cameras

1. Insert the Reflex adapter ① into the relay tube to the microscope ②.
 2. Screw the "T2" ring ③ (not provided) to the reflex adapter.
 3. Connect the Reflex camera ④ to the "T2" just installed (Fig. 39).
 4. Mount the other end of the relay tube ② into the empty hole of the trinocular port, then tighten the clamping screw. (Fig. 37)
- "T2" ring is not provided with the microscope, but is commercially available.
 - While shooting dark specimens, darken eyepieces and viewfinder with a dark cloth to minimize the diffused light.
 - To calculate the magnification of the camera: objective magnification * camera magnification * lens magnification.
 - **If using an SLR camera, mirror movement may cause the camera to vibrate.**
 - **We suggest lifting the mirror, using long exposure times and a remote cord.**

14. Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

15. Troubleshooting


Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
LED operates, but field of view remains dark.	Power supply is unplugged.	Connect
	Brightness is too low	Set brightness to a proper level
Field of view is obscured or not evenly illuminated	Revolving nosepiece is not correctly engaged.	Make sure that the revolving nosepiece clicks properly into place.
	The turret of the phase contrast condenser is in an incorrect position	Move the turret to a click stop
Dirt or dust is visible in the field of view.	Dirt/dust on the specimen	Clean the specimen
	Dirt/dust on the eyepieces	Clean the eyepieces
Image looks double	Aperture iris diaphragm is stopped down too far.	Open aperture iris diaphragm.
	The condenser is not well centered or it is in a wrong height	Set the condenser according to Koehler settings.
Visibility is poor. • Image is not poor. • Contrast is poor. • Details are indistinct. • Phase contrast is poor	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture iris diaphragm is too closed or too open.	Adjust aperture iris diaphragm.
	Dust or dirt on lenses (condenser, objectives, eyepieces and slide)	Clean thoroughly.
	For transmitted light observation, the coverglass thickness must not exceed 0.17mm	Use a coverglass with thickness 0.17mm
	For phase contrast observation, a brightfield objective is used instead a phase contrast one	Use a phase contrast objective
	Phase rings of objective and condenser are not well centered	Operate on centering screws
	Objective in use is not compatible with condenser phase ring	Use a compatible objective
	Focus is not even	Slide holder is not flat. Move the specimen to a flat position.
One side of the image is unfocused	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Slide is mounted not in a flat position (tilted)	Place the specimen in a flat position on the stage
	Poor quality of the glass slide	Use a glass slide with higher quality
II. Mechanical Section:		
Coarse focus knob is hard to turn	Tension adjustment ring is too tight	Loosen tension adjustment ring
Focus is unstable	Tension adjustment ring is too loose	Tighten tension adjustment ring
III. Electrical Section:		
LED doesn't turn on.	Power supply not connected	Check for proper connection
Brightness is not enough	Brightness setting is too low	Adjust brightness
Light blinks	Power supply not well connected	Check for proper connection

IV. Observation tube:		
Field of view of one eye does not match that of the other.	Interpupillary distance is incorrect.	Adjust interpupillary distance.
	Incorrect diopter adjustment.	Adjust diopter.
	Your view is not accustomed to microscope observation.	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope.
V. Microphotography:		
Image edge is unfocused	To a certain extent it is due to achromatic objectives features	To minimize the problem, set the aperture diaphragm in a proper position
Bright spots appear on the image	Stray light entering in the microscope through eyepieces or camera viewfinder	Cover eyepieces and viewfinder with a dark cloth

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.