**B-810 Series**

**INSTRUCTION MANUAL**

<b>Model</b>
B-810

v 1.0 2018

---

## Table of Contents

1. Warning
  2. Symbols and conventions
  3. Safety Information
  4. Overview
  5. Unpacking
  6. Assembling
  7. Summary of brightfield observation procedures (B-810BF)
  8. Using the microscope
  9. Use of universal condenser for brightfield / darkfield / phase contrast
  10. Use of slide condenser for light field / phase contrast
  11. Microphotography
  12. Maintenance
  13. Troubleshooting
- Equipment disposal

---

## 1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

## 2. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

### **CAUTION**

This symbol indicates a potential risk and alerts you to proceed with caution.

### **ELECTRICAL SHOCK**

This symbol indicates a risk of electrical shock.


## 3. Safety Information

### **Avoiding Electrical Shock**

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

## 4. Overview

## 5. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

## 6. Assembling

Once you open the box, these are the microscope's components:

- ① Condenser
- ② Optical head
- ③ Main body
- ④ Eyepieces
- ⑤ Power supply

- ⑥ Stage
- ⑦ Objectives

Put the microscope stand on a solid table. Insert the optical head above the attachment, using the **2mm** Allen wrench to tighten the screw. (Fig.1)

Fig.1

Insert both eyepieces into the tubes of the optical head. (Fig.2)

Fig.2

Insert the condenser under the stage: position until it is well inserted into its holder (under the condenser there is a pin that must fully enter the guide of the holder). (Fig. 3)


Fig.3

Screw each objective into the thread of the nosepiece, in order of magnification. (Fig.4)


Fig.4

Insert the power supply jack  
on the connector at the rear.  
(Fig.5)


## 7. Summary of brightfield observation procedures (B-510BF)

## 8. Using the microscope

### LED settings

In order to activate the transmitted light illuminator, insert the plug of the external power supply into the socket and turn on the main switch on the side of the main body. Turn the brightness adjustment knob to a brightness suitable for observation. (Fig.6)

### Adjust the observation head

Loosen the lock-screw, turn the observation head to a comfortable position for observation, and then lock the lock-screw. (Fig.7)

### Dioptric adjustment

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the dioptric adjustment ring ① to compensate. (Fig.8)

► **The adjustment range is  $\pm 5$  diopter. The number indicated on the adjustment ring graduation should correspond to the operator's dioptric correction.**

## Adjusting the interpupillary distance

Observing with both eyes, hold the two eyepiece prism assemblies. Rotate them around their common axis until the fields of view coincide. (Fig.9)

## Use of eye shields

- **Use with eyeglasses (Fig. 10)**  
Fold rubber eyeshields with both hands. Folded eyeshields avoid scratching the lenses of eyeglasses.
- **Use without eyeglasses (Fig. 11)**  
Raise eye shields and observe at the microscope placing eyes to the shields, avoiding external light to disturb the observation.

## Stage

Stage accepts standard slides 26 x 76 mm, thickness 1,2 mm with coverside 0,17mm. It is possible to place two slides side by side on the stage.

- **Open the spring arm of the slide holder and place frontally the slides on the stage.**
- **Gently release the spring arm of the slide holder. (Fig. 12)**
- ▶ **A sudden release of the the spring arm could cause the falling of the slide.**

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
### Centering the condenser (Fig.13)

1. Place the specimen on the stage, insert 10x objective into the light path and focus.
2. Insert the front lens of the swing-out condenser ①.
3. Rotate the field diaphragm ring ② in the direction showed by the arrow, to fully close the diaphragm.
4. Rotate the condenser height adjustment knob ③ to focus the edges of the diaphragm.
5. Rotate the two centering screws ④ to bring the bright spot in the center of the field of view.
6. Gradually open the diaphragm. The condenser is centered when the diaphragm image is symmetrical to the field of view.
7. In normal use, open the diaphragm until it circoscribes the field of view.

Fig.13

### CENTERING THE CONDENSER

### Effects of the field diaphragm

Field diaphragm adjusts the illuminated area to obtain a high contrast image. Set the diaphragm according to the objective in use until it circoscribes the field of view, in order to eliminate unnecessary light to eyepieces.

### Aperture diaphragm (Fig. 13-14)

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image.
- With low contrast specimens set the numerical aperture value ① (printed on the condenser ring) to about 70%-80% of the objective's N.A. If necessary, remove on eyepiece and, looking into empty sleeve, adjust the condenser's ring in order to obtain an image like the one in fig. 15.

Fig.14

**Example: with objective PLAN 40x / 0,65 set the scale to  $0.65 \times 0.8 = 0,52$**


Fig.15

## Use of oil immersion objective (Fig. 16)

1. Focus the specimen with a low power objective.
  2. Lower the stage (remembering to lock the coarse upper limit knob).
- ▶ **Put a drop of oil (provided) on the area of the specimen to be observed. Make sure that there are no oil bubbles. Air bubbles in the oil damage the image quality.**
- To check for bubbles: remove an eyepiece, fully open the aperture diaphragm and observe the objective exit pupil. (The pupil must be round and bright).
  - To remove the bubbles, gently move the nosepiece to the right and left to move the immersion objective a few times and allow the air bubbles to move.
3. Insert immersion objective.
  4. Return the stage to the upper focusing point and obtain an optimal focus using the fine focus knob.
  5. After use, gently remove the oil with a soft paper towel or a lightly moistened optic paper with a mixture of ethyl ether (70%) and absolute ethyl alcohol (30%).
- ▶ **The immersion oil, if not immediately cleaned, could crystallize creating a glass-like layer. In this situation the observation of the specimen would be difficult if not impossible due to the presence of an additional thickness on the objective.**

## Focus tension adjustment

Turn the tension-adjust knob to get a suitable tension for the focus system. (Fig.17)

**NOTE:** if the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the friction.


Fig.17

## Focus-stop knob

Loosen the focus-stop knob, rotate the coarse focusing knob to bring the slide into focus with a 4X or 10X objective. Then rotate the focus-stop in order to block the height of the stage. This simplifies the next focusing operations. The focus-stop knob is also useful to avoid accidental contacts between objective and specimen. (Fig.18)



Fig.18

## 9. Use of universal condenser for brightfield/darkfield/phase contrast

Universal condenser provided with B-810PH allows observation in brightfield, darkfield and phase contrast.

OBSERVATION MODE	Condenser Turret position
Brightfield	BF (Fig. 17)
Darkfield	DF (Fig. 18)
Phase contrast (10x)	10/20 (Fig.19)
Phase contrast (20x)	10/20 (Fig.19)
Phase contrast (40x)	40 (Fig. 20)
Phase contrast (100x)	100 (Fig. 21)

### 1) BRIGHTFIELD OBSERVATION (BF)

Rotate the condenser turret to insert the “BF” position. Now repeat the steps described in the procedure “Summary of brightfield observation procedures” at pag.47.

### 2) DARKFIELD OBSERVATION (DF)



1. Rotate the condenser turret to insert the “DF” position.
  2. Open the aperture diaphragm.
  3. Place a specimen on the stage and focus.
  4. Observing into eyepieces raise or lower the condenser until a homogeneous illumination of the specimen can be achieved, thus obtaining a proper darkfield effect.
- **Darkfield requires a huge amount of light. Switching from darkfield to brightfield, one could be dazzled. Do not keep your eyes on the eyepieces when moving the condenser turret from DF to BF.**
  - **“Dry” darkfield observation, that is, without the use of oil, is only possible with objectives with N.A. lower than 0,7.**
  - **Observing in darkfield, it may be necessary to raise the condenser from the normal position to obtain a more homogeneous illumination. This is not a defect.**

### 3) PHASE CONTRAST OBSERVATION (PH)

1. Center the condenser as already described at pag. 50.
  2. Rotate the condenser turret to insert the "10/20" position.
  3. Insert 10x objective into the light path.
  4. Open aperture diaphragm.
  5. Place a specimen on the stage and focus.
  6. Remove one eyepiece and insert the centering telescope. (Fig.24)
  7. Rotate the upper part of the centering telescope until the two phase rings (one dark and one bright) visible in the telescope are in focus. (Fig. 25)
  8. Using centering screws on the condenser ①, center the phase rings to make the bright ring ② be concentric to the dark ring ③.
  9. Insert 20x objective (do not rotate the condenser turret) and check the centering of the two rings. (Fig. 26)
  10. Repeat the same operation with other objectives to check the ring centering: 40x objective – turret position "40", 100x objective – turret position "100".
  11. At the end remove the centering telescope, reinstall the eyepiece and begin observation.
- **With 40x and 100x objectives it may be useful to slightly raise the condenser, to obtain a better projection of the phase rings. This is not a defect.**
  - **With the 4X objective, the condenser could have a dark halo at the periphery of the field of view. This is not to be considered a defect.**

#### Use of the green filter (Fig. 28)

- The green filter is used to increase the contrast of the image during phase contrast observation.
- Place the filter on the field lens of the microscope (Fig. 27) and begin the observation.
- For observation in brightfield or darkfield it is advisable to remove the filter from the optical path.

## 10. Use of slide condenser for light field / phase contrast

Universal condenser provided with B-810PH allows observation in brightfield, darkfield and phase contrast.

Fig.29

Fig.30

Fig.31



OBSERVATION MODE	Condenser Turret position
Brightfield	BF (Fig. 29)
Phase contrast (10x)	10 (Fig.30)
Phase contrast (40x)	40 (Fig.31)

### 1) BRIGHTFIELD OBSERVATION (BF)

Rotate the condenser turret to insert the "BF" position. Now repeat the steps described in the procedure "Summary of brightfield observation procedures".

## 2) OBSERVATION IN PHASE CONTRAST (PH)

1. Center the condenser as described on page 25.
  2. Move the condenser slide to insert the "10" position
  3. Insert the 10x objective into the optical path.
  4. Open the aperture diaphragm.
  5. Place a sample on the table and focus.
  6. Remove an eyepiece and insert the centering telescope (Fig.32)
  7. Rotate the top of the telescope to focus on the rings (one light and one dark) visible in the telescope. (Fig. 33)
  8. Using the centering screws on the slide ① center the rings so that the clear ring ② is concentric to the dark ring ③.
  9. Repeat the operation with the 40x objective: 40x objective - slide position "40".
  10. At the end remove the centering telescope, reposition the eyepiece and start the observation.
- **With the 40x objective, it may be useful to slightly raise the condenser, to obtain a better projection of the phase rings. This is not a defect.**

## 11. Microphotography

### Installing the C-mount adapter

1. Loosen the clamping screw ① on the trinocular port and remove the dust cap ②. (Fig.36)
2. Screw the C-mount adapter ③ to the camera ④ and insert the round dovetail of the C-mount into the empty hole of the trinocular port, then tighten the clamping screw ①. (Fig.37)

Fig.36

Fig.37

### Use of Reflex camera

1. Insert the Reflex adapter ① into the relay tube to the microscope ②.
  2. Screw the "T2" ring ③ (not provided) to the reflex adapter.
  3. Connect the Reflex camera ④ to the "T2" just installed (Fig. 38).
- "T2" ring is not provided along with the microscope, but is commercially available.
  - While shooting dark specimens, darken eyepieces and viewfinder with a dark cloth to minimize the diffused light.
  - To calculate the magnification of the camera: objective magnification \* camera magnification \* lens magnification.
- **If using an SLR camera, mirror movement may cause the camera to vibrate.**  
We suggest lifting the mirror, using long exposure times and a remote cord.

Fig.38

---

## 12. Maintenance

### Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 75 % (non condensing). Use a dehumidifier if needed.

### To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

### Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

### Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.  
Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

**For the best results, use the OPTIKA cleaning kit (see catalogue).**

If you need to send the microscope to Optika for maintenance, please use the original packaging.

## 13. Troubleshooting

Review the information in the table below to troubleshoot operating problems.


PROBLEM	CAUSE	SOLUTION
<b>1. Optical System</b>		
LED does not light.	Power cord is unplugged.	Plug power cord into the power outlet.
LED operates, but field of view remains dark.	Aperture and field iris diaphragms are not opened wide enough.	Adjust them to proper sizes.
	Condenser is lowered too much.	Adjust the condenser height position.
	Light path selector knob is set to the camera position.	Move the knob to the eye position.
Field of view is obscured or not evenly illuminated	Light path selector knob is in an intermediate position.	Set the knob according to the observation method.
	Revolving nosepiece is not correctly engaged.	Make sure that the revolving nosepiece clicks properly into place.
	Condenser is not attached properly.	Re-attach it.
	An objective that falls outside of the condenser's illumination range is used.	Use a condenser to match the purpose.
	Condenser is not properly centered.	Center the condenser.
	Field iris diaphragm is stopped down too far.	Open the field iris diaphragm until it circumscribes the field.
Dirt or dust is visible in the field of view.	Dirt/dust on the eyepieces	Clean thoroughly.
	Dirt or the top surface of the condenser	
	Dirt/dust on the specimen	
Visibility is poor. · Image is not clear. · Contrast is poor. · Details are indistinct. · Image glares.	Condenser is lowered too far.	Adjust the condenser height position.
	Aperture iris diaphragm is stopped down too far.	Open aperture iris diaphragm.
	Front lens of objective is dirty.	Clean objective.
	Immersion oil is not being used with an oil immersion objective.	Use immersion oil.
	Immersion oil contains bubbles.	Remove the bubbles.
	Recommended immersion oil is not used.	Use the provided immersion oil.
	Dirt/dust on specimen.	Clean it.
	Dirt/dust on condenser	
	Inappropriate object side or cover glass thickness.	Replace with glass of recommended thickness.

<b>PROBLEM</b>	<b>CAUSE</b>	<b>SOLUTION</b>
One side of image is blurred.	Objective is not correctly engaged in light path.	Make sure that revolving nosepiece clicks into place correctly.
	Specimen is not correctly mounted on stage.	Place specimen correctly on to of stage and secure it with slide holder.
Image appears to waver.	Objective is not correctly engaged in light path.	Make sure that revolving nosepiece clicks into place correctly.
	Condenser is not properly centered.	Center the condenser.
Field of view becomes only slightly brighter when the voltage is raised.	Condenser is not properly centered.	Center the condenser.
	Condenser is lowered too far.	Adjust the condenser height position.
<b>2. Coarse/Fine Adjustment</b>		
Coarse adjustment knob is hard to turn.	Tension adjustment ring is tightened excessively.	Loose ring.
	You are trying to raise stage while focus-lock lever is kept locked.	Unlock focus-lock lever
Stage drifts down by itself or focus is lost during observation.	Tension adjustment ring is too loose.	Tighten ring.
Coarse adjustment will not go all the way up.	Focus-lock lever is locked at a too low height.	Unlock focus-lock lever.
Coarse adjustment will not go all the way down.	Condenser holder is too low.	Raise condenser holder.
Objective makes contact with specimen before focus is obtained.	Specimen is mounted upside down.	Mount specimen correctly.
<b>3. Observation Tube</b>		
Field of view of one eye does not match that of the other.	Interpupillary distance is incorrect.	Adjust interpupillary distance.
	Incorrect diopter adjustment.	Adjust diopter.
	Different eyepieces are used on left and right.	Change on eyepiece to match the other so that both sides are the same type.
	Your view is not accustomed to microscope observation.	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope.
<b>4. Stage</b>		
Specimen stops midway on the Xaxis traverse.	Specimen is not correctly positioned.	Place specimen correctly.
X- and Y-axis knobs are too tight or too loose.	Tension of X- and Y-axis knobs is too high or too low.	Adjust tension.

---

## Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection.

The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment.

Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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**OPTIKA® S.r.l.**