ECOVISION Series

INSTRUCTION MANUAL

| Model |
|-----------|
| B-20R |
| B-20CR |
| M-100FX |
| M-100FLed |

Ver. 2.3 2023

Table of contents

| | | |
|-------|---|----|
| 1. | Warning | 3 |
| 2. | Safety Information | 3 |
| 3. | Unpacking | 3 |
| 4. | Intended use | 3 |
| 5. | Symbols and conventions | 3 |
| 6. | Instrument description | 4 |
| 6.1 | B-20R | 4 |
| 6.2 | B-20CR | 5 |
| 6.3 | M-100FX | 6 |
| 6.4 | M-100FLed | 7 |
| 7. | Assembling | 8 |
| 7.1 | Microscope assembling | 8 |
| 7.1.1 | Installing the batteries (M-100FLed) | 8 |
| 7.1.2 | Installing the batteries (B-20R / B-20CR) | 8 |
| 8. | Use of the microscope | 9 |
| 8.1 | Switching on the microscope | 9 |
| 8.1.1 | B-20R / B-20CR / M-100FX | 9 |
| 8.1.2 | M-100FLed | 9 |
| 8.2 | Light intensity adjustment | 9 |
| 8.2.1 | M-100FX | 9 |
| 8.2.2 | B-20R / B-20CR / M-100FX | 9 |
| 8.2.3 | M-100FLed | 9 |
| 8.3 | Placing the specimen | 10 |
| 8.3.1 | B-20R / M-100FX / M-100FLed | 10 |
| 8.3.2 | B-20CR | 10 |
| 8.4 | Use with batteries | 10 |
| 8.4.1 | M-100FLed | 10 |
| 8.4.2 | B-20R / B-20CR | 10 |
| 8.5 | Aperture diaphragm | 11 |
| 9. | Maintenance | 12 |
| 10. | Troubleshooting | 13 |
| | Equipment disposal | 14 |

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

2. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

3. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

4. Intended use

Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

5. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.


ELECTRICAL SHOCK


This symbol indicates a risk of electrical shock.

6. Instrument description


6.1 B-20R


6.2 B-20CR

6.3 M-100FX

6.4 M-100FLed


7. Assembling

7.1 Microscope assembling

The microscope comes already assembled from the factory. No assembling procedure is requested.

7.1.1 Installing the batteries (M-100FLed)

- **Batteries are not provided. It is necessary to use 3x AA rechargeable NiMh 1.2V.**
1. Remove the cover (Fig. 1) and insert the batteries, respecting the polarity. (Fig. 2)
 2. Put back the battery compartment cover.

7.1.2 Installing the batteries (B-20R / B-20CR)

- **Batteries are not provided. It is necessary to use 3x AAA rechargeable NiMh 1.2V.**
1. Remove the cover and insert the batteries, respecting the polarity. (Fig. 3)
 2. Put back the battery compartment cover.

8. Use of the microscope

8.1 Switching on the microscope

8.1.1 B-20R / B-20CR / M-100FX

1. Operate on the main switch ① placed in the rear side of the microscope, moving the selector on "I". (Fig. 4)

8.1.2 M-100FLed

1. Operate on the main switch ② placed in front of the microscope. (Fig. 5)

8.2 Light intensity adjustment

8.2.1 M-100FX

- No light intensity adjustment can be done on this model.

8.2.2 B-20R / B-20CR / M-100FX

1. Turn the light intensity dial ③ placed on the right side of the microscope. (Fig. 6)

8.2.3 M-100FLed

1. Operate on the main switch ② placed in front of the microscope. (Fig. 5)

8.3 Placing the specimen

8.3.1 B-20R / M-100FX / M-100FLed


1. Place the specimen on the stage and secure it using the stage clips ①. (Fig. 7).

8.3.2 B-20CR

Stage accepts standard slides 26 x 76 mm, thickness 1,2 mm with coverslide 0,17 mm. (Fig. 8)


1. Open the spring arm of the slide holder ② and place the slide from the front on the stage.
 2. Gently release the spring arm of the slide holder.
- **A sudden release of the spring arm could cause the falling of the slide.**

8.4 Use with batteries

8.4.1 M-100FLed

- The microscope is provided with an ON-OFF switch ③ (Fig. 9) placed on the bottom side of the microscope. The use of the switch is as follows:
1. **Switch placed in “OFF” position:** use of the microscope is possible only with external power supply. Battery use and recharge is disabled.
 2. **Switch placed in “ON” position:** use of the microscope is possible either with external power supply or with batteries. In this case the battery recharge is enabled.
- **Once the batteries are charged, it is recommended to disconnect the power supply, or set the switch to “OFF”.**

8.4.2 B-20R / B-20CR

- The batteries recharge when the power supply is connected to the power outlet. Unplugging the power supply starts the microscope's battery operation.
- **Once the batteries are charged, it is recommended to disconnect the power supply.**

8.5 Aperture diaphragm

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image. Move the diaphragm ring ① (Fig. 10) on the value corresponding to the objective in use. In this case the optimal setting of the condenser is achieved.
- With low contrast specimens set the numerical aperture to about 70%-80% of the objective's N.A. If necessary, remove on eyepiece and, looking into empty sleeve, adjust the condenser's diaphragm in order to obtain an image like the one in Fig. 11.


Fig. 10


Fig. 11

9. Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°- 40°C and a maximum relative humidity of 85 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off- position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.


10. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

| PROBLEM | CAUSE | SOLUTION |
|--|---|--|
| I. Optical Section: | | |
| LED operates, but field of view remains dark. | Power supply is unplugged. | Connect |
| | Brightness is too low | Set brightness to a proper level |
| | Batteries discharged | Charge the batteries |
| Dirt or dust is visible in the field of view. | Dirt/dust on the specimen | Clean the specimen |
| | Dirt/dust on the eyepieces | Clean the eyepieces |
| Image looks double | Aperture diaphragm is stopped down too far | Open aperture diaphragm |
| Visibility is poor. <ul style="list-style-type: none"> • Image is not good. • Contrast is poor. • Details are indistinct. • Image glares | Revolving nosepiece is in an incorrect position | Move the nosepiece to a click stop |
| | Aperture diaphragm is too closed or too open | Adjust aperture diaphragm |
| | Dust or dirt on lenses (condenser, objectives, eyepieces and slide) | Clean thoroughly |
| | For transmitted light observation, the coverglass thickness must not exceed 0.17 mm | Use a coverglass with thickness 0.17 mm |
| | Focus is not even | Slide holder is not flat. Move the specimen to a flat position |
| One side of the image is out of focus. | The nosepiece is not in the center of the light path | Turn the nosepiece to a click stop |
| | The specimen is out of place (tilted) | Place the specimen flat on the stage. |
| | The optical performance of the sample cover glass is poor | Use a cover glass of better quality |
| II. Electric Section: | | |
| The LED doesn't turn on. | No power supply | Check the power cord connection |
| The brightness is not enough | The brightness adjustment is low | Adjust the brightness |
| The light blinks | The power cord is poorly connected | Check the power cord |

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.