Digital Refractometer

INSTRUCTION MANUAL

Model
HRD-300N
HRD-400N
HRD-500N

Ver. 1.3 2023

Table of Contents

1.	Warning	3
2.	Unpacking	3
3.	Intended use	3
4.	Package content	4
5.	Instrument description	5
6.	Technical specifications	6
7.	Display and buttons	7
7.1	Display	7
7.2	Buttons	7
8.	Use of the refractometer	8
8.1	Battery installation	8
8.2	Turning on the instrument	8
8.3	Calibration	8
8.4	Measurement	9
9.	Scales conversion and temperature systems conversion	10
9.1	Scales conversion	10
9.2	Temperature system conversion	10
10.	System shutdown	11
11.	Troubleshooting guide	11
12.	Maintenance	11
	Equipment disposal	12

1. Warning

This device is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

2. Unpacking

The device is stored in a molded plastic container. Carefully open the lid and remove the device from its housing. Be careful not to drop the optical parts when taking the device out of the box to prevent damage.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

3. Intended use

Standard models

For research and teaching use only.

4. Package content


① Digital refractometer

② Dropper

③ Screwdriver

5. Instrument description

6. Technical specifications

HRD-300N

MEASURING RANGE	0-50% Brix 1.3330 - 1.4200nD RI
MIN. DIVISION	0.1 Brix 0.0001nD
ACCURACY	$\pm 0.2\%$ Brix ± 0.0003 nD

HRD-400N

MEASURING RANGE	0-50% Brix 0.0 - 28.0% Salinity 1.3330 - 1.4200nD RI
MIN. DIVISION	0.1 Brix 0.1 Salinity 0.0001nD RI
ACCURACY	± 0.2 Brix ± 0.2 Salinity ± 0.0003 nD RI


HRD-500N

MEASURING RANGE	Urine SP. G 1.000-1.050 Serum P. 0-12 g/dl 1.3330 - 1.3900nD RI
MIN. DIVISION	Urine SP 0.001 Serum 0.1 0.0001nD RI
ACCURACY	± 0.002 Urine ± 0.2 Serum ± 0.0003 nD RI

7. Display and buttons

7.1 Display

The LCD screen has three main display areas, namely: instrument display area, temperature display area, multi-function display area, and others.

- ① Battery level

② Reading area

③ % or ‰ unit
- ④ Refractive index (RI)

⑤ Temperature


⑥ Multi-function

Regarding the battery level, the symbols shown on the display vary according to the charge % of the battery.

Battery level	Symbol
100% - 80%	
80% - 50%	
50% - 20%	
20%	flashing

7.2 Buttons

On the instrument three buttons are available (Fig. 3):


1. **READ:** booting / measurement.
2. **CAL:** for entering calibration mode “ZERO POINT”
3. **SCALE:** for converting different scales / converting temperature system between Celsius and Fahrenheit.

8. Use of the refractometer

8.1 Battery installation

1. Open the battery compartment by pushing the cover's lock along the direction of the arrow. (Fig. 1)
2. Insert the battery, respecting polarity.
3. Close the battery compartment.


8.2 Turning on the instrument

Press **“READ”** button for 1 second. The instrument turns on and starts boot. (Fig. 2)

NOTE:

- The Multi-function display shows the current numerical scale. E.g.: S01 indicates the first scale.
- When used outdoors, it is recommended to avoid direct sunlight in order to minimize the impact on measurement accuracy. Alternatively, install the removable cover and close the lid to better resist the impact of light.
- Before placing the liquid sample, clean the lens with a clean soft cloth or soft paper.
- Keep the instrument in a stable and firm position.

Make sure the instrument, surroundings and samples are at the same temperature before taking the measurement.

8.3 Calibration

1. Drop 4 - 5 drops of distilled water in on the prism.
2. Press **“CAL”** button for 2-3 seconds till see the 'CAL' flashing. (Fig. 3)

3. Press **“CAL”** button once again while 'CAL' is flashing and look at the display (Fig. 4): when **“- End”** appears calibration is finished.

4. The value shown is 0.0%. (Fig. 5)
- If no operation is performed within 10 seconds, the instrument will return to the boot position.

Fig. 5

- If the calibration is not completed, the multi-function display will show an Error code. (Fig. 6)

NOTE:

- If the multi-function area shows code A01, it means the calibration temperature exceeds the limit; other error codes are listed in Chapter 10 "Troubleshooting Guide".
- The instrument only accepts calibration with distilled water.

Fig. 6

8.4 Measurement

1. After calibration, wipe distilled water and dry the sample plate.
2. Drop 4-5 drops of the sample material.
3. Quickly press the "READ" button for 1 second.
4. The instrument, accordingly, will show the actual value. (Fig. 7)
- If the measurement range value is exceeded "HHH" or "LLL" appears on the display. (Fig. 8)


Fig. 7

Fig. 8

9. Scales conversion and temperature systems conversion

9.1 Scales conversion

This instrument offers maximum 10 scale marks measure scale.

1. Press **"SCALE"** button each second can convert the scales and the values. (Fig. 9)


Fig. 9

9.2 Temperature system conversion

This instrument offers two temperature units: Celsius (0.0 ~ 50.0 °C) and Fahrenheit (32.0 ~ 122.0 °F).

1. Press **"SCALE"** button for 2 seconds: temperature unit will be converted. (Fig.10)

Fig. 10

- If the higher or lower limit temperature is exceeded, **"HHH"** or **"LLL"** will appear. (Fig. 11)

Fig. 11

10. System shutdown

If the instrument is idle for 1 minute, it will automatically shut down.

11. Troubleshooting guide

ERROR CODE	MEANING
A01	The system exceeds the temperature limits for calibration: 0.0°C ~ 40.0°C
A02	During calibration, no solution or wrong solution
A03	Hardware failure

12. Maintenance

Device environment

This device is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 75 % (non condensing). Use a dehumidifier if needed.

Cleaning the optics


- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the device to Optika for maintenance, please use the original packaging.

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.