IM Series

INSTRUCTION MANUAL

Model
IM-3
IM-3F
IM-3FL4
IM-3LD2
IM-3LD4
IM-3LD4D

Ver. 5.6 2023

Table of Contents

1.	Warning	4
2.	Safety Information	4
3.	Package content	5
3.1	IM-3	5
3.2	IM-3F / IM-3FL4	6
3.3	IM-3LD2	7
3.4	IM-3LD4	8
3.5	IM-3LD4D	9
4.	Unpacking	10
5.	Intended use	10
6.	Symbols and conventions	10
7.	Instrument description	11
7.1	IM-3	11
7.2	IM-3F	12
7.3	IM-3FL4	13
7.4	IM-3LD2	14
7.5	IM-3LD4	15
7.6	IM-3LD4D	16
8.	Assembling	18
8.1	Installing the objectives	18
8.2	Installing stage extension or mechanical stage	18
8.3	Installing the stage insert	19
8.4	Installing the eyepieces	19
8.5	Installing the condenser	19
8.6	Installing the LED housing	20
8.7	Installing color filters	20
8.8	Installing the PC screen (IM-3LD4D)	21
8.9	Cable connection (IM-3LD4D)	21
8.10	Connecting the power supply	22
8.11	Installing fluorescence	23
9.	Brightfield observation (transmitted light)	27
10.	Use of the microscope in brightfield (transmitted light)	28
10.1	Turning on the microscope	28
10.2	Adjusting the light intensity	28
10.3	Adjusting the coarse focus tension	28
10.4	Diopter adjustment	28
10.5	Adjusting interpupillary distance	29
10.6	Use of eyeshields	29
10.7	Selecting the light path	29
10.8	Stage and stage inserts	30
10.8.1	Installing stage inserts	31
10.9	Aperture diaphragm	31
10.10	Using color filters	32
11.	Use of the microscope in phase contrast (transmitted light)	33
11.1	Installing the phase contrast slider	33
11.2	Phase contrast slider	33
11.3	Centering the phase ring	33
12.	Use of the microscope in RPC (transmitted light)	35
12.1	Installing the RPC slider	35
12.2	RPC slider	35
12.3	RPC observation	36
13.	Fluorescence observation (IM-3F / IM-3FL4)	37
14.	Fluorescence observation (IM-3LD2 / IM-3LD4 / IM-3LD4D)	37
15.	Use of the microscope in fluorescence (reflected light)	38
15.1	Centering the mercury bulb	38
15.2	Centering the field diaphragm	40
15.3	Turning on fluorescence LED	40
15.3.1	Switching fluorescence filter cubes	40
15.3.2	Available fluorescence filter cubes	41
15.4	Installing fluorescence filter	43
15.5	Use of the Anti-glow cap	43
15.6	Filter holder / Shutter	44
15.6.1	Inserting a ND filter	44
15.6.2	Using the filter slider	44

16. Simultaneous observation in Phase Contrast + Fluorescence	45
16.1 IM-3F / IM-3FL4	45
16.2 IM-3LD2 / IM-3LD4 / IM-3LD4D	45
17. Microphotography	46
17.1 Use of C-mount cameras	46
17.2 Use of Reflex cameras	46
18. Maintenance	47
19. Troubleshooting	48
Equipment disposal	50

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

2. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

3. Package content

3.1 IM-3

- ① Microscope body
- ② Condenser
- ③ Filter holder
- ④ LED illuminator
- ⑤ Phase contrast slider
- ⑥ Eyepieces

- ⑦ Centering telescope
- ⑧ Metal insert for stage
- ⑨ Glass insert for stage
- ⑩ Objectives
- ⑪ Power supply
- ⑫ Green filter IF550

3.2 IM-3F / IM-3FL4

- | | |
|-------------------------------------|--|
| ① Microscope body | ⑫ UV shield |
| ② Condenser | ⑬ Fluorescence illuminator |
| ③ LED illuminator | ⑭ HBO lamp housing |
| ④ Fluorescence power supply + cable | ⑮ Microscope power supply |
| ⑤ Fluorescence power cable | ⑯ Phase contrast slider |
| ⑥ Filter holder | <ul style="list-style-type: none"> • Provided with IM-3F • Optional with IM-3FL4 |
| ⑦ Metal insert for stage | ⑰ Centering telescope |
| ⑧ Glass insert for stage | <ul style="list-style-type: none"> • Provided with IM-3F • Optional with IM-3FL4 |
| ⑨ Objectives | ⑱ Anti-glow cap |
| ⑩ Eyepieces | |
| ⑪ Green filter IF550 | |

3.3 IM-3LD2

- ① Microscope body
- ② Condenser
- ③ LED illuminator
- ④ Anti-glow cap
- ⑤ Phase contrast slider
- ⑥ Glass insert for stage
- ⑦ Metal insert for stage
- ⑧ Objectives

- ⑨ Eyepieces
- ⑩ Centering telescope
- ⑪ Green filter IF550
- ⑬ UV shield
- ⑮ Power supply + power cord

3.4 IM-3LD4

- ① Microscope body
- ② Condenser
- ③ LED illuminator
- ④ Objectives
- ⑤ Eyepieces
- ⑥ Anti-glow cap
- ⑦ Glass insert for stage

- ⑧ Metal insert for stage
- ⑨ Mechanical stage
- ⑩ UV shield
- ⑪ Stage insert
- ⑫ Lens cleaning tissue
- ⑬ Dust cover
- ⑭ Power supply + power cord

3.5 IM-3LD4D

- | | |
|--------------------------|-----------------------------|
| ① Microscope body | ⑪ Camera |
| ② Condenser | ⑫ "C" mount |
| ③ LED illuminator | ⑬ Mouse |
| ④ Anti-glow cap | ⑭ Keyboard |
| ⑤ Objectives | ⑮ Monitor |
| ⑥ Eyepieces | ⑯ UV shield |
| ⑦ Glass insert for stage | ⑰ Power supply + power cord |
| ⑧ Metal insert for stage | ⑱ Lens cleaning tissue |
| ⑨ Mechanical stage | ⑳ Dust cover |
| ⑩ Stage inserts | |

4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

5. Intended use

Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.


ELECTRICAL SHOCK


This symbol indicates a risk of electrical shock.

7. Instrument description


7.1 IM-3


7.2 IM-3F


7.3 IM-3FL4


7.4 IM-3LD2


7.5 IM-3LD4

7.6 IM-3LD4D


IM-3LD4D (back side)


8. Assembling

8.1 Installing the objectives

1. Turn the coarse focusing knob ① until the nosepiece reaches its lowest position.
- **For a safe transport, the nosepiece is placed in the lowest position and the tension adjustment collar ② is adjusted to the proper tension when the microscope leaves the factory. (Fig.1)**

2. Screw the lowest magnification objective on the nosepiece from the right side, then turn the nosepiece clockwise. Mount the other objectives in the same way, following the sequence from low to high.
- **Note: the objectives can also be installed through the stage opening. (Fig. 2)**
- Clean the objectives regularly. In inverted microscopes, the objectives are very sensitive to dust.
- To prevent dust and contamination from entering the microscope, cover all the unused holes with dust caps ③. (Fig. 3)
- When operating, use the low power objective (10X) to search and focus the specimen, then switch to higher magnifications.
- When switching between objectives, slowly turn the nosepiece until it clicks. The click means that the objective is in the right position, in the center of the light path.

8.2 Installing stage extension or mechanical stage

- **Stage extension and mechanical stage are optional accessories for some models.**
 - Stage extension can be installed on either side of the stage to enlarge the working surface.
 - **Mechanical stage can only be installed on the right side.**
1. Installing the units: screw the bolts in the fixing holes of the stage, then mount the unit from **below the stage**. (Fig. 4)
- **NOTE: The stage has a series of holes in the underside. To install the mechanical stage it is necessary, starting counting from the front of the microscope, to use the third and fifth holes. By using a different set of holes, the mechanical stage will not be installed correctly.**

8.3 Installing the stage insert

1. When using the glass stage, make sure that the insert is horizontal.
2. Install the stage insert in the stage opening. (Fig. 5)

8.4 Installing the eyepieces

Insert both eyepieces into the tubes of the optical head. (Fig. 6)

8.5 Installing the condenser

1. Insert the condenser illumination unit into the bracket. (Fig. 7)

2. Turn the condenser illumination unit clockwise about 90°, with the "AS" mark of filter holder facing front. Align the screw of the condenser illumination unit and the hole of the holder, then screw the bolt in the hole with the supplied Allen wrench. (Fig. 8)

8.6 Installing the LED housing

- **All models except IM-3LD4D**

1. Insert the connector plug into the connector jack.
2. Push the LED housing gently into the holes of the illumination unit. (Fig. 9)

- **IM-3LD4D**


The LED of the IM-3LD4D is already pre-installed on the illuminator. No operation is therefore necessary. (Fig. 10)

8.7 Installing color filters

- **Remove the filter holder, then install the color filter you need. (Fig. 11 - 12)**

1. Mount the color filter flat as shown in ①, checking that they are not tilted.
- **If the color filter is tilted or otherwise out of place ②, it may fall.**


- Color filters can be stacked in the holder. This allows to install as many filters as needed, as long as the whole thickness is less than 11 mm.

8.8 Installing the PC screen (IM-3LD4D)

- For the camera installation, please refer to chapter 17.1.
- Align the notch ① with the bracket ② on the illuminator. (Fig. 13)

- Once the monitor is fully inserted, tighten the wheel ③ to lock the monitor. (Fig.14)

8.9 Cable connection (IM-3LD4D)

- Connect the USB mouse / keyboard receiver in the slot "A" of the PC stick.
- Plug the USB3.0 cable in the slot "B" of the PC stick.
- Plug the PC stick power supply cable in the slot "C".
- Plug the HDMI cable in the slot "D" of the PC stick. (Fig. 15)

- Plug the microscope LED power supply in the connector "E" of the microscope illuminator. (Fig. 16)

6. Plug the monitor power supply in the connector “F” on the monitor. (Fig. 17)
7. Plug the other side of the HDMI cable in the connector “G” on the monitor. (Fig. 17)

8. Plug the USB3.0 cable in the back side of the camera. (Fig. 18)

8.10 Connecting the power supply

1. Turn the main switch to “O” (OFF) before connecting the power supply.
 2. Insert the jack of the power supply into the power socket of the microscope. (Fig.19)
 3. Plug the power supply into the mains socket. Check for a safe connection.
- **Please use the provided power supply. If lost or damaged, please refer to qualified service.**
 - **Connect the power supply to a grounded (earthed) wall outlet only.**

8.11 Installing fluorescence

(IM-3F and IM3FL4 only)

- Disconnect all electrical cables before installing or replacing the bulb.
- The bulb has an anode and a cathode of different sizes. Respect the polarity during assembly, respecting the bulb dimensions.
- Do not touch the bulb of the lamp with bare hands to leave no traces of grease on the bulb. If this happens, clean the bulb with a soft cloth before turning on the lamp.
- The lamp has an average life of about 200-250 hours: a time counter and a voltage indicator are shown on the lamp power supply. Replace the lamp when the hour count exceeds 250 or if the voltage drops below 4.5A.
- During use, the lamp, the lamp housing and the surrounding environment become hot.
- Before replacing the lamp, switch off the power supply, disconnect all cables and wait for the lamp and the lamp housing to cool.
- After switching on the lamp, wait at least 10-15 minutes before switching it off.
- After switching off the lamp, wait for 5-10 minutes before switching it on again so that the mercury vapors have time to condense.
- The lamp contains ultraviolet radiation that could be harmful to eyes and skin. Always look at the lamp arc through the provided orange screen.
- Fluorescence filters are installed before shipment from the factory. Therefore, no user intervention is required

1. Remove the black plastic cover from the microscope rear.
2. Insert the fluorescence illuminator from the back. In order to ease the insertion, just tilt the assembly at about 45° and move it forward. Fix it using the 3 provided Allen screws. (Fig. 20)

Fig. 20

3. Insert the lamp house and fix it with the Allen screw ①. (Fig. 21)
- Please use the provided power supply. If lost or damaged, please refer to qualified service.
 - Connect the power supply to a grounded wall outlet only.

Fig. 21

4. Remove one of the knurled knobs from the filter holder and slide it into the slot placed in the back side of the microscope. (Fig. 22)
5. Once inserted, reinstall the knurled knob.


- **IM-3F**
6. Screw the terminal with the etched **G** on top of the filter lever. Repeat the same steps for the right side, mounting the **B** filterset. (Fig. 23)


- **IM-3FL4**
7. Screw the filter lever on the left side of the microscope. (Fig. 24)


8. In order to prevent possible damages from UV radiation, mount the orange protection screen as shown. (Fig. 25)


9. Open the lamp housing using the door lock screw ① and remove the lamp holder. (Fig. 26)

10. Remove the plastic block ② from the lamp holder (or the exhausted lamp in case of replacement) by loosening the two locking screws ③. (Fig. 27)

11. Insert the mercury bulb ④ (respect the polarity of the bulb), tighten the locking screws and refit the lamp holder inside the lamp housing. (Fig. 28)
- **Do not touch the bulb with bare hands.**

12. Connect the cable from the external power supply to the lamp housing. (Fig. 29)

13. Connect the power cable to the external power supply. (Fig. 30)

- The input voltage for fluorescence power supply is 110-240Vac.
- Please use the standard power cable provided. Select suitable one when missing or damaged.


- Connect the power supply correctly, be sure to have a good earth connection.
- Before connecting the power cord, secure the lamp housing cable to the power supply.
- If the power cord is connected before, there may be a risk of electrical shock.

- Disconnect all electrical cables before installing or replacing the bulb.
- The bulb has an anode and a cathode of different sizes. Respect the polarity during assembly.

9. Brightfield observation (transmitted light)

10. Use of the microscope in brightfield (transmitted light)

10.1 Turning on the microscope

Move the main switch ① in the "I" (ON) position. (Fig. 31)

Fig. 31

10.2 Adjusting the light intensity

Turn the brightness adjustment knob ② placed on the left side of the microscope to increase and decrease the brightness. (Fig. 32)

Fig. 32

10.3 Adjusting the coarse focus tension

- The coarse focusing knob ④ is pre adjusted to a tight tension upon leaving the factory.

If the nosepiece drops down by itself, or the specimen defocuses while adjusting the fine focus knob ⑤, the coarse focus knob is too loose. Turning the tension adjustment collar ④ in clockwise direction tightens the coarse focus tension ③. Rotate in the opposite direction to decrease the tension. (Fig. 33)


Fig. 33

10.4 Diopter adjustment

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the diopter adjustment ring ⑥ to compensate. (Fig. 34)

- The adjustment range is ± 5 diopter. The number indicated on the adjustment ring graduation should correspond to the operator's diopter correction.

Fig. 34

10.5 Adjusting interpupillary distance

Observing with both eyes, hold the two eyepiece prism assemblies. Rotate them around their common axis until the fields of view match.

- **The graduation on the interpupillary distance indicator ①, pointed by the spot “.” on the eyepiece holder, shows the distance between the operator’s eyes. (Fig. 35)**

The range of the interpupillary distance is 48-75mm.

10.6 Use of eyeshields

- **Use with eyeglasses**

Fold rubber eyeshields with both hands. Folded eyeshields avoid scratching the lenses of eyeglasses. (Fig. 36)

- **Use without eyeglasses**

Raise eye shields and observe at the microscope placing eyes to the shields, avoiding external light to disturb the observation. (Fig. 37)

10.7 Selecting the light path

Pull the light path selector lever ② sideways using your thumb, selecting the light path you need. (Fig. 38)

- The tables on the next page show the light distribution depending on the microscope model.

• **IM-3**

LIGHT PATH LEVER	BRIGHTNESS	APPLICATION
In	50% used for binocular observation / 50% used for video or photography	Binocular observation and micrography or video can be operated simultaneously
Out	100% used for binocular observation / 0% used for video or photography	Binocular observation

• **IM-3F / IM-3FL4 / IM-3LD2 / IM-3LD4 / IM-3LD4D**

LIGHT PATH LEVER	BRIGHTNESS	APPLICATION
In	0% used for binocular observation / 100% used for video or photography	Microphotography or video acquisition
Out	100% used for binocular observation / 0% used for video or photography	Binocular observation

10.8 Stage and stage inserts


- **For the best image quality, use flasks, Petri dishes and slides with a 1.2 mm thickness.**



1. Place the proper insert for your specimen (according to the table on the right) on the stage, and fix it with the stage clip.
2. Turning the X and Y knobs, move the specimen to the required position. (Movement Range: 120 (width) × 78 (length) mm).

Moving the specimen

Move the specimen to the desired position by freehand or by turning the knobs ① of the mechanical stage. (Fig. 39)


- **When switching objectives, take care not to touch the holder plates with the objectives, as their weight may damage the front lens.**


	M-793.1 Holder for Petri diameter 38 mm (holder for Terasaki needed)
	M-793.2 Holder for Terasaki and Petri diameter 65 mm
	M-793.3 Holder for slide and Petri diameter 54 mm
	M-793.4 Holder for 2+2 slides
	M-793.6 Holder for Utermöhl-Chamber (holder for Petri diameter 54 mm needed)
	M-793.7 Load-bearing side extension

10.8.1 Installing stage inserts

1. Install the holder in the mechanical stage. (Fig. 40)

2. Multi well plates can be directly inserted in the mechanical stage. (Fig. 41)


10.9 Aperture diaphragm

The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image.

With low contrast specimens move the Aperture Diaphragm lever (AS) ① to set the numerical aperture to about 70%-80% of the objective's N.A. (Fig. 42)

If necessary, remove one eyepiece and, looking into empty sleeve, adjust the condenser's ring in order to obtain an image like the one in fig. 43.

10.10 Using color filters

Selecting the appropriate color filter according your need. (Fig. 44)

FILTER	USE
Green (IF550)	Phase contrast microscopy

11. Use of the microscope in phase contrast (transmitted light)

11.1 Installing the phase contrast slider


1. Insert the slider into the illumination system, printed face up. (Fig. 45)
2. Pull the slider into the desired position, until it arrives to the click stop.
3. When in phase contrast observation, keep the aperture diaphragm adjustment lever ① on the "O" (open) position.

11.2 Phase contrast slider

- The phase ring is pre-centered when the microscope leaves the factory. It should therefore need no further adjustment. Should a re-centering is needed, it can be performed via the two side bolts (see chapter 10.3).
- The 4x/10x position ② must be used with 4x and 10x phase contrast objectives, the 20x/40x position ③ with the 20x and 40x and the SL position ④ is used for brightfield. (Fig. 46)

SLIDER POSITION	MEANING	APPLICATION
SL	empty hole	brightfield observation
4x/10x	phase ring 4x/10x	phase contrast observation with 4x and 10x objectives
20x/40x	phase ring 20x/40x	phase contrast observation with 20x and 40x objectives



11.3 Centering the phase ring

- **Usually this operation is not needed. If necessary, please proceed with the following steps:**
1. Place a specimen on the stage and focus it.
 2. Take out the eyepiece from the eyepiece tube without the diopter adjustment, and replace it with the centering telescope (CT). (Fig. 47)
 3. Check that the phase ring and the objective match, and that both are steadily set on a click stop.

4. Use the CT to focus the condenser phase ring (bright) ① and the objective phase ring (dark) ② image. If the bright phase ring's image is not sharp, adjust the CT head until you can see a clear image of the phase ring. (Fig. 48)
5. Adjust the bolts of the two centering holes in the phase contrast slider using a screwdriver until the bright ring and the dark ring match. (Fig. 49)
6. The 4x and the 10X phase contrast objectives use the same ring on the phase contrast slider. The coincidence of the phase ring center and the phase contrast center must be verified with both objectives. (Fig. 50)
 - **If the phase ring is incorrectly centered, the contrast will be severely impaired.**
 - **The phase ring may need recentering during and after observation of very thick specimens.**
 - **The phase ring may show an apparent misalignment if the specimen is not flat.**

12. Use of the microscope in RPC (transmitted light)

Relief phase contrast (RPC) is a modification of conventional phase contrast that leads to visible improvements in image quality in optical microscopy. Specifically, the following parameters can be improved: contrast, focal depth, sharpness, three-dimensionality, flatness, and halo artifacts. These effects can be achieved when the phase rings of the condenser are replaced by slit rings.

Similar to phase contrast observation, RPC observation requires the use of a slider containing slit phase rings and dedicated RPC objectives.

The use of the slider and objective are identical to those for phase contrast.

12.1 Installing the RPC slider

1. Insert the slider into the illumination system, printed face up. (Fig. 51)
2. Pull the slider into the desired position, until it arrives to the click stop.
3. When in RPC observation, keep the aperture diaphragm adjustment lever ① on the "O" (open) position.

Fig. 51

12.2 RPC slider

- Two sliders are available for the use with different objectives.
- One slider is dedicated to 4X objective (Fig. 52) and another is for 10X/20X/40X objectives. (Fig. 53)
- Both have an empty hole and a RPC ring.

SLIDER POSITION	MEANING	APPLICATION
EMPTY	empty hole ②	brightfield observation
4x	RPC ring 4x ③	RPC observation with 4x objective
10x/20x/40x	RPC ring 10x/20x/40x ④	RPC observation with 10x, 20x and 40x objectives


Fig. 52



Fig. 53

12.3 RPC observation


- **RPC rings don't need a centering.**
 1. Place a specimen on the stage and focus it.
 2. Check that the RPC ring and the objective match, and that both are steadily set on a click stop.
 3. While observing in the eyepiece, modulate the contrast of the sample by turning the ring nut mounted on the slider. (Fig. 54)
- The image will take on a different three-dimensional effect depending on the position of the slit.
- **When using the RPC with the 4x objective, the use of the diffuser filter is recommended (see section 10.8.8).**

13. Fluorescence observation (IM-3F / IM-3FL4)

14. Fluorescence observation (IM-3LD2 / IM-3LD4 / IM-3LD4D)

15. Use of the microscope in fluorescence (reflected light)

15.1 Centering the mercury bulb


Only IM-3F / IM-3FL4

- **Wait around 5 minutes before proceeding with this operation to allow the bulb to properly warm up.**


1. Turn on the HBO mercury bulb by operating the power supply switch ①. (Fig. 55)
2. Turn the nosepiece into an empty position (without objectives) and remove the protective cap, or remove an objective from the nosepiece.


3. Place a piece of white paper on the stage and insert the filter cube "B" into the optical path. (Fig. 56)



4. Acting on the focus screw of the collector lens ② and on the centering screws ③ try to obtain the light spot of the bulb's arc. (Fig. 57-58)


5. Using the focus screw of the collector lens ②, put the image of the arc projected onto the paper. The light spot must be brighter and sharper as possible. (Fig. 58)

6. Using the centering screws ③ on the side of the lamp housing, center the image of the arc. (Fig. 59-60)

7. Using the focusing screw of the collector lens ② enlarge the image until a homogeneous illumination is achieved. (Fig. 61). At this point, insert an objective into the optical path and, looking into the eyepieces, optimize the illumination always using the screws ② and ③.


8. After replacing the exhausted bulb, reset the time counter on the power supply by pressing the "Reset" button ①. (Fig. 62)

15.2 Centering the field diaphragm


Only IM-3F / IM-3FL4

1. Place the specimen on the stage, insert 10x objective into the light path and focus.
2. Rotate the field diaphragm lever ①, to fully close the diaphragm. (Fig. 63)
3. Rotate the two centering screws ② to bring the bright spot in the center of the field of view.
4. Gradually open the diaphragm. The field diaphragm is centered when the diaphragm image is symmetrical to the field of view.
5. In normal use, open the diaphragm until it circoscribes the field of view.

Effects of the field diaphragm

Field diaphragm adjusts the illuminated area to obtain a high contrast image. Set the diaphragm according to the objective in use until it circoscribes the field of view, in order to eliminate unnecessary light to eyepieces. (Fig. 64)

15.3 Turning on fluorescence LED

• IM-3F / IM-3FL4

Turn on the power supply for the mercury bulb and wait 5 minutes for the arc to stabilize. (Fig. 55)

• IM-3LD / IM-3LD4 / IM-3LD4D


Move the main switch in the position "I" to turn on the microscope, then rotate the adjustment knob ③ (Fig. 65) to adjust the fluorescence light intensity.

15.3.1 Switching fluorescence filter cubes

• IM-3F

Move the filter levers (placed on the right or on the left sides of the microscope) ④ to insert the desired filter block (B or G). (Fig. 66)

• **IM-3FL4**

Move the filter lever (placed on the left side of the microscope) ⑤ to insert the desired filter cube: B, G (V and UV - optional). (Fig. 67)

• **IM-3LD2 / IM-3LD4 / IM-3LD4D**

Move the filter lever (placed on the left side of the microscope) ⑥ to insert the desired filter cube (see tables below). (Fig. 68)

15.3.2 Available fluorescence filter cubes

• **IM-3F**

FILTER NAME	EXCITATION FILTER	DICHROIC MIRROR	EMISSION FILTER	APPLICATIONS
B	460-495 nm	505 nm	510LP nm	<ul style="list-style-type: none"> FITC: fluorescent antibodies Acridine orange: DNA, RNA Auramine
G	527-553 nm	565 nm	575LP nm	<ul style="list-style-type: none"> Rhodamine, TRITC: fluorescent antibodies Propidium iodide: DNA, RNA RFP

• **IM-3LD2**

FILTER NAME	EXCITATION FILTER	DICHROIC MIRROR	EMISSION FILTER	APPLICATIONS
B	455-495 nm	500 nm	510LP nm	<ul style="list-style-type: none"> FITC: fluorescent antibodies Acridine orange: DNA, RNA Auramine
G	510-550 nm	570 nm	575LP nm	<ul style="list-style-type: none"> Rhodamine, TRITC: fluorescent antibodies Propidium iodide: DNA, RNA RFP

• **IM-3FL4**

FILTER NAME	EXCITATION FILTER	DICHROIC MIRROR	EMISSION FILTER	APPLICATIONS
B	460-495 nm	505 nm	510LP nm	<ul style="list-style-type: none"> • FITC: fluorescent antibodies • Acridine orange: DNA, RNA • Auramine
G	527-553 nm	565 nm	575LP nm	<ul style="list-style-type: none"> • Rhodamine, TRITC: fluorescent antibodies • Propidium iodide: DNA, RNA • RFP
UV	325-375 nm	415 nm	435LP nm	<ul style="list-style-type: none"> • Nuclear counterstaining
V	390-420 nm	440 nm	455LP nm	<ul style="list-style-type: none"> • Acridine orange: DNA, RNA



• **IM-3LD4 / IM-3LD4D**

FILTER NAME	EXCITATION FILTER	DICHROIC MIRROR	EMISSION FILTER	APPLICATIONS
M-1233	325-375 nm	415 nm	435LP nm	<ul style="list-style-type: none"> • DAPI • Hoechst
M-1233.1	340-390 nm	405 nm	420-470 nm	<ul style="list-style-type: none"> • DAPI • Hoechst
M-1232	390-420 nm	440 nm	450LP nm	<ul style="list-style-type: none"> • Coumarin • Pacific Blue
M-1230	455-495 nm	500 nm	510LP nm	<ul style="list-style-type: none"> • FITC • GFP
M-1230.1	455-495 nm	500 nm	518-542 nm	<ul style="list-style-type: none"> • FITC • GFP
M-1231	510-550 nm	570 nm	575LP nm	<ul style="list-style-type: none"> • Rhodamine • TRITC • Propidium Iodide
M-1231.1	510-550 nm	570 nm	585-625 nm	<ul style="list-style-type: none"> • Rhodamine • TRITC • Propidium Iodide
M-1238	582-603 nm	610 nm	615-645 nm	<ul style="list-style-type: none"> • Red#2 FISH
M-1234	590-650 nm	660 nm	665LP nm	<ul style="list-style-type: none"> • Cy5
M-1235	595-645 nm	655 nm	665-715 nm	<ul style="list-style-type: none"> • Cy5 • AlexaFluor 647 • Draq5
M-1236	623-678 nm	685 nm	690-750 nm	<ul style="list-style-type: none"> • Cy5.5
M-1237	720-760 nm	770 nm	780LP nm	<ul style="list-style-type: none"> • Li_Cor IR dye 800

15.4 Installing fluorescence filter

Only IM-3LD4 / IM-3LD4D

- **This procedure can be needed when installing new fluorescence cubes or when replacing an existing cube with a different one.**
1. Disconnect the power supply plug from the microscope.
 2. Unscrew the filter selector lever.
 3. Open the side cover of the illuminator, by unscrewing the side screws ①. (Fig. 69)
 - The cubes are mounted on the opposite side of the cover: opening the left cover acts on the right side of the slider and vice versa.
-
4. Remove the filter slider from the microscope and put it on a table.
 5. Using the screws provided with the fluorescence cube, fix the new cube on the slider. (Fig. 70)
 6. Reinstall the filter slider on the microscope.
 7. Close the side cover.
-
8. Apply the adhesive marker ② for the fluorescence cube on the side cover. (Fig. 71)
 9. Connect the power supply.
 10. Start working.

15.5 Use of the Anti-glow cap


Use the anti-glow cap to prevent flare coming from the condenser front lens. (Fig. 72)

15.6 Filter holder / Shutter

Only IM-3F / IM-3FL4

- The microscope is equipped with a filter holder / shutter located on the back side of the microscope. (Fig. 73)
- The filter holder has three positions: ① filter holder for ND filters, ② empty hole, ③ shutter.


15.6.1 Inserting a ND filter

The bleaching of the fluorescent specimen caused by excessive energy of the HBO bulb can be delayed by inserting an ND (Neutral Density) filter in the optical path.

1. Slide the filter into the hole and place it. (Fig. 74)

FILTER	USE
ND25	Transmission of 25% of the total amount of light coming from the HBO bulb
ND50	Transmission of 50% of the total amount of light coming from the HBO bulb

15.6.2 Using the filter slider

- Moving the slider toward left or right to insert the desired position. A click sound will confirm the proper positioning of the slider.
1. When the position ① is inserted (provided that a ND filter is installed) the light coming from the lamp housing will be decreased, according to the filter type installed.
 2. When the position ② is inserted, no filters are engaged into the light path.
 3. When the position ③ is inserted, the shutter is engaged in the light path and therefore no light will come.
- We suggest to close the shutter by interrupting the observation for a limited time and not subjecting the sample to unnecessary lighting in the period in which it is not observed. (Switching off and on frequently the HBO lamp considerably reduces its duration).

16. Simultaneous observation in Phase Contrast + Fluorescence

- **Fluorescence models allow observation in transmitted light Phase contrast in combination with reflected light Fluorescence. Samples with rapid decay must first be observed in Fluorescence and then in Phase Contrast. The combined observation allows you to easily identify some areas of the sample that emit fluorescence.**

16.1 IM-3F / IM-3FL4

1. Turn on the power supply for the HBO fluorescent bulb and wait 5 minutes before the arc stabilizes.
2. Move the filter selector to an empty position or, if the filter holder is full, to the position containing the UV filter.
3. Insert the desired PH objective and move the phase contrast slider to the position containing the corresponding phase ring.
4. Focus the sample.
5. Adjust the light intensity of the transmitted light.
6. Move the fluorescence filter selector to the desired position.
7. To obtain the proper observation of the sample, adjust the light intensity of the transmitted light to modulate the intensity of the fluorescence with the one of the phase contrast.


16.2 IM-3LD2 / IM-3LD4 / IM-3LD4D

1. Turn on the microscope main switch.
2. Move the filter selector to an empty position or, if the filter holder is full, to the position containing the UV filter.
3. Insert the desired PH objective and move the phase contrast slider to the position containing the corresponding phase ring.
4. Focus the sample.
5. Adjust the light intensity of the transmitted light.
6. Move the fluorescence filter selector to the desired position.
7. Adjust the light intensity of the reflected light.
8. To obtain the proper observation of the sample, adjust the light intensity of the transmitted light to modulate the intensity of the fluorescence with the one of the phase contrast.

17. Microphotography

17.1 Use of C-mount cameras

1. Loosen the clamping screw ① on the trinocular port and remove the dust cap ②. (Fig. 75)
2. Screw the C-mount adapter ③ to the camera ④ and insert the round dovetail of the C-mount into the empty hole of the trinocular port, then tighten the clamping screw ①. (Fig. 76)

17.2 Use of Reflex cameras

1. Insert the Reflex adapter ② into the relay tube to the microscope ①.
 2. Screw the "T2" ring ③ (not provided) to the reflex adapter.
 3. Connect the Reflex camera ④ to the "T2" ring just installed (Fig. 77).
 4. Mount the other end of the relay tube ① into the empty hole of the trinocular port, then tighten the clamping screw. (Fig. 75)
- "T2" ring is not provided along with the microscope, but is commercially available.
 - While shooting dark specimens, darken eyepieces and viewfinder with a dark cloth to minimize the diffused light.
 - To calculate the magnification of the camera: objective magnification * camera magnification * lens magnification.
 - **If using an SLR camera, mirror movement may cause the camera to vibrate.**
 - **We suggest lifting the mirror, using long exposure times and a remote cord.**

18. Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 75 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off- position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send the microscope to Optika for maintenance, please use the original packaging.

19. Troubleshooting


Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
I. Optical Section:		
The illumination is on, but the field of view is dark	The plug of the LED holder is not connected to the illumination set	Connect them
	The brightness is too low	Adjust to a proper setting
	Too many colour filters have been stacked	Minimize the number of the filters
	Fluorescence filter selector is not in a click stop	Move the selector to a click stop
	Fluorescence shutter is closed	Open the shutter
	Fluorescence filter is not suitable for the specimen	Use a suitable filter
The edge of the field of view is vignett-ed or the brightness is asymmetric	The nosepiece is not in the correct position	Turn the nosepiece to a click stop
	The color filter is partially inserted	Insert the filter to full depth
	The phase contrast slider is not in the proper position	Move the slider to a click stop
Dust and stains can be seen in the field of view	There are stains and dust on the specimen	Clean the specimen
	There are stains and dust on the eyepiece	Clean the eyepiece
Image looks double	The size of the aperture diaphragm is too small	Open the aperture diaphragm
Poor image quality: <ul style="list-style-type: none"> • The image is not sharp • The contrast is not high • The details are not clear • The phase contrast is low 	The nosepiece is not in the center of the light path	Turn the nosepiece to a click stop
	The aperture diaphragm in the view of field is opened too much or too little	Adjust the aperture diaphragm
	The lenses (condenser, objective, eyepieces are culture dish) is dirty	Thoroughly clean all the optical system
	In phase contrast observation, the bottom thickness of the sample is more than 1.2mm	Use a sample holder whose bottom thickness is 1.2mm
	A bright field objective is used for phase contrast observation	Switch to a phase contrast objective
	The condenser ring is not aligned with the objective phase ring	Adjust the condenser ring to match the objective phase ring
	The light ring and/or the phase contrast ring is not centered	Adjust the bolts to center them
	The objective used is not compatible with the phase ring	Please use a compatible objective
	The phase contrast depends on the sample position	The sample holder is not flat. Move the sample around until a compatible area is found.
One side of the image is out of focus	The nosepiece is not in the center of the light path	Turn the nosepiece to a click stop
	The specimen is out of place (tilted)	Place the specimen flat on the stage
II. Mechanical Section:		
The coarse focus knob is hard to turn	The tension adjustment collar is too tight	Loosen the tension adjustment collar
The focus is unstable	The tension adjustment collar is too loose	Tighten the tension adjustment collar

III. Electrical Section:		
The LED doesn't turn on.	No power supply	Check the power cord connection
The brightness is not enough	The brightness adjustment is low	Adjust the brightness
The light blinks	The power cord is poorly connected	Check the power cord
IV. Observation tube:		
The field of view of the two eyes is different	The interpupillar distance is not correct	Adjust the interpupillar distance
	The dioptic correction is not right	Adjust the dioptic correction
	The viewing technique is not correct, and the operator is straining the eyesight	When look into the objective, do not stare at the specimen but look at the whole field of view. Periodically, move the eyes away to look at a distant object, then back into the objective
V. Microphotography and video:		
The image is unfocused	Incorrect focussing	Adjusting the focus system as in the present manual
The edge of the image is unfocussed	To some degree, it is inherent to the nature of achromatic objectives	The problem can be minimized by a correct setting of the aperture diaphragm
Bright patches appear on the image	Stray light is entering the microscope through the eyepieces and through the camera viewfinder	Cover the eyepieces and the viewfinder with a dark cloth

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. “According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal.”

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.