---

## INSTRUCTION MANUAL

Model
M-666
ST-666
M-1144

Version: 2  
Issued: 04, 01, 2017

---

**Application:**

This equipment offers temperature control for microscope stages, it is mainly applied in the biological field in case of special request of constant temperature (inspection of living cells).

**Features:**

- It can be installed on several microscope models.
- Quick reaching of constant temperature: after plugging in the power and setting the temperature, it takes 1-2 minutes to heat the stage and reach the constant temperature in the desired range of temperature

**Specification:**

- Input: selectable 110V or 220V AC (50/60HZ)
- Power: 50W
- Temperature controlling range: from 20° C up to 50° C
- Temperature setting division: 1°C
- Temperature displaying division: 0,1°C
- Accuracy (between actual stage temperature and displayed temperature): +/- 1°C

**Standard outfit:**

- Temperature Controller: 2 LED digital displays and two buttons for temperature adjusting, ON/OFF switch at the rear.
- Temperature Controlled stage: built-in heater and heat sensor with size 150mm x 130mn x 5.8mm
- Electric cable: 1 pc
- Spare fuses: 2 pcs

**Installation:**

- If it is not pre-installed on stage, unscrew the specimen clamp and fix the temperature-controlled surface on the microscope stage.
- Plug the cable from the temperature-controlled surface to the connector on the rear of controller.
- Check that the input voltage marked on the rear of controller is in accordance with the local voltage. 110V or 220V can be changed with appropriate selector if necessary.
- Plug the controller to voltage line with electric cable. Turn on the switch at the rear of controller. The stage is ready to work.

**Use:**

In order to set the desired temperature on stage surface (suitable to the specimen in use), set the working temperature by using the “**Temp Setting**” buttons on the controller panel.

You can set any temperature from 20°C up to 50°C, then wait a few minutes until it displays a stable temperature value.

**Warning:**

- Check the input voltage: the input voltage marked on the rear of controller should be in accordance with the local voltage

## Equipment disposal

Art.13 DLsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection.

The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment.

Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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## **OPTIKA® S.r.l.**