ACCESSORIES Series

# INSTRUCTION MANUAL

Model
M-905

Ver. 1.0 2024

---


## Table of contents

1.	Warning	3
2.	Package content	3
3.	Unpacking	3
4.	Intended use	3
5.	Assembling	4
6.	Use of the system	6
6.1	Light extinction	6

## 1. Warning

This system is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

## 2. Package content

① Side covers

② Rotatable analyzer

③ Polarizer

④ Allen wrench

## 3. Unpacking

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

## 4. Intended use

### Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

### IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

## 5. Assembling

- **NOTE: This accessory can only be used with the IM-3, IM-300 and IM-300D models.**

1. Using the Allen wrench provided, unscrew the screw ① securing the side covers. (Fig. 1)

2. Remove the left side cover mounted on the microscope and insert the side cover ② with slot. (Fig. 2)  
3. Secure the side cover using the screw disassembled earlier.  
4. Repeat the operation with the cover placed on the right side of the microscope.


5. Insert the slider with the rotatable analyzer from the left side. (Fig. 3)

6. Insert the polarizing filter into the supplied filter slider. (Fig. 4)

- 
7. Insert the filter slider into the slot placed in the condenser.  
(Fig. 5)

## 6. Use of the system

### 6.1 Light extinction

1. Remove the slide from the light path and insert 10X objective.
2. Insert rotatable analyzer in the light path. (Fig. 6)
  - When the analyzer is in the correct position, a “click” is heard, indicating that the analyzer is perfectly positioned in the center of the optical path.
3. Insert the filter slider containing the polarizer into the optical path. (Fig. 5)
4. While looking inside the eyepieces, rotate the analyzer dial ① to obtain total extinction (total dark in the eyepieces).
5. Once the dark is achieved (“extinction” or “Crossed Nicol” position) it is possible to begin the observation.
  - **NOTE: In this case, since the polarizer is positioned within the optical path and cannot be adjusted in any way, the graduation shown on the analyzer dial should need not be taken into account. The dial is only used to rotate the analyzer to achieve light extinction.**


Fig. 6

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**OPTIKA® S.r.l.**