SZO Series

INSTRUCTION MANUAL

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1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

2. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

3. Safety Information

Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

4. Intended use


For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

5. Overview


5.1 SZO-B / SZO-T

5.2 SZO-1 / SZO-2


5.3 SZO-3 / SZO-4


5.4 SZO-5 / SZO-6

5.5 SZO-7 / SZO-8

5.6 SZO-9 / SZO-10

6. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

7. Assembly

Once opened the microscope box, the microscope parts are the following:

7.1 SZO-B / SZO-T

- ① Microscope body
(SZO-B Binocular) / SZO-T Trinocular)
- ② Eyepieces

- ③ Allen wrench (only for SZO-T)
- ④ Dust cover
- ⑤ Photo tube (only for SZO-T)


7.2 SZO-1 / SZO-2

- ① Microscope body
(SZO-1 Binocular) / SZO-2 Trinocular)
- ② Eyepieces
- ③ Microscope stand
- ④ Stage clips (one pair)

- ⑤ Photo tube (only for SZO-2)
- ⑥ Black/white plate
- ⑦ Dust cover
- ⑧ Allen wrench (only for SZO-2)

7.3 SZO-3 / SZO-4

- ① Microscope body
(SZO-3 Binocular) / SZO-4 Trinocular)
- ② Eyepieces
- ③ Microscope stand
- ④ Stage clips (one pair)
- ⑤ Photo tube (only SZO-4)
- ⑥ Transparent plate
- ⑦ Incident light LED illuminator
- ⑧ Black/white plate
- ⑨ Dust Cover
- ⑩ Allen Wrench (only SZO-4)
- ⑪ Power supply

7.4 SZO-5 / SZO-6

- ① Microscope body
(SZO-5 Binocular) / SZO-6 Trinocular)
- ② Eyepieces
- ③ Microscope stand with LED flexible arm
- ④ Stage clips (one pair)
- ⑤ Photo tube (only SZO-6)

- ⑥ Transparent plate
- ⑦ Black/white plate
- ⑧ Dust Cover
- ⑨ Allen Wrench (only SZO-6)
- ⑩ Power supply

7.5 SZO-7 / SZO-8

- ① Microscope body
(SZO-7 Binocular) / SZO-8 Trinocular)
- ② Eyepieces
- ③ Microscope stand SZ-STL1

- ④ Photo tube (only SZO-8)
- ⑤ Dust Cover
- ⑥ Allen Wrench (only SZO-8)

7.6 SZO-9 / SZO-10

- ① Microscope body
(SZO-9 Binocular) / SZO-10 Trinocular)
- ② Eyepieces
- ③ Microscope stand SZ-STL2

- ④ Photo tube (only SZO-10)
- ⑤ Dust Cover
- ⑥ Allen Wrench (only SZO-10)

7.7 Assembly procedure (SZO-B / SZO-T)


1. Remove the eyepieces dust caps and insert the eyepieces in the empty eyepiece sleeve. (Fig. 1)

2. Lock the eyepieces using the locking knob ①. (Fig. 2)

3. Insert the microscope into the head holder and lock the clamping knob ②. (Fig. 3)

7.8 Assembly procedure (SZO-1 / SZO-10)


1. Insert the microscope into the head holder and lock the clamping knob ②. (Fig. 3)
2. Remove the eyepieces dust caps and insert the eyepieces in the empty eyepiece sleeve. (Fig. 1)
3. Lock the eyepieces using the locking knob ①. (Fig. 2)

Only for SZO-3/4/5/6

4. Connect the jack of the power supply to the socket in the rear part of the microscope base. (Fig. 4)

7.9 Assembly procedure (SZO-7 / SZO-8)

1. Screw the pillar on the base. (Fig. 5)

2. Tighten the screw to lock the pillar. (Fig. 6)

3. Insert the drop preventing ring. (Fig. 7)

4. Insert the horizontal arm and lock it with the locking knob. (Fig. 8 - 9)

5. Install the focus holder. Unscrew the prevention knob ④, and insert the head holder arm ③ from above in the hole of the horizontal arm. (Fig. 10-11)


6. Lock the fixing knob ⑤ (Fig. 11)

7. Install the head holder, lock the fixing knob ⑥ and screw the prevention knob ④ from below. (Fig. 12-13)
8. Insert the head into the head holder and lock the clamping knob ②. (Fig. 3)

7.10 Assembly procedure (SZO-9 / SZO-10)

1. Screw the pillar on the base. (Fig. 14)


2. Tighten the screw to lock the pillar. (Fig. 15)


3. Insert the drop preventing ring. (Fig. 16)

4. Insert the horizontal arm and lock it with the locking knob ① (Fig. 17). For a bigger safety, lock the fixing screw ② with the provided allen wrench. (Fig. 18)

5. Install the head holder. Unscrew the prevention knob ③ and insert the head holder from below in the vertical arm. Lock the fixing knob ④. At the end screw again the prevention knob ③. (Fig. 19-20)
6. Insert the head into the head holder and lock the clamping knob. (Fig. 3)

8. Using the microscope


8.1 Adjusting interpupillary distance

Hold the right and left eyepiece tube with both hands and adjust the interpupillary distance by moving the two parts until one circle of light can be seen. If two circles appear, the interpupillary distance is too big, and if two overlapped circles appear, the interpupillary distance is too small. (Fig. 21)

8.2 Focus


Put the sample to be observed on the stage plate and focus the sample using the focusing knobs of the stand ①. If necessary adjust the height of the microscope head along the vertical stand. Remember to lock the lock screw ② after aligning the height of the microscope. (Fig. 22)

8.3 Adjusting the tension of the focus knob

- **This adjustment allows to increase or decrease the tension of the knob by avoiding an involuntary descent of the microscope body under its own weight. Adjust the tension just above the point where the focus is stable.**


Grab the knobs ③ with both hands and, holding right knob rotate the left one. Tension increases or decreases according the sense of rotation of the left focus knob.

8.4 Dioptic compensation

This compensation makes possible for people wearing glasses to adjust the microscope to their eyes and use the microscope without glasses.

1. Put the zoom down to the lowest magnification ④ and focus the specimen with the focusing knobs ①. (Fig. 24)
2. Put the zoom to the maximum magnification and repeat the focusing.
3. Return to the lowest magnification: the specimen will be out of focus.
4. Adjust the diopter compensation ring of the right eyepiece until the image of the right eyepiece is clear and sharp. Repeat the procedure for the left eyepiece.
5. Check the focus of the image for the whole zoom range. It should be perfectly parfocal (focus is maintained during the change of magnification).

8.5 Magnification

Select the desired magnification by adjusting the zoom knob. (Fig. 25)

Change the eyepieces and/or add an appropriate additional lens if necessary

The microscope body is equipped with a “click-stop” function that allows to obtain a precise setting of the desired magnification.

Click stop can be activated or deactivated by operating with the provided allen wrench in the hole placed in the right side of the microscope body.

Fig. 25

8.6 Use of additional lens

Screw the desired additional lens on the microscope body. (Fig. 26)

Each additional lens has a specific Working Distance. The stroke of the focus adapter cannot compensate the different working distances of the several additional lenses.

Should the focus adapter not be able to focus properly the specimen, move the whole microscope body accordingly.

Total magnification used can be calculated as:
 Eyepiece magnification * Zoom magnification * Objective lens magnification.

Fig. 26

| Eyepiece | 10x | | 15x | | 20x | | 25x | |
|-------------------|--------------|--------------|----------------|------------|--------------|------------|----------------|------------|
| Field number (mm) | 23 | | 16 | | 12 | | 9 | |
| Objective | Total Mag | F.O.V (mm) | Total Mag | F.O.V (mm) | Total Mag | F.O.V (mm) | Total Mag | F.O.V (mm) |
| 0.3X | 2.01X-13.5X | 114.42-17.04 | 3.02X-20.25X | 52.98-7.90 | 4.02X-27X | 29.85-4.44 | 5.02X-33.75X | 17.93-2.67 |
| 0.5X | 3.35X-22.5X | 68.66-10.22 | 5.02X-33.75X | 31.87-4.74 | 6.7X-45X | 17.91-2.67 | 8.37X-56.25X | 10.75-1.6 |
| 0.75X | 5.02X-33.75X | 45.82-6.81 | 7.54X- 0.62X | 21.22-3.16 | 10.05X-67.5X | 11.94-1.78 | 12.56X-84.37X | 7.17-1.07 |
| 1X | 6.7X-45X | 34.32-5.11 | 10.05X-67.5X | 15.24-2.37 | 13.4X-90X | 8.96-1.33 | 16.75X-112.5X | 5.37-0.80 |
| 1.5X | 10.05X-67.5X | 22.89-3.41 | 15.07X-101.25X | 10.62-1.58 | 20.1X-135X | 5.97-0.89 | 25.12X-168.75X | 3.58-0.53 |
| 2X | 13.4X-90X | 7.16-2.56 | 20.1X-135X | 7.96-1.19 | 26.8X-180X | 4.48-0.67 | 33.5X-225X | 2.69-0.40 |

8.7 Contrast disc (SZO-1/2/3/4/5/6)

You can use the black/white disc in order to enhance the contrast of the image when using the incident illumination. (Fig. 27).

If you want to look at bright opaque objects, place the disc with the black side up.

Fig. 27

8.8 Video capturing (SZO-2/4/6/8/10)

Two kinds of observation heads, binocular or trinocular, are available for the SZO series. Trinocular head can be connected to cameras via an adapter. (Fig. 28)

Please refer to the adapter and camera manuals for further details.


8.9 Use of the illumination (SZO-3/4/5/6)

Turn the incident light knob ① in order to turn ON/OFF or to change the intensity of the incident light LED.

Turn the transmitted light knob ② in order to turn ON/OFF or to change the intensity of the transmitted light LED.

It is possible to use both illuminations at the same time.

8.10 Use of the overhanging stand (SZO-7/8)

Moving the horizontal arm

Unlock the knob on the right side of the horizontal arm ③.

The arm can be extended or shortened according to specific needs. (Fig. 31)

Swivel the head

Unlock the fixing knob ④ and rotate the head to the desired swivel angle (left or right), then fix the knob again. (Fig. 32)

Rotating the horizontal arm

Loosen the horizontal arm fixing knob ⑤ and rotate the arm, then fix again the fixing knob. (Fig. 33)

- **NOTE: 180° rotation of the microscope with respect to the base could cause a rollover of the entire system.**

8.11 Use of the overhanging stand (SZO-9/10)


Moving the horizontal arm

Unlock the knob on the right side of the horizontal arm ⑥. (Fig. 34)
The arm can be extended or shortened according to specific needs. (Fig. 35)


Tilting the head holder

Loosen the knob ⑦ (Fig. 36) on the right side of the head holder. The tilting vertical arm can be moved. Once the desired angle is achieved, tighten the knob again.

Swivel the head

Unlock the fixing knob ⑧ and rotate the head to the desired swivel angle (left or right), then fix the knob again. (Fig. 37 - 38)

Rotating the horizontal arm

Loosen the horizontal arm fixing knob ⑨ and rotate the arm, then fix again the knob. (Fig. 39)

- **NOTE: 180° rotation of the microscope with respect to the base could cause a rollover of the entire system.**

9. Maintenance

Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 75 % (non condensing). Use a dehumidifier if needed.

To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

For the best results, use the OPTIKA cleaning kit (see catalogue).

If you need to send back the microscope to Optika for maintenance, please use the original packaging.

10. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

| PROBLEM | CAUSE | SOLUTION |
|---|--|--|
| I. Optical Section: | | |
| The illumination is ON, but the field of view is dark. | The plug is not connected to the illumination | Connect the cable |
| | The brightness is too low | Adjust to a proper setting |
| The edge of the field of view is vignetted or the brightness is asymmetric. | The incident illuminator is not correctly oriented | Change the angle of the incident illuminator |
| Dust and stains can be seen in the field of view. | There are stains and dust on the specimen | Clean the specimen |
| | There are stains and dust on the eyepiece | Clean the eyepiece |
| Poor image quality. <ul style="list-style-type: none"> The image is not sharp. The contrast is not high. The details are not clear Image flares | The lenses (additional lens, objective, eyepieces) are dirt | Thoroughly clean all the optical system |
| One side of the image is out of focus. | The specimen is out of place (tilted) | Place the specimen flat on the stage. |
| II. Mechanical Section: | | |
| The focus knob is hard to turn. | The tension adjustment collar is too tight | Loosen the tension adjustment collar |
| The focus is unstable. | The tension adjustment collar is too loose | Tighten the tension adjustment collar |
| III. Electric section: | | |
| The LED doesn't turn on. | No power supply | Check the power cord connection |
| The brightness is not enough | The brightness adjustment is low | Adjust the brightness |
| The light blinks | The power cord is poorly connected | Check the power cord |
| IV. Viewing tube assembly: | | |
| The field of view of the two eyes is different | The interpupillar distance is not correct | Adjust the interpupillar distance |
| | The dioptric correction is not right | Adjust the dioptric correction |
| | The viewing technique is not correct, and the operator is straining the eyesight | When look into the eyepieces, do not stare at the specimen but look at the whole field of view. Periodically, move the eyes away to look at a distant object, then back into the eyepieces |
| V. Microphotography and video: | | |
| The image is unfocussed | Incorrect focussing | Adjusting the focus system as in the present manual |
| Bright patches appear on the image | Stray light is entering the microscope through the eyepieces and through the camera viewfinder | Cover the eyepieces and the viewfinder with a dark cloth |

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.