SZP Series

INSTRUCTION MANUAL

Model
SZP-FL

Ver. 1.1 2020

Contents

1. Warning	2
2. Symbols and conventions	2
3. Safety Information	2
4. Intended use	2
5. Overview	3
6. Unpacking	4
7. Assembling	4
7.1 Assembling procedure	5
8. Use of the microscope	7
8.1 HBO bulb switch on	7
8.2 Centering of HBO bulb	8
Equipment disposal	9

1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use uses that does not comply with this manual.

2. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.

3. Safety Information

Avoiding Electrical Shock


Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

4. Intended use

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

5. Overview

This manual is only for fluorescence attachment.
For information on the microscope, consult the manual provided with the microscope.

6. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

7. Assembling

Once opened the box, the microscope parts are the following:

- ① Fluorescence attachment
- ② HBO lamp housing
- ③ Fluorescence power supply + power cord
- ④ HBO burner


- ⑤ Centering slide
- ⑥ Trino port tube
- ⑦ Filter replacement tool
- ⑧ Allen wrenches

7.1 Assembling procedure

1. Install the fluorescence attachment on the microscope body, aligning the guide under the attachment with the pin on the microscope body.
2. Tight the locking screw ① using the provided Allen wrench. (Fig. 1)


3. Loosen the two fixing screws on the lamp housing port and remove the dust cover.
4. Install the lamp housing and tight the fixing screws. (Fig. 2)

5. Open the lamp housing using the door lock screw ③ and remove the lamp holder. (Fig. 3)

6. Remove the plastic block ④ from the lamp holder (or the exhausted lamp in case of replacement) by loosening the two locking screws ⑤. (Fig. 4)


7. Insert the mercury bulb ⑥ (respect the polarity of the bulb), tighten the locking screws and refit the lamp holder inside the lamp. (Fig. 5)

8. Connect the cables to the external fluorescence power supply. (Fig. 6-7)

9. Install the observation head aligning the guide under the head with the pin on the microscope body. (Fig. 8)

8. Use of the microscope

8.1 HBO bulb switch on

1. Turn on the power supply using the main switch. Wait until the current display shows about 4,5A. (Fig. 9). If the current falls below 4 A, replace the lamp. You should now wait for at least 10 minutes before aligning and using the bulb.

2. Move the shutter pulling toward the user to open the light path. (Fig. 10)


3. Select the filter type moving the filter selector to the desired position. (Fig. 11)


- Fluorescence illuminator has 4 positions. Two are used for the preinstalled filters. It is possible to install additional filter dedicated to specific needs. Contact OPTIKA for filter selection and installation.

8.2 Centering of HBO bulb

1. Put the provided centering slide under the microscope.
2. Operating on the focusing screw of the collector lens ① and on the centering screws ② try to obtain the bright spot of bulb arc. (Fig. 12)

3. Using the focusing screw of the collector lens ① project the arc image on the centering slide. (Fig. 13 A).
 4. With centering screws ② place the arc image at the center of the centering slide. (Fig. 13 B-C)
 5. Using the focusing screw of the collector lens ① focus the arc image. (Fig. 13 D)
- In caso of bulb replacement, the centering procedure must be repeated.

6. After replacing the exhausted bulb, reset the time counter on the power supply by pressing the "Reset" button ③. (Fig. 14)

Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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OPTIKA® S.r.l.