OPTISTAR By OPTIKA Series

# INSTRUCTION MANUAL

Model
OS-A1
OS-A2
OS-A3
OS-B2
OS-B3
OS-C3

Ver. 1.0 2021

---

## Table of contents

1.	Warning	3
2.	Safety Information	3
3.	Unpacking	3
4.	Intended use	3
5.	Symbols and conventions	3
6.	Instrument description	4
6.1	OS-A1	4
6.2	OS-A2	5
6.3	OS-A3	6
6.4	OS-B2	7
6.5	OS-B3	8
6.6	OS-C3	9
7.	Assembling	10
7.1	All models	10
7.1.1	OS-A / OS-C Series	10
7.1.2	OS-B Series	10
7.2	OS-B3	10
7.3	OS-C3	10
8.1	Switching on the microscope	11
8.1.1	OS-A1	11
8.1.2	OS-A2 / OS-A3	11
8.1.3	OS-B2 / OS-B3	11
8.1.4	OS-C3	11
8.2	Light intensity adjustment	12
8.2.1	OS-A1	12
8.2.2	OS-A2 / OS-A3	12
8.2.3	OS-B2 / OS-B3 / OS-C3	12
8.3	Coarse focus tension adjustment	12
8.3.1	OS-B2 / OS-B3	12
8.3.2	OS-C3	13
8.4	Stage	13
8.4.1	OS-A1	13
8.4.2	All other models	13
8.5	Adjust the interpupillary distance	13
8.5.1	OS-A3	13
8.5.2	OS-B / OS-C	14
8.6	Diopter adjustment	14
8.6.1	OS-B Series	14
8.6.2	OS-C Series	14
8.7	Use of oil immersion objective	15
8.8	Condenser centering (OS-B Series only)	15
8.9	Aperture diaphragm	15
8.9	Use of color filter (OS-B Series only)	16
9.	Microphotography	17
10.	Maintenance	18
11.	Troubleshooting	19
	Equipment disposal	20

## 1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

## 2. Safety Information

### Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

## 3. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

## 4. Intended use

### Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

## 5. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

### CAUTION

This symbol indicates a potential risk and alerts you to proceed with caution.

### ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.


---

## 6. Instrument description


### 6.1 OS-A1


6.2 OS-A2

6.3 OS-A3


6.4 OS-B2


6.5 OS-B3

6.6 OS-C3

## 7. Assembling

### 7.1 All models

1. Remove the dust protection cap from the eyepiece sleeve and insert the eyepiece into the empty eyepiece sleeve. (Fig. 1)


#### 7.1.1 OS-A / OS-C Series

Insert the power supply jack in the socket ① placed at the rear side of the microscope. (Fig. 2)

#### 7.1.2 OS-B Series

Insert the power cord in the socket placed at the rear side of the microscope.

### 7.2 OS-B3

1. Remove the dust protection cap of the photo tube.
  2. Insert the "Eyepiece Linker" ②. (Fig. 3)
  3. Insert the projection lens in the empty hole of the Eyepiece Linker.
- **The camera will not be in focus if the Eyepiece Linker is not used.**

### 7.3 OS-C3

1. Unscrew the locking screw ③ on the right side of the photo port and remove the dust protection cap. (Fig. 4)
2. Insert the photo tube ④ in the empty hole and lock the screw.
3. Insert the projection lens in the empty hole of the photo tube.

## 8. Use of the microscope

### 8.1 Switching on the microscope

#### 8.1.1 OS-A1

Operate on the main switch ① placed in the rear side of the microscope, moving the selector on "I". (Fig. 5)

Fig. 5

#### 8.1.2 OS-A2 / OS-A3

1. Operate on the main switch/light intensity dial ② placed on the left side of the microscope, to turn on transmitted light. (Fig. 6)
2. Operate on the main switch/light intensity dial ③ placed on the right side of the microscope, to turn on incident light.

Fig. 6

#### 8.1.3 OS-B2 / OS-B3

Operate on the main switch ④ placed in the rear side of the microscope, moving the selector on "I". (Fig. 7)

Fig. 7

#### 8.1.4 OS-C3

Operate on the main switch ⑤ placed on the right side of the microscope, moving the selector on "I". (Fig. 8)

Fig. 8

## 8.2 Light intensity adjustment

### 8.2.1 OS-A1

Operate on the light intensity dial ① to increase or decrease the illumination intensity. (Fig. 9)

Fig. 9

### 8.2.2 OS-A2 / OS-A3

1. Operate on the main switch/light intensity dial ② placed on the left side of the microscope, to increase or decrease the transmitted light illumination intensity. (Fig. 10)
2. Operate on the main switch/light intensity dial ③ placed on the right side of the microscope, to increase or decrease the incident light illumination intensity.

Fig. 10

### 8.2.3 OS-B2 / OS-B3 / OS-C3

Operate on the light intensity dial ④ to increase or decrease the illumination intensity (Fig. 11)

Fig. 11

## 8.3 Coarse focus tension adjustment

### 8.3.1 OS-B2 / OS-B3

The coarse knob tension is preset in the factory.

1. To modify the tension according to personal's needs, rotate the ring ①. (Fig. 12)
- Clockwise rotation increases the tension.
  - If the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.

Fig. 12

### 8.3.2 OS-C3

The coarse knob tension is preset in the factory.

1. To modify the tension according to personal's needs, rotate the ring ① using a flat screwdriver. (Fig. 13)
- Clockwise rotation increases the tension.
  - If the tension is too loose, the stage could go lower by itself or the focus easily lost after fine adjustment. In this case, rotate the knob in order to increase the tension.


Fig. 13

## 8.4 Stage

### 8.4.1 OS-A1

Put the slide on the stage and hold it using the stage clips ②. (Fig. 14)

Fig. 14

### 8.4.2 All other models

Stage accepts standard slides 26 x 76 mm, thickness 1,2 mm with coverslide 0,17 mm. (Fig. 15)

1. Open the spring arm of the slide holder ③ and place the slide from the front on the stage.
  2. Gently release the spring arm of the slide holder.
- **A sudden release of the spring arm could cause the falling of the slide.**

Fig. 15

## 8.5 Adjust the interpupillary distance

### 8.5.1 OS-A3

Rotate toward left or right the interpupillary distance knob ④ until one circle of light can be seen. (Fig. 16)

The range of the interpupillary distance is 55- 75 mm.

Fig. 16

### 8.5.2 OS-B / OS-C

1. Hold the right and left parts of the observation head using both hands and adjust the interpupillary distance by turning the two parts until one circle of light can be seen. (Fig. 17)
  - The graduation on the interpupillary distance indicator ①, pointed by the spot “.” on the eyepiece holder, shows the distance between the operator’s eyes.

The range of the interpupillary distance is 55- 75 mm.

## 8.6 Diopter adjustment

### 8.6.1 OS-B Series

1. Look into the left eyepiece with your left eye only, and focus on the specimen.
  2. Look into the right eyepiece with your right eye only. If the image is not sharp, use the diopter adjustment ring ② to compensate. (Fig. 18)
- **The adjustment range is  $\pm 5$  diopter. The number indicated on the adjustment ring graduation should correspond to the operator’s diopter correction.**

### 8.6.2 OS-C Series

1. Look into the right eyepiece with your right eye only, and focus on the specimen.
  2. Look into the left eyepiece with your left eye only. If the image is not sharp, use the diopter adjustment ring ③ to compensate. (Fig. 19)
- **The adjustment range is  $\pm 5$  diopter. The number indicated on the adjustment ring graduation should correspond to the operator’s diopter correction.**

## 8.7 Use of oil immersion objective

- **OS-B / OS-C Series**
- 1. Focus the specimen with a low power objective.
- 2. Lower the stage.
- 3. Put a drop of oil (provided) on the area of the specimen to be observed. (Fig. 20)
- **Make sure that there are no oil bubbles. Air bubbles in the oil damage the image quality.**
- To check for bubbles: remove an eyepiece, fully open the aperture diaphragm and observe the objective exit pupil. (The pupil must be circular and bright).
- 4. Insert immersion objective.
- 5. Return the stage to the upper focusing point and obtain an optimal focus using the fine focus knob.
- To remove the air bubbles, gently move the nosepiece to the right and left to move the immersion objective a few times and allow the air bubbles to move away.
- 6. After use, gently remove the oil with a soft paper towel or a lightly moistened optic paper with a mixture of ethyl ether (70%) and absolute ethyl alcohol (30%).
- **The immersion oil, if not immediately cleaned, could crystallize creating a glass-like layer.**
- **In this situation the observation of the specimen would be difficult (even not impossible) due to the presence of an additional thickness on the objective.**


## 8.8 Condenser centering (OS-B Series only)

1. Focus the specimen with 4x objective.
2. Fully close the aperture diaphragm by rotating the ring ①. (Fig. 21)
3. If the condenser is not perfectly centered, move the screws ② until the diaphragm image is in the center of the field of view.
4. Fully open the diaphragm.

## 8.9 Aperture diaphragm

- The Numerical Aperture (N.A.) value of the aperture diaphragm affects the image contrast. Increasing or reducing this value one can vary resolution, contrast and depth of focus of the image. Rotate the diaphragm ring toward left or right until a good contrast is achieved.
- With low contrast specimens set the numerical aperture to about 70%-80% of the objective's N.A. If necessary, remove on eyepiece and, looking into empty sleeve, adjust the condenser's diaphragm in order to obtain an image like the one in Fig. 22.

### 8.9 Use of color filter (OS-B Series only)

Rotate the filter holder ① (Fig. 23) to insert a suitable color filter for the specimen under observation.

---

## 9. Microphotography

- Refer to specific camera Instruction manual for detailed explanation about microphotography procedure.
- Instruction manual can be downloaded from this link: <https://www.optikamicroscopes.com/optikamicroscopes/download/optistar-by-optika-digital-cameras/>

---

## 10. Maintenance

### Microscopy environment

This microscope is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 75 % (non condensing). Use a dehumidifier if needed.

### To think about when and after using the microscope

- The microscope should always be kept vertically when moving it and be careful so that no moving parts, such as the eyepieces, fall out.
- Never mishandle or impose unnecessary force on the microscope.
- Never attempt to service the microscope yourself.
- After use, turn off the light immediately, cover the microscope with the provided dust-cover, and keep it in a dry and clean place.

### Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

### Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Do not disassemble objectives or eyepieces in attempt to clean them.

**For the best results, use the dedicated cleaning kit.**

## 11. Troubleshooting

Review the information in the table below to troubleshoot operating problems.

PROBLEM	CAUSE	SOLUTION
<b>I. Optical Section:</b>		
LED operates, but field of view remains dark.	Power supply is unplugged.	Connect
	Brightness is too low	Set brightness to a proper level
Dirt or dust is visible in the field of view.	Dirt/dust on the specimen	Clean the specimen
	Dirt/dust on the eyepieces	Clean the eyepieces
Image looks double	Aperture diaphragm is stopped down too far	Open aperture diaphragm
Visibility is poor. <ul style="list-style-type: none"> <li>Image is not good.</li> <li>Contrast is poor.</li> <li>Details are unsharp.</li> <li>Image glares</li> </ul>	Revolving nosepiece is in an incorrect position	Move the nosepiece to a click stop
	Aperture diaphragm is too closed or too open	Adjust aperture diaphragm
	Dust or dirt on lenses (condenser, objectives, eyepieces and slide)	Clean thoroughly
	For transmitted light observation, the coverglass thickness must not exceed 0.17mm	Use a coverglass with thickness 0.17mm
	Focus is not even	Slide holder is not flat. Move the specimen to a flat position
One side of the image is out of focus.	The nosepiece is not in the center of the light path	Turn the nosepiece to a click stop
	The specimen is out of place (tilted)	Place the specimen flat on the stage.
	The optical performance of the sample cover glass is poor	Use a cover glass of better quality
<b>II. Mechanical Section:</b>		
The coarse focus knob is hard to turn.	The tension adjustment collar is too tight	Loosen the tension adjustment collar
The focus is unstable.	The tension adjustment collar is too loose	Tighten the tension adjustment collar
<b>III. Electric section:</b>		
The LED doesn't turn on.	No power supply	Check the power cord connection
The brightness is not enough	The brightness adjustment is low	Adjust the brightness
The light blinks	The power cord is poorly connected	Check the power cord
<b>IV. Observation tube:</b>		
Field of view of one eye does not match that of the other.	Interpupillary distance is incorrect.	Adjust interpupillary distance.
	Incorrect diopter adjustment.	Adjust diopter.
	Your view is not accustomed to microscope observation.	Upon looking into eyepieces, try looking at overall field before concentrating on specimen range. You may also find it helpful to look up and into distance for a moment before looking back into microscope.
<b>V. Microphotography:</b>		
Image edge is unfocused	To a certain extent it is due to achromatic objectives features	To minimize the problem, set the aperture diaphragm in a proper position
Bright spots appear on the image	Stray light entering in the microscope through eyepieces or camera viewfinder.	Cover eyepieces and viewfinder with a dark cloth

---

## Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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