DIGITAL IMAGING UNIT

INSTRUCTION MANUAL

Model
TB-3AW
TB-5W

Ver. 2.9     2026

---

# Table of Contents

1.	Warning	3
2.	Safety Information	3
3.	Package content	4
4.	Unpacking	5
5.	Intended use	5
6.	Symbols and conventions	5
7.	Instrument description	6
8.	Assembling	7
8.1	System assembling	7
9.	Use of camera and software	9
10.	Micrometric Slide M-005	9
11.	Maintenance	10
	Equipment disposal	11

---

## 1. Warning

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use. We remind you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users. We decline any responsibility deriving from incorrect instrument use that does not comply with this manual.

## 2. Safety Information

### Avoiding Electrical Shock

Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off position. Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users have full responsibility to use this equipment safely. Please follow the guidelines below, and read this manual in its entirety to ensure safe operation of the unit.

### 3. Package content

① Camera

- TB-3AW: 3MP
- TB-5W: 5MP

② Projection lens

③ Digital imaging unit

④ Digital imaging unit power supply

⑤ USB cable 0,5 m

⑥ Micrometric slide

**NOTE: OPTIKA reserves the right to make corrections, modifications, enhancements, improvements and other changes to its products at any time without notice.**

---

## 4. Unpacking

The microscope is housed in a moulded Styrofoam container. Remove the tape from the edge of the container and lift the top half of the container. Take some care to avoid that the optical items (objectives and eyepieces) fall out and get damaged. Using both hands (one around the arm and one around the base), lift the microscope from the container and put it on a stable desk.

Do not touch with bare hands optical surfaces such as lenses, filters or glasses. Traces of grease or other residuals may deteriorate the final image quality and corrode the optics surface in a short time.

## 5. Intended use

### Standard models

For research and teaching use only. Not intended for any animal or human therapeutic or diagnostic use.

### IVD Models

Also for diagnostic use, aimed at obtaining information on the physiological or pathological situation of the subject.

## 6. Symbols and conventions

The following chart is an illustrated glossary of the symbols that are used in this manual.

### CAUTION


This symbol indicates a potential risk and alerts you to proceed with caution.


### ELECTRICAL SHOCK

This symbol indicates a risk of electrical shock.


7. Instrument description

Left side


Right side

## 8. Assembling

### 8.1 System assembling

1. Remove the dust cover from camera and from projection lens.
2. Screw the projection lens to the camera's thread. (Fig. 1)

3. Insert the bottom part of the projection lens into the empty hole of the photo/tv exit. (Fig. 2)

4. Mount the digital imaging unit. Insert the pins ① placed on the camera bracket into the holes ② placed on the digital imaging unit bracket and gently push down until the connectors are properly locked. (Fig. 3-4-5)
- You can select one of three sets of holes to position the digital imaging unit at the desired height.
  - To remove the digital imaging unit from the bracket, gently push up and pull the connectors out of the bracket.

5. Plug one side of the USB cable ① to the digital head and the other side to the digital imaging unit PC using the connector ②. (Fig. 6-7)
6. Plug the power supply cable to the digital imaging unit PC for battery recharge using the USB-C connector ③. (Fig. 8)

- **NOTE: The tablet comes with a 12V/2A power adapter featuring a European-style plug and a USB-C plug. In regions where a different type of wall outlet is used, the user must use a USB-A/USB-C cable and a Power Delivery adapter with the following maximum specifications:**  
**Input 100/240V – 50/60Hz – 1.7A**  
**Output 5/20V – 6.2/3.25A (Max 67W)**

- The digital imaging unit has been set with the Rotation function disabled: this prevents any flipping of the Live View in order to get a continuous and as large as possible view of your slide also when the digital imaging unit is removed from the holder.
- To enable this function again: you can activate the Rotation by swiping the screen on his bottom right side and selecting Settings + Screen. Anyway, it's not suggested to activate the function when the camera is in Live View mode as it may give troubles when the camera runs at high resolutions.

Fig. 5

Fig. 6


Fig. 7


Fig. 8

## 9. Use of camera and software

1. Tap the icon for the ProView software (or double-click the icon with the mouse). The software starts.
  2. In the "Camera List" panel, the "C-B3A" or "C-B5" items are shown.
  3. Tap on the "C-B3A" or "C-B5" item (or click with the mouse): the live image will be displayed on the main window of the software.
  4. Adjust the camera parameters by acting on the exposure time ("Exposure and Gain" panel) and the white balance ("White Balance" panel).
  5. Once the first adjustments have been made, you can operate normally.
- The software's User Manual is available in PDF format within the application itself and can be opened using the "F1" function key. The manual contains all the operating instructions for using the camera and for the various functions of the software.
  - You must have Acrobat Reader installed to view the manual.


## 10. Micrometric Slide M-005

**Micrometric slide, 26x76mm, with 2 scales**  
**(1mm/100div. for biological microscopes / 10mm/100div. for stereo microscopes)**

1 DIV=0.01mm

**For biological microscopes calibration**

1 DIV=0.1mm

**For stereo microscopes calibration**

---

## 11. Maintenance

### Operative environment

This system is recommended to be used in a clean, dry and shock free environment with a temperature of 5°-40°C and a maximum relative humidity of 75 % (non condensing). Use a dehumidifier if needed.

### To think about when and after using the system

- The system should always be kept vertically when moving it.
- Never mishandle or impose unnecessary force on the system.
- Never attempt to service the system yourself.
- After use, turn off the light immediately, cover the system with the provided dust cover, and keep it in a dry and clean place.

### Electrical safety precautions

- Before plugging in the power supply, make sure that the supplying voltage of your region matches with the operation voltage of the equipment and that the lamp switch is in off-position.
- Users should observe all safety regulations of the region. The equipment has acquired the CE safety label. However, users do have full responsibility to use this equipment safely.

### Cleaning the optics

- If the optical parts need to be cleaned try first to: use compressed air.
- If that is not sufficient: use a soft lint-free piece of cloth with water and a mild detergent.
- And as a final option: use the piece of cloth moistened with a 3:7 mixture of ethanol and ether.
- **Note: ethanol and ether are highly flammable liquids. Do not use them near a heat source, near sparks or near electric equipment. Use these chemicals in a well ventilated room.**
- Remember to never wipe the surface of any optical items with your hands.
- Fingerprints can damage the optics.


**For the best results, use the OPTIKA cleaning kit (see catalogue).**

If you need to send the systema to Optika for maintenance, please use the original packaging.

---

## Equipment disposal

Art.13 Dlsg 25 July 2005 N°151. "According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal."

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste. The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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**OPTIKA® S.r.l.**