
OPTIKA

M I C R O S C O P E S
I T A L Y

Ver. 3.0.0

XZ-2

OPERATION MANUAL

GUIDA UTENTE

MANUAL DE INSTRUCCIONES

OPTIKA MICROSCOPES - ITALY

www.optikamicroscopes.com - info@optikamicroscopes.com

1.0 DESCRIPTION	page 4
2.0 INTRODUCTION	page 5
3.0 UNPACKING AND ASSEMBLY	page 5
4.0 USING THE MICROSCOPE	page 15
5.0 MAINTENANCE	page 19
6.0 ELECTRICAL SPECIFICATIONS	page 19
7.0 RECOVERY AND RECYCLING	page 20

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily use.

Optika reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

Optika declines any responsibility deriving from instrument uses that do not comply with this manual.

Safety guidelines

This manual contains important information and warnings regarding safety about installation, use and maintenance of the microscope. Please read this manual carefully before using the equipment. To ensure safe use, the user must read and follow all instructions in this manual. OPTIKA products are designed for safe use in normal operating conditions. The equipment and accessories described in the manual are manufactured and tested according to industry standards for safety instrumentation laboratory. Misuse can cause personal injury or damage to the instrument. Keep this manual at hand close to the instrument, for an easy consultation.

Electrical safety

Before connecting the power cord to wall outlet, ensure that your mains voltage for your region corresponds to the voltage supply of the instrument, and that the illuminator's switch is in position OFF. The user must observe the safety regulations in force in his region. The instrument is equipped with CE safety marking, in any case the user has full responsibility concerning the safe use of that instrument.

Warning/Caution symbols used in this manual

The user should be aware of safety aspects when using the instrument. Warning or hazard symbols are shown below. These symbols are used in this manual.

**DANGER**

The instructions on this symbol to avoid possible severe personal injuries.

**WARNING**

Warning of use; the incorrect operation on the instrument can cause damages to the person or instrument.

**WARNING**

Possibility of electric shock.

**HOT!**


Attention: high temperature surfaces. Avoid direct contact.

**NOTE**

Technical notes or usage tips.

1.0 DESCRIPTION

This microscope is a scientific precision instrument designed to last for many years with a minimum of maintenance. It is built to high optical and mechanical standards and to withstand daily laboratory use.

Optika reminds you that this manual contains important information on safety and maintenance, and that it must therefore be made accessible to the instrument users.

Optika declines any responsibility deriving from instrument uses that do not comply with this manual.

Connect the mains plug into the socket at the base

Make sure, before you turn the illumination on, that the voltage selector is set to the mains voltage for your region.

The power cord should be used only on network sockets equipped with adequate grounding. Contact a technician to check the state of your electrical system. If there is no need to install additional accessories, the instrument is now ready for use. Once positioned and installed with the necessary components, the microscope is ready to be used. Your microscope is a laboratory instrument designed to last. Handle it always carefully and avoid abrupt vibrations or shocks. Always disconnect the power cable from the microscope when not in use for long time, while you clean it or when you perform any maintenance.


AVOID DISASSEMBLING THE INSTRUMENT

Do not disassemble the instrument. This entails the cancellation of the warranty and may cause malfunction.

3.0 UNPACKING AND ASSEMBLY

Components of the microscope

Unpack the components of the microscope from the box and put them on a stable surface. In the packaging you find:

1 - PILLAR BASE

2 - HEAD HOLDER WITH FOCUSING MECHANISM

3 - SUPPORT ROD FOR LCD


4 - OPTICAL HEAD WITH BUILT-IN DIGITAL CAMERA

5 - 8" TFT LCD SCREEN

6 - EYEPIECES WF10X/18MM


7 - LED RING ILLUMINATOR WITH BRIGHTNESS CONTROL

8 - LED RING POWER SUPPLY (OUT: 12VDC 0.5A)

9 - DIGITAL CAMERA POWER SUPPLY (OUT: 5,5VDC 200mA)

10 - LCD SCREEN POWER SUPPLY (OUT: 12VDC 2,6A)


Assembling the microscope


1. Put the optical head in the head holder and then slide it on the pillar. Tighten the screw on the back at about half height of the pillar (you can move it later when you focus your sample).

2. Put the LCD screen on a table with its rear up-facing, then insert the slider of the support rod all the way up and tighten the two screws.

SCREWS


3.0 UNPACKING AND ASSEMBLY

2. Put the support rod with the screen on the pillar base, locking the screw.


USE THIS SCREW TO ADJUST
LCD INCLINATION

LOCK THIS SCREW

4. Put the LED ring illuminator at the bottom of the optical head.

LOCK THIS SCREW


5. Take the cable connected to the digital camera, and plug it to the power supply (5,5V).

DIGITAL CAMERA CABLE

JACK FROM 5,5V POWER SUPPLY

6. Plug the other jack of the digital camera cable into the AV input on the rear of the LCD screen.


3.0 UNPACKING AND ASSEMBLY

7. Plug the jack from LED ring power supply (12Vdc 0,5A) into the brightness control.

JACK FROM 12V-0.5A POWER SUPPLY

8. Plug the jack from the LCD screen power supply (12Vdc 2,6A) into the DC12V input on the rear of the LCD screen.

Now the XZ-2 is ready to be used.

- 4.1** Switch on the power supplies and press the POWER button on the LCD screen.
- 4.2** Use the brightness control to adjust the LED light:
- position 0: LED OFF
 - position I: LED ON with adjustable brightness
 - position II: LED ON at maximum brightness
- 4.3** Put your sample on the base under the objective and rotate the focusing knob to bring the object into focus on the LCD screen.
- 4.4** Rotate the zoom adjustment ring at the desired magnification (zoom range: 0,7x – 4,5x). The ring has click-stop positions at specific zoom values for precise magnifications.
- 4.5** **Focus tension adjustment:** The tension of the focusing knob can be adjusted by rotating the right focusing knob keeping still the left focusing knob.
- 4.6** **Interpupillary distance of the optical head:** move the two eyepiece tubes until only one circular field can be seen through the two eyepieces. If two circles appear, the interpupillary distance is too big, and if 2 overlapped circles appear, the interpupillary distance is too small.
- 4.7** **Dioptric compensation of the optical head:** this compensation makes it possible for people with glasses to adjust the microscope to their eyes and use the microscope without glasses. Turn the zoom down to the lowest magnification. Adjust the diopter compensation ring of the right eyepiece tube until the image of the right eyepiece is clear and sharp. Repeat the procedure for the left eyepiece. Then, check the focus of the image for the whole zoom range. It should now be perfectly parfocal (focus is always maintained during the change of magnification).

- 4.8 Orientation of the digital camera:** in order to have the image on the LCD aligned with the object under observation, you can rotate the digital camera on the top of the optical head. Note that the image in the eyepieces is reversed compared to the LCD image.
- 4.9 Setting the LCD screen.** On the front panel of LCD you see the buttons for LCD management:

- **POWER** (button at right): press this button to switch ON / OFF the screen.

- Brightness adjustment:

Press buttons:

MENU → ► → ► or ◀

To exit press:

MENU

- Contrast adjustment:

Press buttons:

MENU → ► → ▼ → ► or ◀

To exit press:

MENU

- Saturation adjustment:

Press buttons:

MENU → ► → ▼ → ▼ → ► or ◀

To exit press:

MENU

- Language selection:

Press buttons:


MENU → ▼ → ► → ► → ▲ or ▼

To save press:

To exit press:

MENU


- Cross display:

Press buttons:

MENU → ▼ → ▼ → ► → ► → ▲ → ►
→ ▲ or ▼

To save press: ►

To exit press: **MENU**

- Analog camera input settings:

Press button:

to select "AV" for the XZ-2 camera.

4.0 USING THE MICROSCOPE

- Reset the image settings (all values set to 50):

Press buttons:


MENU → ▲ → ▲ → ► → ► → ► → ▲ *or* ▼

To save press:

To exit press:

MENU

- 4.10 Analog input.** Connect an AV coaxial cable to TV connector on the rear, in order to see an external source.

5.1 Always think about

- The following environment is required: Indoor temperature: 0-40°C, Maximum relative humidity: 85 % (non condensing).
- Keep the microscope away from dust and shocks while in use.
- Turn off the light immediately after use.
- Use a soft lens tissue to clean the optics after use.
- Only if needed, use a cloth moistened with water and a mild detergent, rinsing with water and drying immediately with a lint-free cloth.
- After use, cover the microscope with the included dust-cover, and keep it in a dry and clean place.

5.2 Do not!

- Wipe the surface of any optical items with your hands. Fingerprints can damage the optics.
- Use solvents, neither on the microscope, nor on the optics.
- Disassemble objective or eyepieces to attempt to clean them.
- Mishandle or impose unnecessary force on the microscope.
- Clean the unit with volatile solvents or abrasive cleaners.
- Attempt to service the microscope yourself.

- 5.3** If you need to send the microscope to OPTIKA for maintenance, please use the original packaging.

6.0 ELECTRICAL SPECIFICATIONS

Power Supply:

INPUT:	100-230 Vac, 50/60 Hz
OUTPUTS:	LED ring: 12Vdc 0.5A Digital camera: 5,5Vdc 200mA LCD screen: 12Vdc 2,6A

Art.13 Dlsg 25 july 2005 N°151. “According to directives 2002/95/EC, 2002/96/EC and 2003/108/EC relating to the reduction in the use of hazardous substances in electrical and electronic equipment and waste disposal.”

The basket symbol on equipment or on its box indicates that the product at the end of its useful life should be collected separately from other waste.

The separate collection of this equipment at the end of its lifetime is organized and managed by the producer. The user will have to contact the manufacturer and follow the rules that he adopted for end-of-life equipment collection. The collection of the equipment for recycling, treatment and environmentally compatible disposal, helps to prevent possible adverse effects on the environment and health and promotes reuse and/or recycling of materials of the equipment. Improper disposal of the product involves the application of administrative penalties as provided by the laws in force.

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