



144-LED ring illuminator, EU plug

Chromaticity and LED	Coordinates (x,y)	(0.32, 0.32)
	Color Temperature (K)	6,300
	LED number	144
	Brightness control	Manual
	Quadrant ON/OFF control	Yes
	Lifetime (hours)	> 65,000
	Max. required power (W)	2
Optical power	Luminous flux (lumen)	120
	Illuminance (at 10 cm distance, lux)	6,000 (on 10 cm diameter)
Power Supply for Illumination	Type	Internal
	Power plug type	EU plug
	Input voltage	100/240 Vac, 50/60 Hz, 150mA
	Output voltage	10-12Vdc
Accessories Included	Fixing ring	Yes