



LED ring illuminator, EU plug

Chromaticity and LED	Coordinates (x,y)	(0.32, 0.32)
	Color Temperature (K)	6,300
	Brightness control	Manual
	Quadrant ON/OFF control	Yes
	Lifetime (hours)	> 65,000
	Max. required power (W)	2
Optical power	Luminous flux (lumen)	240
	Illuminance (at 10 cm distance, lux)	8,000 (on 3 cm diameter)
Power Supply for Illumination	Type	Internal
	Power plug type	EU plug
	Input voltage	100/240 Vac, 50/60 Hz, 150mA
	Output voltage	10-12Vdc
Accessories Included	Fixing ring	Yes