

**Model:**

Fluorescence attachment with mercury lamp

Description:

Fluorescence with 6-position filterwheel and HBO mercury lamp.

| | |
|---------------------|---|
| Source | HBO mercury vapour lamp. Power: 100 W. Lifetime: about 300 hours. Knobs for lamp and rear mirror centering. |
| Filter Set | 6 positions fluorescence filter holder: Standard mounted filtersets: - B: EX 460-490, DM 500, EM 510LP; - G: EX 510-550, DM 570, EM 590LP; Optional filtersets: - UV: EX: 330-385 , DM 400 , EM 420LP; - V: EX 400-410, DM 455, EM 455LP. Custom filtersets available on request. |
| Diaphragms | Aperture and field diaphragms. |
| Power supply | Included. Voltage: 100/240 Vac, 50/60 Hz Fuse: T3,15 A 250 V. Max power required: 130 W. |
| Dimensions | 620 mm overall length. |