

**Model:**

Rackless mechanical stage

**Description:**

*Double-layer mechanical stage with coaxial translation knobs, with belt-drive movement.*

<b>Dimensions</b>	220x150 mm.
<b>Translation</b>	X: 80 mm. Y: 50 mm. Coaxial drive control length: 137 mm.
<b>Slide holder</b>	For 2 standard slides.
<b>Graduated axis</b>	On both X and Y: 0.1 mm vernier scale.
<b>Notes</b>	Internal belt-drive system on X axis (no protruding parts during movement).