



Model:

STX

Typology:

EDUCATIONAL STEREO MICROSCOPE

Description:

Entry-level stereomicroscope and educational field.

Head	Binocular, vertical, fixed.
Eyepieces	Wide Field, WF 10x/16.
Objective	2x fixed. Working distance: 67 mm.
Stand	Fixed with focus mechanism.
Focusing	By means of a rack and pinion controlled by a pair of knobs placed on both sides of the microscope.
Stage	Fitted with specimen clamps.
Illumination	Incident illuminator fitted with high-efficiency 0.1 W LED light source. Maximum power consumption: 0.15 W. Powered by 2 AA batteries (not included).
Dimensions	HEIGHT: 230 mm WIDTH: 100 mm DEPTH: 135 mm WEIGHT: 487 gr
Accessories included	Manual and dust cover.