



*Fluorescence attachment for SZP series,  
Blue & Green filtersets*

<b>Observation Method - Incident Light</b>	Fluorescence	Yes
<b>Accessories Included</b>	Dust cover	Yes
	Allen wrench	Yes
	User Manual	Digital version (downloadable)
<b>Product Dimensions</b>	Total width (mm)	520
	Total depth (mm)	200
<b>Product Weight</b>	(kg)	6.4 (including FLUO power supply)
<b>Fluorescence Attachment</b>	Number of positions	4
	Blue filter set (included) specs	Excitation: 460 - 500 nm; Dichroic: 505 nm; Emission: 510LP nm
	Green filter set (included) specs	Excitation: 510 - 550 nm; Dichroic: 570 nm; Emission: 575LP nm
	Filter dimensions	Excitation: 18 mm diam.; Dichroic: 27 mm x 18 mm; Emission: 18 mm diam.
	Filter set selection	Manual
	Shutter	Yes
<b>Fluorescence Light Source</b>	Light source	HBO
	Watt (W)	100
	Lifetime (hours)	300
<b>Fluorescence Power Supply</b>	Type	External
	Microscope connector	C13 plug
	Power plug type	Multi-plug (EU, UK, US)
	Input voltage	100/240 Vac, 50/60 Hz
	Max. power required (W)	130
	LED indicators	Yes