

Large stand, fixed arm type

Focusing System	Type	Coarse
	Coarse total travel (mm)	100
	Adjustable tension	Yes
	Head holder internal diameter (mm)	76
Stand	Type	Fixed arm
Incident Illumination	Type	No
	Annular	No
	N° of LED	No
	X-LED type	No
	Light source power (W)	No
	Illuminance (lux)	No
	Brightness control	No
	Lifetime (hours)	No
	Temperature (K)	No
	Max. required power (W)	No
Accessories Included	Dust cover	Yes
	User Manual	Digital version (downloadable)
Product Dimensions	Total height (mm)	320
	Base height (mm)	50
	Base width (mm)	380
	Base depth (mm)	260
Product Weight	(kg)	4,5